

Product datasheet for **SC327081**

PPP2R5C (NM_001161725) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPP2R5C (NM_001161725) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPP2R5C
Synonyms:	B56G; B56gamma; PR61G
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327081 representing NM_001161725.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCGAATAAAAAACAAGAAGGAGAAAAGAATCACAAAAGCAGGGAAGAGTGGAAAAAGTTCAAAAGAA
GGACAAGACACAGTAGAATCAGAGGGAACAAGTCTGAAGAACCCTCAAGCCCTAAAGTACCACCTCCA
CTTCTTCCTGAGTTGCTGGTCTTATTTTTGGTGGGCTTCAAGGAAGAGATGTTCCCTCCTGCTGATCAA
GAGAAGCTTTTTATCCAGAAGTTACGTCAGTGTTGCGTCTCTTTGACTTTGTTTCTGATCCACTAAGT
GACCTAAAGTGGAAAGGAAGTAAAACGAGCTGCTTAAAGTAAATGGTAGAATATATCACCCATAATCGG
AATGTGATCACAGAGCCTATTTACCCAGAAGTAGTCCATATGTTTGCAGTTAACATGTTTCGAACATTA
CCACCTTCCTCCAATCTACGGGAGCGGAATTTGACCCGGAGGAAGATGAACCAACGTTAGAAGCAGCC
TGGCCTCATCTACAGCTTGTATGAATTTTTCTTAAGATTTTTAGAGTCTCCAGATTTCCAACCTAAT
ATAGCGAAGAAAATATATTGATCAGAAGTTGTATTGCAGCTTTAGAGTCTTTGACAGTGAAGATCCT
CGGGAGAGAGATTTTCTAAAACCACCTTCACAGAATCTATGGGAAATTCCTAGGCTTGAGAGCTTAC
ATCAGAAAACAGATAAATAATATATTTTATAGGTTTATTTATGAAACAGAGCATCATAATGGCATAGCA
GAGTTACTGGAAATATTGGGAAGTATAATTAATGGATTTGCCTTACCCTAAAAGAGAGACAAGATT
TTCTTATTGAAGGTGTTACTACCTTTGCACAAAGTAAATCTCTGAGTGTCTACCATCCCCAGCTGGCA
TACTGTGTAGTGCAGTTTTTAGAAAAGGACAGCACCCCTCACGGAACCAAGTGGTATGGCACTTCTCAA
TACTGGCCAAAGACTCACAGTCCAAAAGAAGTAATGTTCTTAAACGAATTAGAAGAGATTTTAGATGTC
ATTGAACCATCAGAAATTTGTGAAGATCATGGAACCCCTTCCGGCAGTTGGCCAAATGTGTCTCCAGC
CCACACTTCCAGGTGGCAGAGCGAGCTCTCTATTACTGGAATAATGAATACATCATGAGTTAATCAGT
GACAACGCAGCGAAGATTCTGCCATCATGTTTCTTGTACCGCAACTCAAAGACCCATTGGAAC
AAGACAATACATGGCTTGATATACAACGCCCTGAAGCTTTCATGGAGATGAACAAAAGCTATTTGAT
GACTGTACACAACAGTTCAAAGCAGAGAACTAAAAGAGAAGCTAAAAATGAAAGAACGGGAAGAAGCA
TGGGTTAAAATAGAAAATCTAGCCAAAGCCAATCCCCAGTACACAGTGTATAGTCAAGCCAGCACCATG
AGCATTCCGGTTGCAATGGAGACAGATGGGCTTTATTTGAAGATGTGCAGATGCTGAGAAAAGACAGTG
AAGGACGAGGCTCATCAGGCACAGAAAGATCCGAAGAAGGACCGTCTCTTGCACGCCGCAAGTCCGAG
CTGCCTCAGGACCCACACCAAGAAAGCCTTGAAGCTCACTGCAGGGCCGATGAGCTGGCCTCCCAG
GACGGCCGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001161725

Insert Size: 1668 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001161725.1](#)

RefSeq Size: 4494 bp

RefSeq ORF: 1668 bp

Locus ID: 5527

UniProt ID: [Q13362](#)

Cytogenetics: 14q32.31

Protein Families: Druggable Genome, Phosphatase

Protein Pathways: Oocyte meiosis, Wnt signaling pathway

MW: 64.3 kDa

Gene Summary: The product of this gene belongs to the phosphatase 2A regulatory subunit B family. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The B regulatory subunit might modulate substrate selectivity and catalytic activity. This gene encodes a gamma isoform of the regulatory subunit B56 subfamily. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (5) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (e) has a distinct N-terminus and is longer than isoform a.