

Product datasheet for **SC327085**

WDR46 (NM_001164267) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WDR46 (NM_001164267) Human Untagged Clone
Tag:	Tag Free
Symbol:	WDR46
Synonyms:	BING4; C6orf11; FP221; UTP7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327085 representing NM_001164267.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGACAGCCCCAAGCCGGGCAAGGATGTCCCGCCCAAGAAAGACAAACTTCAGACCAAGAGAAAG
AAACCCCGAGAATGGAAGAACCCGGAGTCCCAGCGCGGCTTGCCGGGACCCAAGATCCATTCCCAGGC
CCGCCCCCGTCCCTGTGGAAGTGGTCCAGAAGTTCTGTGCGATTGACAAATCCCAGAAAGCTACCACAT
TCTAAAGCCAAAACCTCGAAGCCGACTTGAGGTGGCTGAAGCTGAGGAAGAGGAAACAAAGTATCAAAGCT
GCTCGTTCTGAGCTGCTGCTTGTGAAGAACCTGGGTTTCTGGAAGGGGAGGATGGGGAAGACACAGCA
AAGATATGCCAGGCTGACATTGTGGAGGCTGTGGACATTGCAAGTGCAGCCAAGCACTTTGACTTGAAT
CTGCGGCAGTTTGGACCCTACAGACTAACTACTCTGAACTGGAAGACACCTGGCTTTTGGAGGGCGC
CGAGGTATGTGGCTGCCCTTATTGGGTAAACAAAGAAGCTTATGTGCGAGATCAACGTCATGGAGGCG
GTGCGGGACATCCGTTTCTCCATTCTGAGGCACTGCTTGTCTGCTCAGAACCCTGGCTCCACATC
TATGACAATCAGGGCATTGAGCTCCACTGTATCCGCCGCTGTGACCGAGTAACACGGCTTGAGTTCTTG
CCCTTCCACTTCCTCCTGCTACAGCTTCAAGAACAGGGTTTCTAACCTACCTGGATGTGTCAAGTGGG
AAGATTGTGGCAGCTCTGAATGCTCGAGCTGGGCGGCTCGATGTTATGAGTCAGAACCTTACAATGCC
GTCATCCATCTCGGACACAGCAATGGTACTGTGCTTTATGGAGTCCAGCTATGAAGGAGCCACTGGCA
AAGATTCTCTGTCATCGTGGTGGGTCGGGCTGTGGCAGTAGATTCTACAGGCACGTACATGGCCACC
TCTGGCCTAGACCACAGCTGAAGATCTTTGACTTGCGAGGGACGTACCAGCCTCTGAGCACTCGGACC
CTGCCCATGGAGCAGGGCACCTGGCCTTCTCCAGAGGGGACTGCTGGTGGCGGAATGGGTGACGTT
GTCAACATCTGGCAGGGCAGGGCAAGGCCAGCCACCTCCCTTGAACAGCCCTACCTACCCACCGG
CTCTCAGGCCCTGTGCATGGCCTTCAGTTCTGCCCTTTGAAGATGTGCTGGGGTGGGGCACACTGGG
GGCATCACAGCATGCTGGTCCCTGGGGCCGGTGAGCCCACTTCGATGGCTGGAGAGTAATCCATAC
AGAAGCCGGAAGCAGCGCCAGGAGTGGGAGGTGAAGGCCCTGCTAGAGAAGGTACCTGCAGAGCTTATT
TGTCTGGACCCACGAGCCCTGGCCGAGGTGGATGTCTCTCCCTGGAGCAGGGAAGAAGGAGCAGATA
GAGAGGCTGGGCTATGACCCGAGGCTAAGGCTCCCTTCCAGCCAAAGCCAAAGCAGAAGGGCCGACG
TCCACGGCAAGCCTGGTGAAGAGGAAGGAAGGTATGGATGAGGAACACAGGGACAAGGTCCGGCAG
AGCCTTACGAGCAGCATATAAGGAGGCAAGGCCAAGCCACGGGGGCCCGGCATCTGCCCTGGAC
AGATTTGTGCGCTGA
ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001164267

Insert Size: 1671 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001164267.1</u>
RefSeq Size:	2201 bp
RefSeq ORF:	1671 bp
Locus ID:	9277
UniProt ID:	<u>O15213</u>
Cytogenetics:	6p21.32
Protein Families:	Druggable Genome
MW:	61.8 kDa
Gene Summary:	Scaffold component of the nucleolar structure. Required for localization of DDX21 and NCL to the granular compartment of the nucleolus.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) uses an alternate in-frame splice site in the 5' coding region, compared to variant 1. This results in a shorter protein (isoform 2), compared to isoform 1.