

Product datasheet for **SC327086**

PLBD2 (NM_001159727) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PLBD2 (NM_001159727) Human Untagged Clone
Tag:	Tag Free
Symbol:	PLBD2
Synonyms:	P76
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327086 representing NM_001159727.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTGGGCCAGATGACTGCTACCCCGCAGCCACCTGGCCCGGGCGCTGACGCGGGCGCTGGCGCTG
GCCCTGGTGTGGCCCTGCTGGTCGGGCGTTCTTGAGCGGCCTGGCGGGGCGATCCCAGCGCCGGGG
GGCCGCTGGGCGCGGATGGGCAGGTCCCTCCAGCCTCCCGCAGCCGCTCGGTGCTCCTGGACGCTCG
GCGGGCCAGCTGCTTATGGTGGACGACGCCACCTGACGCCGTGGCCTGGGCCAACCTACCAACGCC
ATCCGCGAGACTGGTGGGCTTCTGGAGCTGGGCACAAGTGGCCAATAACAATGACAGCTTGCAGGCC
TATGCAGCCGGTGTGGTGGAGGCTGCTGTGTCGGAGGAGCTCATCTACATGCACTGGATGAACACGGTG
GTGAATTACTGCGCCCTTCGAGTATGAAGTCGGCTACTGCGAGAGGCTGAAGAGCTTCTGGAGGCC
AACCTAGAGTGGATGCAGGAAGAGATGGAGTCAAACCCAGACTCACCTTACTGGCACCAGGTGCGGCTG
ACCCTCTGCAGCTGAAAGGCCCTGGAGGACAGCTACGAAGCCGTGTGAGCTTCCAGCTGGGAAGTTC
ACCATCAAACCTTGGGGTTCTCTGCTGCAGCTCTCTGGGACCTGGAAGACCTGGAGCTGGCCCTG
AAACAAGACCAAGATCAAACCTTCTCTGGGCTCTGGCTCCTGTTCTGCCCTCATCAAGCTGCTCCCTGGC
CAGAGTGACCTCCTGGTTGCCCAACACCTGGAACAACCTACCAGCACATGCTGCGTGTCAAGAAAG
TACTGGCTCCAGTTCGGGAAGGCCCTGGGGGACTACCCGCTGGTTCCCGGCAACAAGCTGGTCTTC
TCCTCCTACCCCGCACCATCTTCTCTGCGACGACTTCTACATCCTGGGAGTGGGCTGGTGACTG
GAGACCACCATTTGGCAACAAGAACCAGCCCTGTGGAAGTATGTGCGGCCAGGGGCTGTGTGCTGGAG
TGGGTACGCAACATCGTGGCCAACCGCTGGCCTCGGATGGGGCCACCTGGGCAGACATTTCAAGAGG
TTCAACAGCGGCACCGGCATGGTGGTGGCTGACAAGACCTCGGAGCTCTACCAGAAGACCTACTGG
GCCAGCTACAACATACCGTCTTCGAGACTGTGTTCAATGCCAGTGGGCTGCAGGCCCTAGTGGCCAG
TATGGGGACTGGTTTTCTTATGACGGGAGCCCCGGGCCAGATCTTCCGGCGGAACCAAGTCACTGGTA
CAAGACATGGACTCCATGGTCAGGCTGATGAGGTACAATGACTTCTCCATGACCTCTGTCACTGTGC
AAAGCCTGCAACCCAGCCCAATGGGGAGAATGCTATCTCCGCCGCTCCGACCTCAACCCGGCCAAT
GGCTCCTACCCCTTCCAGGCCCTGCGTCAGCGCTCCCATGGGGTATCGATGTGAAGGTGACCAGCATG
TCACTGGCCAGGATCCTGAGCCTGCTGGCGCCAGCGGTCCACGTGGGACCAGGTGCCCCGTTCCAG
TGGAGCACCTCGCCCTCAGCGGCTGCTGCACATGGGCCAGCCAGACCTCTGGAAGTTCGCGCCTGTC
AAGGTTTCATGGGACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001159727

Insert Size: 1674 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001159727.1](#)

RefSeq Size: 2736 bp

RefSeq ORF: 1674 bp

Locus ID: 196463

UniProt ID: [Q8NHP8](#)

Cytogenetics: 12q24.13

Protein Families: Transmembrane

MW: 61.9 kDa

Gene Summary: Putative phospholipase.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) lacks an in-frame exon in the coding region, compared to variant 1. The encoded isoform (2) is shorter than isoform 1. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on the chimeric transcripts BC002929.1 and BC071832.1.