

Product datasheet for **SC327114**

NRG1 (NM_001160001) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NRG1 (NM_001160001) Human Untagged Clone
Tag:	Tag Free
Symbol:	NRG1
Synonyms:	ARIA; GGF; GGF2; HGL; HRG; HRG1; HRGA; MST131; MSTP131; NDF; NRG1-IT2; SMDF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327114 representing NM_001160001.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGAAAGGACGCGCGGGCCGAGTTGGCACCACAGCCTTGCCCTCCCGATTGAAAGAGATGAAAAGC
CAGGAATCGGCTGCAGGTTCCAAACTAGTCCTTCGGTGTGAAACCAGTTCTGAATACTCCTCTCAGA
TTCAAGTGGTTCAAGAATGGGAATGAATTGAATCGAAAAACAAACCACAAAATATCAAGATACAAAA
AAGCCAGGGAAGTCAGAACTTCGCATTAACAAAGCATCACTGGCTGATTCTGGAGAGTATATGTGCAAA
GTGATCAGCAAATTAGGAAATGACAGTGCCCTGCGCAATATCACCATCGTGAATCAAACGCTACATCT
ACATCCACCCTGGGACAAGCCATCTTGTAAAATGTGCGGAGAAGGAGAAAACTTCTGTGTGAATGGA
GGGAGTGCTTCATGGTAAAGACCTTTCAAACCCCTCGAGATACTTGTGCAAGTGCCCAATGAGTTT
ACTGGTATCGTGCCAAAACACGTAATGGCCAGCTTCTACAAGCATCTGGGATTGAATTTATGGAG
GCGGAGGAGCTGTACCAGAAGAGAGTGTGACCATAACCCGCATCTGCATCGCCCTCCTTGTGGTCGGC
ATCATGTGTGGTGGCCTACTGCAAAACCAAGAAACAGCGGAAAAAGCTGCATGACCGTCTTCGGCAG
AGCCTTCGGTCTGAACGAAACAATATGATGAACATTGCCAATGGGCCCTACCATCCTAACCCACCCCC
GAGAATGTCCAGCTGGTGAATCAATACGTATCTAAAAACGTCATCTCCAGTGAGCATATTGTTGAGAGA
GAAGCAGAGACATCCTTTCCACCAGTCACTATACTTCCACAGCCATCACTCCACTACTGTCAACCAG
ACTCCTAGCCACAGCTGGAGCAACGGACACACTGAAAGCATCCTTTCCGAAAGCCACTCTGTAATCGTG
ATGTCATCCGTAGAAAACAGTAGGCACAGCAGCCAACTGGGGGCCAAGAGGACGCTTAAATGGCACA
GGAGGCCCTCGTGAATGTACAGCTTCTCAGGCATGCCAGAGAAACCCCTGATTCTACCGAGACTCT
CCTCATAGTGAAGGTATGTGTCAGCCATGACCACCCCGGCTCGTATGTACCTGTAGATTTCCACAGC
CCAAGCTCCCCAAATCGCCCCCTTCGAAATGTCTCCACCCGTGTCCAGCATGACGGTGTCCATGCCT
TCCATGGCGGTCAGCCCTTTCATGGAAGAAGAGAGACCTCTACTTCTCGTGACACCACCAAGGCTGCGG
GAGAAGAAGTTTGACCATCACCCCTCAGCAGTTCAGCTCCTTCCACCACAACCCCGCGCATGACAGTAAC
AGCCTCCCTGCTAGCCCTTGAGGATAGTGAGGATGAGGAGTATGAAACGACCCAAGAGTACGAGCCA
GCCAAGAGCCTGTTAAGAACTCGCCAATAGCCGGCGGGCCAAAAGAACCAAGCCCAATGGCCACATT
GCTAACAGATTGGAAGTGGACAGCAACACAAGCTCCCAGAGCAGTAACTCAGAGAGTGAACAGAAGAT
GAAAGAGTAGGTGAAGATACGCCTTCTGGGCATACAGAACCCCTGGCAGCCAGTCTTGAGGCAACA
CCTGCCTTCGCCTGGCTGACAGCAGGACTAACCCAGCAGGCCGCTTCTCGACACAGGAAGAAATCCAG
GCCAGGCTGTCTAGTGAATTGCTAACCAAGACCTATTGCTGTATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001160001

Insert Size: 1773 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001160001.2](#)

RefSeq Size: 5518 bp

RefSeq ORF: 1773 bp

Locus ID: 3084

UniProt ID: [Q02297](#)

Cytogenetics: 8p12

Protein Families:	Druggable Genome, Secreted Protein, Transcription Factors, Transmembrane
Protein Pathways:	ErbB signaling pathway
MW:	65.5 kDa
Gene Summary:	<p>The protein encoded by this gene is a membrane glycoprotein that mediates cell-cell signaling and plays a critical role in the growth and development of multiple organ systems. An extraordinary variety of different isoforms are produced from this gene through alternative promoter usage and splicing. These isoforms are expressed in a tissue-specific manner and differ significantly in their structure, and are classified as types I, II, III, IV, V and VI. Dysregulation of this gene has been linked to diseases such as cancer, schizophrenia, and bipolar disorder (BPD). [provided by RefSeq, Apr 2016]</p> <p>Transcript Variant: This variant (HRG-beta1d), which uses the type IV promoter, contains an alternate 5' terminal exon and thus differs in the 5' UTR and 5' coding region, and it also lacks two in-frame exons in the central coding region, compared to variant HRG-beta1. The resulting isoform (HRG-beta1d, also known as IV-beta 1a or type IV fetal A beta 1a) has a distinct N-terminus, lacks an internal segment, and is shorter than isoform HRG-beta1.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>