

Product datasheet for **SC327129**

AIFM3 (NM_001146288) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AIFM3 (NM_001146288) Human Untagged Clone
Tag:	Tag Free
Symbol:	AIFM3
Synonyms:	AIFL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327129 representing NM_001146288.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGGCGGCTGCTTCTCCAAACCCAAACCAGGCGCTGCCTGCCACAGTGGAGCTCAAGATCGAGGTG
GTGCTGCCTGAGAAGGAGCGAGGCAAGGAGGAGCTGTCGGCCAGTGGGAAGGGCAGCCCCGGGCTAC
CAGGGCAATGGCACGGCCCGCCACTTCCACACGGAGGAGCGCCTGTCCACCCCTACCCCTACCCAGC
CCTCAGGATTGCGTGGAGGCTGCTGTCTGCCACGTCAAGGACCTCGAGAATGGCCAGATGCGGGAAGTG
GAGCTGGGCTGGGGGAAGGTGTTGCTGGTGAAGGACAATGGGAGTTCACGCCCTGGGCCATAAGTGT
CCGCACTACGGCGCACCCCTGGTAAAGGCGTTCTGTCCCGTGGTGGGTGCGCTGCCCTGGCACGGC
GCCTGTTCAACATCAGCACTGGGGACCTGGAGGACTTCCCTGGCCTGGACAGTCTACACAAGTTCAG
GTGAAGATTGAGAAGGAGAAGGTGTACGTCGGGCCAGCAAGCAGGCCCTACAGTGCAGCGAAGGACC
AAGGTGATGGCCAAGTGTATCTCTCAAGTGTGGGTACAGCAGTAGCACCATGTGCTCATTGTGGGT
GCAGGTGCAGCTGGCCTGGTGTGTGCAGAGACTGCGGCAGGAGGGCTTCTCCGACCGGATCGTCTG
TGCACGCTAGACCGGCACCTTCCCTACGACCGTCCCAAGCTCAGCAAGTCCCTGGACACACAGCCTGAG
CAGCTGGCCCTGAGGCCAAGGAGTTTTTCCGAGCCTATGGCATCGAGGTGCTCACCGAGGCTCAGGTG
GTCACAGTGGACGTGAGAATAAGAAGGTGCTGTTCAAGGATGGCTTCAAGCTGGAGTACAGCAAGCTG
CTGCTGGCACCAGGGAGCAGCCCAAGACTCTGAGCTGCAAAGGCAAAGAAGTGGAGAACGTGTTCACT
ATCCGGACGCCAGAGGATGCCAATCGCGTGGTGGGCTGGCCCGAGGCCCAACGTGGTCTGTTGGGA
GCCGGCTTCTGGGGATGGAGGTGGCCGCTTACCTGACGGAGAAGGCCACTCTGTGTCTGTGGTGGAG
CTGGAGGAGACGCCCTTCCAGGAGTTCCTGGGGAGCGCGTGGTCTGCCCTCATGAAGATGTTTGAG
AACAAACCGGTGAAGTTCTACATGCAGACGGAGTGTCTGAGCTGCGGGCCAGGAGGGAAGCTGAAG
GAGGTTGTGCTGAAGAGCAGCAAGGTGCTGCGGGCTGACGTCTGCGTGGTGGGCATTGGTGCAGTGCCC
GCCACAGGCTTCTGAGGCAAAGCGGCATCGGTTTGGATTCCCGAGGCTTCCCTGTCAACAAGATG
ATGCAGACCAATGTCCCAGGCGTGTGTCAGTGGCGATGCTGTACCTTCCCCCTTGCCTGGAGGAAC
AACCGCAAAGTGAACATTCCACATTGGCAGATGGCTCATGCTCAGGGGCGCGTGGCAGCCAGAACATG
TTGGCGCAGGAGGCGGAGATGAGCACTGTGCCCTACCTCTGGACCGCCATGTTGGCAAGAGCCTGCGC
TACGCGGGCTACGGAGAAGGCTTCGACGACGTCATCATCCAGGGGATCTGGAGGAGCTGAAGTTTGTG
GCTTTTTACTAAAGGCGACGAGGTGATCGCCGTGGCCAGCATGAACTACGATCCATTGTGTCCAAG
GTCGCTGAGGTGCTGGCCTCAGGCCGTGCCATCCGGAAGCGGGAGGTGGAGACTGGCGACATGCTCTGG
CTTACGGGAAAGGATCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001146288

Insert Size: 1815 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146288.1](#)

RefSeq Size: 2439 bp

RefSeq ORF: 1815 bp

Locus ID: 150209

UniProt ID: [Q96NN9](#)

Cytogenetics: 22q11.21

Protein Families:	Druggable Genome
MW:	66.5 kDa
Gene Summary:	Induces apoptosis through a caspase dependent pathway. Reduces mitochondrial membrane potential.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) differs in the 5' UTR and has multiple coding region differences compared to variant 1. The encoded isoform (3) is shorter than isoform 1.