

Product datasheet for **SC327148**

ARHGAP25 (NM_001166276) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARHGAP25 (NM_001166276) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGAP25
Synonyms:	HEL-S-308; KAIA0053
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327148 representing NM_001166276.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCCCTCGGTCAAGGAGTGTGATGACTGGCGAG
CAGATGGCTGCCCTCCATCCATCGTCCACCCCAACCCGCTGGAGAGGCCATCAAGATGGGCTGGCTG
AAGAAGCAGAGGTCCATCGTGAAGAAGTGGCAGCAGAGGTAATTTGTGCTGAGGGCGCAGCAGCTCTAC
TACTACAAGGATGAAGAGGACACGAAGCCCAAGGCTGCATGTATCTACCAGGATGTACAATCAAGGAG
ATCGCCACAAACCCAGAAGAAGTGGGAAGTTTGTCTTTGAAATCATTCCAGCCTCATGGGACCAGAAT
CGCATGGGACAGGACTCCTATGCTCCTATGGCCAGCTCTCAGGCGGAGATGGAGGAGTGGTTAAATTC
CTCAGGAGAGTTGCTGGCACACCCTGTGGAGCAGTGTGGCCAGCGCTTGGATGAGACTGTGGCCTAT
GAACAGAAATTCGGCCCCATCTGGTGGCCATCCTGGTGGAGAAATGTGCAGAGTTCATCCTGGAGCAC
GGCCGGAATGAAGAGGGCATCTCCGTCTGCTGGCAGGACAACCTGGTGAAGCAGCTGAGAGACGCT
TTTGATGCTGGGAGCGCCCTCTTTGACAGAGACACAGATGTGCACACTGTGGCTTCCCTGTTAAAG
CTCTACCTCCGAGACCTCCAGAGCCCGTGGTTCCCTGGAGCCAGTACGAAGGGTTCCTGCTCTGTGGG
CAGCTCACGAATGCGGATGAGGCAAGGCTCAGCAGGAGTTGATGAAGCAGCTCTCCATCCTTCCCTCGT
GACAACTATAGTCTCCTGAGCTACATCTGCAGGTTCTACATGAAATACAGCTGAACTGTGCTGTTAAC
AAGATGAGTGTGGACAACCTGGCTACTGTGATTGGTGTGAATCTCATCAGGTGCAAGGTGCAAGACCCCT
GCCGTGATCATGAGAGGGACTCCTCAGATCCAAGAGTGTGACTATGATGATCAGAGACCATGAAGTC
CTTTCCCAAGTCCAAGGATAACCCCTGTACCCCTGCCAGAAAAATGACCCCAAGAAAGCTCCA
GTGGCCGAAGCTCTGTAGGCTGGGATGCCACTGAAGACCTCCGAATTTCTAGGACAGACAGCTTCAGT
AGCATGACAAGCGACTCTGATACAACCAGCCCAACCCGACAGCAGCCGAGCGATGCGTTTTCCGGAGGAC
AGCAGCAAAGTACCCAGGAAAAGCCAGGAGACTGGAAAATGCAATCTCGTAAAAGGACTCAAACACTC
CCTAACCGGAAATGTTTCTTGACATCAGCTTTTCAGGTTGCCAACAGCAGCAAAATGGAGATCTTTAAA
AATGAATTCTGGTCGCTTCTCAGAGGCTAAGGCAGGGGAAGGGCACAGGAGAACGATGTCTCAAGAC
TTGCGCAACTTTCTGACTCCCAACGGACTTCCACCTACGATAACGTCCCTTCCCTGCCAGGGTCCCTT
GGGAGGAAGCCAGTGCCTCTTCCCAAGCCTGTGACTCCAAGGAGATACTTTGCCAGTCCAAAC
TCTGAAACTGGGCTGGAAAAAAGAACTCTGGAGAAGAGGAAATGATTCTTTGCAGAGGATGGTCCAA
GAGCTACGAAAGGAAATAGAAACACAGAAGCAAAATGTATGAGGAACAGATTAACCTTGAGAAGGAA
AATTATGACGTTTGGGCTAAAGTGGTGAAGCTCAATGAAGAACTGGAGAAGGAAAAGAAGTCTGCA
GCCCTAGAGATCAGCCTCCGCAACATGGAGCGCTCCCGGAGGATGTTGAGAAGGGAACAAGGCCTTG
GAAGAAGAAGTCAAGGAATTTGTCAAATCCATGAAGGAACCCAAGACCGAGGCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001166276

Insert Size: 1920 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001166276.1](#)

RefSeq Size: 2960 bp

RefSeq ORF: 1920 bp

Locus ID: 9938

UniProt ID: [P42331](#)

Cytogenetics: 2p13.3

MW: 72.5 kDa

Gene Summary: ARHGAPs, such as ARHGAP25, encode negative regulators of Rho GTPases (see ARHA; MIM 165390), which are implicated in actin remodeling, cell polarity, and cell migration (Katoh and Katoh, 2004 [PubMed 15254788]).[supplied by OMIM, Mar 2008]
Transcript Variant: This variant (3) differs in the 5' UTR and 5' coding region, compared to variant 1. The resulting isoform (c) is shorter and has a distinct N-terminus, compared to isoform a.