

Product datasheet for **SC327161**

CHFR (NM_001161344) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CHFR (NM_001161344) Human Untagged Clone
Tag:	Tag Free
Symbol:	CHFR
Synonyms:	RNF116; RNF196
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327161 representing NM_001161344.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGCGGCCCGAGGAAGGCAAGCAGTCGCCGCCCGCAGCCCTGGGACGGCTCCTGCGTCTGGGC
GCGGAGGAGGGCGAGCCGACGTCTCTGAGGAAGCGGGAGTGGACCATCGGGCGAGACGAGGTTGC
GACCTTTCTTCCCAGCAATAAACTGGTCTCTGGAGATCACTGTAGAATTGTAGTGGATGAAAAATCA
GGTCAGGTGACTGGAAGATACCAGCACCAGTGAACAGTGATTAACAAGCTGAAGTTGTTAAGAAG
CAGACATGCCCTTACAGACTGGGGATGTCATCTACTTGGTGTACAGGAAGAATGAACCGGAACACAAC
GTGGCATACTCTATGAATCTTTAAGTAAAAGCAAGGCATGACACAAGAATCCTTTGAGCTAACAAAG
GAAAATGTGTTCCATGGGACAAAGATACCTCAGGTGCAGGTGCAGGGCGAGGGGCCGATCCCCGGGTC
CCTCCGTCGTGCGCCGCACTCAGGTGTGCTTTGAGGAACCACAGCCATCAACATCGACGTCAGACCTC
TCCCCACAGCCTCGGCCTTTCCACGGAGCCTTCTCCTGCAGGGCGAGAGCGTTCTCCAGTTGTGGG
TCTGGGGTGGTGGCATCTCCCCTAAAGGAAGTGGTCCCTCTGTGGCAAGTGATGAAGTCTCCAGCTTT
GCCTCAGCTCTCCAGACAGAAAGACTGCGTCTTTTCTGTCGTTGGAACCCAGGATCAGGAGGATTTG
GAGCCCGTGAAGAAGAAAATGAGAGGAGATGGGGACCTTGACCTGAACGGGCAGTTGTTGGTTCGCACAA
CCGCGTAGAAATGCCAAACCGTCCACGAGGACGTGAGCAGCGGCTGGGAAGCCAGACAAGATGGAG
GAGACGCTGACATGCATCATCTGCCAGGACCTGCTGCAGACTGCGTGAGTTTGCAGCCCTGCATGCAC
ACGTTCTGCGCGGCTTGCTACTCGGGCTGGATGGAGCGCTCGTCCCTGTGTCCTACCTGCCGCTGTCCC
GTGGAGCGGATCTGTA AAAAACCACATCTCAACAACCTCGTGAAGCATACTCATCCAGCATCCAGAC
CCCAAAGTCAGGCGGCTTTTTCTGATGAAGAAGGGAGTTCAGAGGACCTGCTGGAGCTGTCAGACGTT
GACAGTGAGTCCTCAGACATTAGCCAGCCATACGTCGTGTGCCGGCAGTGTCTGAGTACAGAAGGCAG
GCGGCGCAGCCTCCCCTGCCCAGCACCAGGGCGAGCCAGGAGCCCCACAGGCCCTGGGGATGCA
CCCTCCACGTCCGTGAGCCTGACGACAGCAGTCCAGGATTACGTGTGCCCTCTGCAAGGAAGCCACGCC
CTGTGCACCTGCTGCTTCCAGCCATGCCGACCCGAGAGCGGAGCGGAGCAGGACCCGCGTGTGCGCC
CCTCAGCAGTGTGCGGTGCTGCTGCAGCCTTTCTGCCACCTGTACTGGGGCTGCACCCGGACCGGCTGC
TACGGCTGCTGGCCCCGTTTTGTGAGCTCAACCTGGGTGACAAGTGTCTGGACGGCGTGTGAACAAC
AACAGCTACGAGTCAGACATCTGAAGAATTACCTGGCAACCAGAGGTTTACATGAAAAACATGTTG
ACCGAGAGCCTCGTGGCTCTCCAGCGGGAGTGTCTGCTGTCTGATTACAGAGTCACGGGAGACACC
GTTCTGTGTTACTGCTGTGGCCTGCGCAGCTTCGGTGAGCTGACCTATCAGTATCGGCAGAACATTCCT
GCTTCCGAGTTGCCAGTGCCGTAACATCCCCTGCTGACTGCTACTGGGGCCGTAACCTGCCGCACTCAG
GTGAAAGCTACCCAGCCATGAAATTCAATCATATCTGTGAACAGACAAGGTTCAAAAACATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001161344

Insert Size: 1995 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001161344.1](#)

RefSeq Size: 3291 bp

RefSeq ORF: 1995 bp

Locus ID: 55743

UniProt ID: [Q96EP1](#)

Cytogenetics: 12q24.33

Protein Families: Druggable Genome

MW: 73.4 kDa

Gene Summary: This gene encodes an E3 ubiquitin-protein ligase required for the maintenance of the antephasic checkpoint that regulates cell cycle entry into mitosis and, therefore, may play a key role in cell cycle progression and tumorigenesis. The encoded protein has an N-terminal forkhead-associated domain, a central RING-finger domain, and a cysteine-rich C-terminal region. Alternatively spliced transcript variants that encode different protein isoforms have been described. [provided by RefSeq, Mar 2014]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.