

Product datasheet for **SC327164**

RED1 (ADARB1) (NM_001160230) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RED1 (ADARB1) (NM_001160230) Human Untagged Clone
Tag:	Tag Free
Symbol:	ADARB1
Synonyms:	ADAR2; DRABA2; DRADA2; NEDHYMS; RED1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327164 representing NM_001160230.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGATATAGAAGATGAAGAAAACATGAGTTCACGACACTGATGTGAAGGAAAACCGCAATCTGGAC
AACGTGTCCCAAGGATGGCAGCACACCTGGCCTGGCAGGGCTCTCAGCTCTCCAATGGGGTGGT
GGTGGCCCGGCAGAAAGCGGCCCTGGAGGAGGCAGCAATGGCCACTCCAAGTACCGCTGAAGAAA
AGGAGGAAAACACCAGGGCCCGTCTCCCAAGAAGCCCTGATGCAGCTGAATGAGATCAAGCCTGGT
TTGCAGTACACTCCTGTCCCAGACTGGGCCCGTGCACGCGCCTTTGTTTGTATGTCTGTGGAGGTG
AATGGCCAGGTTTTGAGGGCTCTGGTCCACAAAAGAAAAGGCAAACTCCATGCTGCTGAGAAGGCC
TTGAGGTCTTTCAGTTTCTAATGCCTCTGAGGCCACCTGGCCATGGGGAGGACCCTGTCTGTC
AACACGGACTTCACATCTGACCAGGCCACTTCCCTGACACGCTCTTCAATGGTTTTGAACTCCTGAC
AAGGCGGAGCCTCCCTTTACGTGGGCTCCAATGGGGATGACTCCTTCAAGTTCCAGCGGGGACCTCAGC
TTGTCTGCTTCCCGGTGCTGCCAGCCTAGCCAGCCTCCTCTCCCTGTCTTACCACCTTCCACCC
CCGAGTGGGAAGAATCCCGTGTGATCTTGAACGAACTGCGCCAGGACTCAAGTATGACTTCTCTCC
GAGAGCGGGGAGAGCCATGCCAAGAGCTTCGTATGTCTGTGGTCTGGTGGTCAAGTCTTTGAAGGC
TCGGGGAGAAAACAAGAAGCTTGCCAAGGCCCGGCTGCGCAGTCTGCCCTGGCCGCCATTTTTAACTTG
CACTTGGATCAGACGCCATCTCGCCAGCCTATCCAGTGGGGTCTTACAGTGCATTTACCGCAGGTT
TTAGCTGACGCTGTCTCACGCTGGTCTGGGTAAGTTTGGTACCTGACCGACAACCTTCTCTCCCT
CACGCTCGCAGAAAAGTCTGGTGGAGTGCATGACAACAGGCACAGATGTTAAAGATGCCAAGGTG
ATAAGTGTCTACAGGAACAAAATGTATTAATGGTGAATACATGAGTGCATGCTGGCCTTGCAATTAAT
GACTGCCATGCAGAAATAATATCTCGGAGATCCTTGTCTCAGATTTCTTTATACACAACCTGAGCTTAC
TTAAATAACAAAGATGATCAAAAAGATCCATCTTTCAGAAATCAGAGCGAGGGGGTTTAGGCTGAAG
GAGAATGTCCAGTTTCATCTGTACATCAGCACCTCTCCCTGTGGAGATGCCAGAATCTTCTCACCACAT
GAGCCAATCCTGGAAGAACCAGCAGATAGACACCCAAATCGTAAAGCAAGAGGACAGCTACGGACAAA
ATAGAGTCTGGTGGGGACGATTCCAGTGCCTCCAATGCGAGCATCAAACGTGGGACGGGGTGTG
CAAGGGGAGCGGCTGCTCACCATGTCTGCAGTACAAGATTGCACGCTGGAACGTGGTGGGATCCAG
GGATCCCTGCTCAGCATTTTCTGGAGCCATTTACTTCTCGAGCATCATCTGGGCAGCCTTTACCAC
GGGACCACCTTTCCAGGGCCATGTACCAGCGGATCTCCAACATAGAGGACCTGCCACCTCTCTACACC
CTCAACAAGCCTTTGCTCAGTGGCATCAGCAATGCAGAAGCACGGCAGCCAGGGAAGGCCCAACTTC
AGTGTCAACTGGACGGTAGGCGACTCCGCTATTGAGGTGATCAACGCCACGACTGGGAAGGATGAGCTG
GGCCGCGCTCCCGCTGTGTAAGCACGCTGTACTGTGCTGGATGCGTGTGCACGGCAAGGTTCCC
TCCCACTTACTACGCTCCAAGATTACCAAGCCCAAGTGTACCATGAGTCCAAGCTGGCGGCAAAGGAG
TACCAGGCCCAAGGTACACTGA
AGCGGACCGACGCTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
```

Restriction Sites: SgfI-RsrII

Plasmid Map:


ACCN: NM_001160230

Insert Size: 2025 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001160230.1](#)

RefSeq Size: 3485 bp

RefSeq ORF: 2025 bp

Locus ID: 104

UniProt ID: [P78563](#)

Cytogenetics: 21q22.3

Protein Families: Druggable Genome

MW: 73.7 kDa

Gene Summary: This gene encodes the enzyme responsible for pre-mRNA editing of the glutamate receptor subunit B by site-specific deamination of adenosines. Studies in rat found that this enzyme acted on its own pre-mRNA molecules to convert an AA dinucleotide to an AI dinucleotide which resulted in a new splice site. Alternative splicing of this gene results in several transcript variants, some of which have been characterized by the presence or absence of an ALU cassette insert and a short or long C-terminal region. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (7), also known as ADAR2d, lacks the ALU cassette insert and an alternate segment in the 3' coding region, compared to variant 2. The resulting isoform (7) lacks an internal segment and has a shorter and distinct C-terminus, compared to isoform 2, and contains the short C-terminal region. Variants 7 and 10 encode the same isoform.