

Product datasheet for **SC327165**

ARHGAP40 (NM_001164431) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARHGAP40 (NM_001164431) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGAP40
Synonyms:	C20orf95; dJ1100H13.4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327165 representing NM_001164431.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCGAGCCTGCCCTCCTCCCGCCGCCAGATGGAGAGCTGGCCCCAGGGCCCTAGCCTCACCG
TGTCTCGGATCCCGCGGCCCGCATTGCCAGGCGCTGCGCCAGCGCTGGGCCGACCTGGGCTGCAGC
TCTGGCCCTCCTCAGGCCGAATGGATCAGCTTCCCAGAAGAATCTCCTGCGACTACACCCTGCTGGC
TCTGCTGGCTGTTCCACTGGAGTGGAGAGCTCCAGCATGGATGGCTTTTGGATGGAGGTGGAACAGATC
CAACAGAGAGATGAGCTGAGGGAAGAGGACAGTGGCGGAATGAAGGCCAGCTTCCAGAGGGGGAAGCT
GAATCCCAGTGGCTGCAGGACACAGGCTGTGGGCTCCTTGGTGGCTGGGCTGGATGGTGATCAC
CAGGAGCTCCTGTCCACCCTGACACAGACCCAGGTGGCCGCTGTGTGCCCGCGCTGGACATCTATGCT
CGCTCAGTGCAGAACACACAAGACACCTGTGAGAGATGTCAGGGATGCTTTGGGGTCTTCAATTCA
GGGAAAATGTCGTGAGAGAAATGGCGACTCCGGTATGAAGGGGGCCAGCTCAGTTCGGGGCCTCTAAG
TTTCTCCAGCAGCAGAGCTGGGGGGTGCAGGAGCAGGCTGGGAGGGAAGAAGCTTCAACATGGAC
TCTGCCTACTCAGAACAGCTGCCGTGCTCCTGCAGAGGAGCAGGCCATCCCGGGGAGGCACCTCTGCC
TGGGGCAAGTGTCCCTGCCGAAGTTTACCGTCCCAAAGGCAGACTTGGCGTGACGAGGATAGGAGAC
TTGTCCCTGCAGGATATGAGGAAGGTACCTTCTGTGCCCTCATAGAGCTGACCGCGCTCTGTGACATC
CTCGGCTTGGACCTGAAGAGGAGCAAGGCGGGGAGTGGAAAGCGGCAGAAACTCGCCTCTTTGGTGTG
CCCCTTGACAGCCTGCTAGAAGCTGACCACAAAGTCTCCCGCAGCACACAGGTCCCGCTGGTCTTCAA
GCCCTGCTGTCTGTTTGGAAAAGAGAGGACTGGACATGGAAGGCATTCTCAGGGTGCCCGGATCCCAG
GCCAGGGTCAAGGGGCTGGAACAGAAAATGGAGAGAGACTTCTATGCTGGCCTTTTGTAGCTGGGACGAG
GTTTCATACAATGACGCCTCTGATTTGCTCAAAGGTTTCATCCGGAAGCTGCCGACACCTTTGCTCAG
GCTGAGTACCTCCCGCCCTTCGCCGTGGTGCCTAACATCCCAACCTGAAGCAGCCCTGCAGGTTCTT
CACCTGCTCATCCTCATCTCCAGAACCAACAGAAATGCCTTAAAGGCACTGCTGGAATTCCTCAGG
AAGGTGGTGGCCCGGGAACAGCACAAAGATGACTCTGAGGAATGTTTCCACCGTATGGCCCTAAC
CTCTTTTGCACCAAGGGCGGCCCCCAAGCTCCCAAAGGCAAGGAGAAGCAGCTGGCAGAAGGGGCA
GCCGAGGTGGTGCAGATAATGGTGCCTACCAGGATCTGCTGTGGACGGTGCCTCTTCTGGTGCCT
CAGGTGCGAAAATGAACGACAGTAGCAGCAGGCGCCCCAGCTCTGCGACGAGGCTCAAGCTTGG
CTGCGGAGGATGCACGACAGGACAAGGCGGGGACGGCCTCGAGGCGACTCCAAGGTGGCAAAG
ATCCAGGTGGTCTGGCCTATCAAGGACCCCTTGAAGGTGCCTCTACCCCGCAGCACCAGTGGCCAC
GTCCTGAGGCAGTTCACAGAGCAGCTCAGCCCTGGTTCCAAGGTCAGGAAGACAGTGGGACATGGAC
AGCCTCCTTACACCATAGGTCATGGAGTACGCAACATCCTCCTCTATGAAGTTGGAGGGAATATC
AATGAGCATCGCCTGGACCCAGATGCCTACCTCTTAGATCTGTACCGTGCCAACCCGCATGGTGAGTGG
GTCCTTAAACAGAACCCACCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAAATCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001164431

Insert Size: 2025 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164431.1](#)

RefSeq Size: 2838 bp

RefSeq ORF: 2025 bp

Locus ID: 343578

Cytogenetics: 20q11.23

MW: 74.4 kDa