

## Product datasheet for **SC327183**

### ZNF611 (NM\_001161500) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF611 (NM_001161500) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF611
Mammalian Cell Selection:	None
Vector:	<u><a href="#">pCMV6-XL5</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001161500, the custom clone sequence may differ by one or more nucleotides

```

ATGTTACGTGAGGAAGCAGCTCAGAAGAGGAAAGGAAAGGAGCCAGGCATGGCTCTTCT
CAGGGACGCTTGACTTCCGGGATGTGGCTATAGAATTCTCATTGGCAGAGTGGAATGC
CTGAACCTTTCACAGAGGGCTTTGTACAGGGAAGTGATGTTGGAGAACTACAGGAACCTG
GAGGCTGTGGATATCTCTTCCAAATGCATGATGAAGGAGTCTTGTCAACAGGGCAAGGC
AATACAGAAGTGATCCACACAGGGACATTGCAAAGACATGAAAGTCATCACATTGGAGAT
TTTTGCTTCCAGGAAATTGAGAAAGAAATTCATGACATTGAGTTTCAGTGTCAAGGAAGAT
GAAAGAAATGGCCTTGAAGCACCCATGACAAAAATAAAAAAGTTGACTGGTAGCACAGAC
CAACATGATCACAGGCATGCTGGAACAAGCCTATTAAAGATCAGCTTGGATCAAGCTTT
TATTCACATCTGCCTGAACTCCACATATTTAGATCAAAGGTGAAATTGGTAATCAACTT
GAGAAGTCTACTAATGATGCTCCCTCAGTTTCAACATTCCAAAGAATTTCTGTAGGCC
CAAACCCAGATTTCTAATAACTATGGGAATAATCCCCTGAATTCTTACTCCACAA
AAACAGGAAGTACACATGAGAGAAAAATCTTCCAATGTAATAAGAGTGGCAAAGCCTTT
AATTGTAGCTCACTCTTAAGGAAACACCAGATACCCCTTTAGGAGACAAACAATATAAA
TGTGATGTATGTGGCAAGCTCTTAATCACGAGCAATACCTTGCCATGCCATGATAGATGT
CACACTGTTGAGAAACCTTACAAGTGTAAAGAGTGTGGCAAGACCTTCAGTCAGGAGTCA
TCCCTTACCTGCCATCGTAGACTTCATACTGGAGTAAAACGTTACAATTGTAATGAGTGT
GGCAAGATCTTTGGTCAAATTCAGCCCTTCTAATTGATAAGGCAATTGATACTGGAGAA
AATCCTTACAAGTGAATGAATGTGACAAGGCTTTTAAATCAACAATCACAACTTTACAT
CATAGAATTCATACTGGAGAGAAACCTTACAAATGTGAAGAATGTGACAAAGTTTTAGT
CGGAAATCAACCATTGAGACACATAAGAGAATTCATACTGGAGAGAAACCATACAGATGT
AAGGTTTGTGACACAGCTTTCACGTGGCATTACAGCTGGCTCGACATAGAAGAATTCAT
ACTGCAAAGAAAACCTTATAAATGTAATGAGTGTGGCAAGACCTTCAGTCACAAGTCATCC
CTTGATGCCATCATAGACTTCATGGTGGAGAGAAATCTTACAAATGTAAGGTTTGTGAC
AAGGCTTTTGTGTGGAGTTCACAACCTGGCAAAACATACTAGAATTGACTGTGGAGAAAA
CCTTACAAGTGAATGAATGTGGGAAGACCTTTGGTCAAATTCAGATCTTTAATTCAT
AAGTCAATTCATACTGGGAGCAACCTTACAAATGTGATGAATGTGAAAAGGTTTTAGT
CGTAAATCAAGCCTTGAACACATAAGATAGGTCATACTGGAGAGAAACCATACAATGT
AAGGTTTGTGACAAGGCTTTTGCCTGTCATTCTATCTGGCAAAACATACTAGAATTCAT
AGTGGAGAGAAACCTTACAAGTGAATGAGTGCAGCAAGACCTTCAGTCACAGGTCATAC
CTTGATGCCATCATAGAGTTCATAGTGGTGGAGAAACCTTACAAGTGAATGAGTGCAGC
AAGACCTTCAGTCGACAGTTCATCCCTTATTGCCATCGTAGACTTCATAGTGGTGGAGAA
CCTTACAAGTGAATGAGTGTGGCAATACCTTCCGTCAGTCTCATCCCTTATATACCAT
CGTAGACTTCATACTGGAGAGAAATCTTACAAATGTACAATTTGTGACAAGGCTTTTCGTG
CGTAATTCACCTCTGTCAAGACATACCAGAATTCACACTGCAGAGAAACCTTACAAGTGT
AATGAATGTGGGAAGGCTTTAATCAACAATCACACCTTTCACGTCATCATAGAATTCAT
ACTGGAGAGAAACCT

```

**Restriction Sites:** Please inquire

**ACCN:** NM\_001161500

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001161500.1, NP_001154972.1</u>
<b>RefSeq Size:</b>	4417 bp
<b>RefSeq ORF:</b>	2118 bp
<b>Locus ID:</b>	81856
<b>UniProt ID:</b>	<u>Q8N823</u>
<b>Cytogenetics:</b>	19q13.41
<b>Protein Families:</b>	Transcription Factors
<b>Gene Summary:</b>	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.