

Product datasheet for **SC327187**

DDHD2 (NM_001164232) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DDHD2 (NM_001164232) Human Untagged Clone
Tag:	Tag Free
Symbol:	DDHD2
Synonyms:	iPLA(1)gamma; iPLA1gamma; SAMWD1; SPG54
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327187 representing NM_001164232.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGTCATCAGTGCAGTCACAACAGGAGCAGTTGTCCAGTCAGATCCATCTCCGTCACAAACTCATGT
AGTTCCTTTGAGCTAATAGACATGGATGCTGGCAGCTTGTATGAACCAAGTTTCTCCCATTTGGTTTTAT
TGTAAGATAAATAGATTCTAAGGAGACATGGATTCTTTCAACTCTGAGGATTCACAGCAGCTGGAAGAG
GCATATAGCTCTGAAAAAGGTTGAATGGGAGAGTTGTTCTACTGATGGGGGCAGATATGATGTTTCAT
TTGGGGGAGAGGATGCGGTATGCTGTATACTGGGATGAAGTGGCATCGGAAGTGAGACGATGTACGTGG
TTTTACAAGGGGACAAAGACAATAAGTATGTTCCCTACTCGGAGAGCTTCAGCCAAGTTTTAGAGGAA
ACTTACATGCTTGCTGTAACCTTTGGATGAATGGAAAAAGAACTGGAATCTCCCAACAGAGAAATTATT
ATTTTACACAATCCAAGCTTATGGTGCATTACCAGCCAGTTGCAGGGTCTGATGATTGGGGTTCAACA
CCCACGGAGCAGGGTCGACCAAGAACTGTGAAGAGAGGAGTTGAGAATCTCTGTTGACATTCATTGT
GGAGAACCTTTACAAATAGATCACTTGGTTTTTGTAGTCCATGGGATTGGACCAGCTTGTGATCTCCGC
TTTCGAAGCATTGTACAGTGTGTTAATGATTTTCGCAAGTGTTCCTTGAAGTTGCTACAGACACATTTT
AAGAAAGCCCAAGAAAATCAGCAGATTGGGAGGGTAGAATTTCTCCAGTCAACTGGCACAGTCCTTTG
CATTCTACTGGTGTGGATGTAGATCTGCAGCGAATAACCTGCCAGCATTAAACGCCTCAGGCACCTTC
ACCAATGACACAATCTGGATGTCTTCTTCTACAATAGTCCACCTACTGTCAGACTATTGTGGACACA
GTTGCTTCTGAAATGAACCGAATATACACATTTTTCTACAGAGGAACCTGATTTCAAAGGGGGTGTA
TCCATTGCTGGTCATAGTTTAGGTTTCGCTTATATTGTTTGATATCCTAACAAATCAGAAAGATTCTTTG
GGGATATTGACAGTGAAAAGGATTCGCTAAATATTGTAATGGATCAAGGAGATACACCTACACTAGAG
GAAGATTTGAAGAACTTCAGCTCTCTGAATCTTTGATATCTTTGAGAAGGAGAAAGTAGATAAGGAA
GCTCTGGCTTTATGTACAGACCGAGATCTTCAGGAAATAGGAATTCCTTTAGGACCAAGAAAGAGATA
TTAAACTATTTACGACACAGAAAAAATCAATGGGATTAAGAGACCAGCCCCGAGCCTGCTTCAGGG
GCAACATCCCCAAGAATCTGAGTTCTGCAGTAGCAGTAATACTAGAAATGGTGACTATCTGGATGTT
GGCATTGGGCAGGTGTCTGTGAAATACCCCGGCTCATCTATAAACCAGAGATATTCTTTGCCTTTGGA
TCTCCCATTTGGAATGTTCTTACTGTCCGAGGACTAAAAAGAATTGATCCCAACTACAGATTTCCAACG
TGCAAAGGTTTCTCAATATTTATCACCTTTTGATCCTGTGGCCTATAGGATTGAACCAATGGTGGTC
CCAGGAGTGGAATTTGAGCCAATGCTGATCCACATCATAAAGGCAGGAAGCGGATGCATTTAGAAGTGA
AGAGAGGGCTTGACCAGGATGAGTATGGACCTTAAGAACAATGCTAGGTTTCGTCGGGATGGCCTGG
AAGTCTTTTACCAGAGCTCCATACCCTGCCTTACAAGCTTCAGAAACACCAGAAGAACTGAAGCAGAA
CCTGAATCAACTTCAGAGAAGCCTAGTGATGTTAACACAGAAGAGACCTCTGTGGCAGTTAAGAAGAA
GTCCTGCCTATCAATGTGGGGATGCTGAATGGAGGCCAACGCATTGACTATGTGCTACAGGAGAAGCCT
ATTGAAAGTTTAAATGAGTATTTATTTGCTTTACAAAGCCATCTATGCTACTGGGAGTCTGAAGATACA
GTATTGCTCGTCCTCAAAGAGATCTACCAAAACCAGGGTATCTTCCTTGATCAGCCTTTACAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI
ACCN: NM_001164232
Insert Size: 2136 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001164232.1</u>
RefSeq Size:	4546 bp
RefSeq ORF:	2136 bp
Locus ID:	23259
UniProt ID:	<u>O94830</u>
Cytogenetics:	8p11.23
MW:	81 kDa
Gene Summary:	This gene encodes a phospholipase enzyme containing sterile-alpha-motif (SAM), WWE, and DDHD domains. This protein participates in membrane trafficking between the endoplasmic reticulum and the Golgi body. Mutations in this gene can cause autosomal recessive spastic paraplegia 54. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013] Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same isoform (1).