

Product datasheet for **SC327192**

ZC3H14 (NM_001160104) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZC3H14 (NM_001160104) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZC3H14
Synonyms:	MRT56; MSUT-2; NY-REN-37; SUT2; UKp68
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327192 representing NM_001160104.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGATCGGCACCGAGATCAGCCGCAAGATCCGGAGTGCCATTAAGGGGAAATTACAAGAATTAGGA
GCTTATGTTGATGAAGAAGCTTCTGATTACATTATGGTGATGGTGCCCAACAAGAAAAGTCAGGACCAA
ATGACAGAGGATCTGTCCCTGTTTCTAGGGAACAACAATTTCGATTCACCGTATGGCTTCATGGTGTA
TTAGATAAACTTCGCTCTGTTACAACCTGAACCCTCTAGTCTGAAGTCTTCTGATACCAACATCTTTGAT
AGTAACGTGCCTTCAAACAAGAGCAATTCAGTCGGGGAGATGAGAGGAGGCATGAAGCTGCAGTGCCA
CCACTTGCCATTCTAGCGCGAGACCTGAAAAAGAGATTCCAGAGTTTCTACAAGTTCGCAGGAGTCA
AAAACCACAAATGTCAGACAGACTTACGATGATGGAGCTGCAACCCGACTAATGTCAACAGTGAAACCT
TTGAGGGAGCCAGCACCCCTCTGAAGATGTGATTGATATTAAGCCAGAACCAGATGATCTCATTGACGAA
GACCTCAACTTTGTGAGGAGAATCCCTTATCTCAGAAAAACCTACAGTGACACTTACATATGGTTCT
TCTCGCCCTTCTATTGAAATTTATCGACCCTGCAAGTAGAAAATGCAGATAGTGGTGTTCAATTTAAC
AGGTTGCAATTTCAACAGCAGCAGAATAGTATTTCATGCTGCCAAGCAGCTTGATATGCAGAGTAGTTGG
GTATATGAAACAGGACGTTTGTGTGAACAGAGGTGCTTAACAGCTTAGAAGAAACGTATAGTCCGTTT
TTAGAAACAACCTCGGAGAAAATGAGTATGGAGGATGAAAACCTTCGGAAGAGAAAGTTGCCTGTGGTA
AGTTTCAGTTGTTAAAGTAAAAAATTCATCATGATGGAGAAGAGGAGGAAGAAGATGATGATTACGGG
TCTCGAACAGGAAGCATCTCCAGCAGTGTGTCTGTGCCTGCAAAGCCTGAAAGGAGACCTTCTCTCCA
CCTTCTAAACAAGCTAACAGAATCTGATTTTGAAGGCTATATCTGAAGCTCAAGAAATCCGTAACAAAA
ACAACAACTACTCTACAGTCCACAGAAACAGACACTCCAGTTGCTCCAGAAGCTCGAACTTCTCAA
GAAGAATTGCTAGCAGAAGTGGTCCAGGGACAAGTAGGACCCCAAGAAATAGTCCCCCATTAAGAA
GAGGAAACAAGGAGATTCTGTAGAAAAAATCAAGGAATCAACAGAGGCAATTATTATCCCGACTG
CAAATCGACCCAGTAATGGCAGAACTCTGCAGATGAGTCAAGATTACTATGACATGGAATCCATGGTC
CATGCAGACACAAGATCATTTTCTGAAGAAGCCAAAGCTGTCTGAGGAAGTAGTAGTGGCACCAAC
CAAGAGTCGGGGATGAAGACTGCAGATTCCTTCGGGACTTTTCAGGACACCTTATGCAGACACGAGAT
CTTGTAACAACAGATAAACCTGCAAGTCCCAAGTTTATAGTGACGCTGGATGGTGTCCCCAGCCCCCA
GGATACATGTCAGATCAAGAGGAGGACATGTGCTTTGAAGGAATGAAACCCGTAACCAAACTGCAGCC
TCAAACAAGGGACTCAGAGGTCTCCTCCACCCACAGCAGTTGCACTTGCTGAGCAGGACGCTTGAGGAC
CCAAATGCTGAGATGAGTGAAGTGTGGCACAGAAACCAGAAAACTTTTGGAGCGCTGCAAGTAC
TGGCCTGCTGTAAAAATGGGGATGAGTGTGCCTACCATACCCCATCTCACCCCTGCAAAGCCTTCCCC
AATTGTAATTTGCTGAAAAATGTTTGTGTTTACCCAAATTGTAATATGATGCAAAGTGTACTAAA
CCAGATTGTCCCTTCACTCATGTGAGTAGAAGAATTCCAGTACTGTCTCCAAAACAGTTGCACCACCA
GCACCACCTTCCAGTAGTCAGCTCTGCCGTTACTTCCCTGCTTGAAGAAGATGGAATGTCCCTTCTAT
CATCCAAAACATTGTAGGTTTAACTCAATGTACAAGACCGGACTGCACATTCTACCATCCCACCATT
AATGTCCCACCACGACATGCCTTGAATGGATTTCGACCTCAAACCAGCGAATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001160104

Insert Size: 2193 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001160104.1](#)

RefSeq Size: 4190 bp

RefSeq ORF: 2193 bp

Locus ID: 79882

UniProt ID: [Q6PJT7](#)

Cytogenetics: 14q31.3

MW: 82.3 kDa

Gene Summary: The protein encoded by this gene is a poly(A)-binding protein that can affect gene expression and poly(A) tail length. The encoded protein may influence mRNA stability, nuclear export, and translation. [provided by RefSeq, May 2016]
Transcript Variant: This variant (6) uses two alternate in-frame splice sites in the 3' coding region, compared to variant 1. The resulting isoform (6) lacks a five-aa and a one-aa segment, compared to isoform 1.