

Product datasheet for **SC327203**

ADAM8 (NM_001164489) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ADAM8 (NM_001164489) Human Untagged Clone
Tag:	Tag Free
Symbol:	ADAM8
Synonyms:	CD156; CD156a; MS2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327203 representing NM_001164489.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCGCGGCCCTCGGGCTCTGGCTGCTGGCGCGATGATGCTGCCTGCGATTGCCCCAGCCGGCCCTGG
GCCCTCATGGAGCAGATAGAGTCTGTTGCCGTGGCGTCTGCCAGGCCCCAGTCCGCCGAGCTCTG
CCCTCCCACTTGGGCCTGCACCCAGAGAGGGTGAAGTACGTCCTTGGGGCCACAGGGCACAACTTACC
CTCCACCTGCGGAAGAACAGGGACCTGCTGGGCTCCGGCTACACAGAGACCTATACGGCTGCCAATGGC
TCCGAGGTGACGGAGCAGCCTCGCGGGCAGGACTGCTTCTACCAGGGCCACGTAGAGGGGTACCCG
GACTCAGCCGCCAGCCTCAGCACCTGTCCGGCCTCAGGGGTTTCTTCCAGGTGGGGTACAGCTGCAC
CTGATCGAGCCCTGGATGAAGTGGCGAGGGCGGACGGCAGCCGTGTACCAGGCTGAGCACCTGCTG
CAGACGGCCGGGACCTGCGGGTTCAGCGACGACGCTGGGCAGCCTCCTGGACCCGGACGGCAGCC
GTCTTCAGGCCTCGGCCGGGACTCTCTGCCATCCCAGAGACCCGCTACGTGGAGCTGTATGTGGTC
GTGGACAATGCAGAGTTCAGATGCTGGGAGCGAAGCAGCCGTGCGTCAATCGGGTCTGGAGGTGGTG
AATCACGTGGACAAGCTATATCAGAAACTCAACTTCCGTGTGGTCTGGTGGGCCTGGAGATTTGGAAT
AGTCAGGACAGGTTCCACGTGAGCCCGACCCAGTGTACACTGGAGAACCTCCTGACCTGGCAGGCA
CGGCAACGGACACGGCGGCACCTGCATGACAACGTACAGCTCATCACGGGTGTCGACTTACCCGGGACT
ACCGTGGGGTTTCCAGGGTGTCCGCCATGTGCTCCCACAGCTCAGGGGCTGTGAACCAGGACCACAGC
AAGAACCCCGTGGGCGTGGCCTGTACCATGGCCATGAGATGGGCCACAACCTGGGCATGGACCATGAT
GAGAACGTCCAGGGTCCCGTCCAGGAACGCTTCGAGGCCGGCCGCTGCATCATGGCGGGCAGCATT
GGCTCCAGTTTCCCAGGATGTTCAAGTACTGCAGCCAGCCCTACCTGGAGAGCTTTTGGAGCGGCCG
CAGTCGGTGTGCCTCGCCAACGCCCTGACCTCAGCCACCTGGTGGCGGCCCGCTGTGGGAACCTG
TTTGTGGAGCGTGGGGAGCAGTGCAGTCCGGCCCGCCCGAGGACTGCCGGAACCGCTGCTGCAACTCT
ACCACCTGCCAGCTGGCTGAGGGGGCCAGTGTGCGCACGGTACCTGCTGCCAGGAGTGAAGTGAAG
CCGGCTGGTGAAGTGTGCGTCCCAAGAAGGACATGTGTGACCTCGAGGAGTTCTGTGACGGCCGGCAC
CCTGAGTGCCCGAAGACGCCTTCCAGGAGAACGGCACGCCCTGCTCCGGGGGCTACTGTACAACGGG
GCCTGTCCACACTGGCCAGCAGTCCAGGCCTTCTGGGGCCAGGTGGGCAGGCTGCCGAGGAGTCC
TGCTTCTCCTATGACATCTACCAGGCTGCAAGGCCAGCCGTTACAGGGCTGACATGTGTGGCCTCTG
CAGTGCAAGGGTGGGCAGCAGCCCTGGGGCGTGCCATCTGCATCGTGGATGTGTGCCACGCGCTCACC
ACAGAGGATGGCACTGCGTATGAACCAGTCCCGAGGGCACCCGTTGTGGACCAGAGAAGTTTGTCTGG
AAAGGACGTTGCCAGGACTTACAGTTCACAGATCCAGCAACTGCTCTGCCAGTCCACAACCATGGG
GTGTGCAACCACAAGCAGGAGTGCCACTGCCACGCGGGTGGGCCCGCCCACTGCGCGAAGCTGCTG
ACTGAGGTGCACGAGGCTGCCAGCCGCTGCCGGCAAGGGCGGGGCTCCAGCCCCATCCAGGGGCC
CCAAGAGCTGGTCCCACCACCCACCCGGGCCAGCCCGCCGACACCCGGCTCCTCGGTGGCTCTGAA
GAGGCCGCCCCCTGCTCCTCCGGTCACTGTGTCCAGCCACCCCTCCCAGTTCCTGTCTACACCCGGCA
GGCACCAAGCAGGTCAAGCCAACGTTCCGACCCCAAGTCCCGCCAGTCAAACCCGGGGTGGTGC
GGCCAACCTGGTCCAGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001164489

Insert Size: 2229 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164489.1](#)

RefSeq Size: 3149 bp

RefSeq ORF: 2229 bp

Locus ID: 101

UniProt ID: [P78325](#)

Cytogenetics: 10q26.3

Protein Families: Druggable Genome, Transmembrane

MW: 80.3 kDa

Gene Summary: This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The protein encoded by this gene may be involved in cell adhesion during neurodegeneration, and it is thought to be a target for allergic respiratory diseases, including asthma. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2009]

Transcript Variant: This variant (2) lacks an alternate exon and uses an alternate splice site, which causes a frameshift in the 3' coding region, compared to variant 1. The resulting isoform (2) has a distinct C-terminus and is shorter than isoform 1.