

Product datasheet for **SC327208**

DISC1 (NM_001164539) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DISC1 (NM_001164539) Human Untagged Clone
Tag:	Tag Free
Symbol:	DISC1
Synonyms:	C1orf136; SCZD9
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001164539, the custom clone sequence may differ by one or more nucleotides

```

ATGCCAGGCGGGGTCCTCAGGGCGCCCCAGCCGCCCGCGGGCGGGCGGTGAGCCAC
CGCGCAGGCAGCCGGGATTGCTTACCACCTGCAGCGTGCCTTCGGAGGCGGGCTGGCA
CGGAGGCCGGGCTACATGAGAAGCTCGACAGGGCTGGATCGGGTTCCTTTCCCGCA
GTGGGCACACTGTTCCGGTTCAGGAGGGGTGTCTGGCGAGGAGTCCCACCACTCGGAG
TCCAGGGCCAGACAGTGTGGCCTTGACTCGAGAGGCCTTTGGTCCGGAGCCCTGTTCC
AAGAGTGCAGCAGCCCTACTGTGACCTGTGTGAGAGGAACCTCGGCGCACTTTGGGATT
CAGCTCAGAGGTGGCACCAGATTGCCTGACAGGCTTAGCTGGCCGTGTGGCCCTGGGAGT
GCTGGGTGGCAGCAAGAGTTTGCAGCCATGGATAGTTCTGAGACCCTGGACGCCAGCTGG
GAGGCAGCTGCAGCATGGAGCAAGGCGTGTCCGGGCAGCAGGCTCTCTGCCATCAGCA
GAGTTGAGTAGCAACAGCTGCAGCCCTGGCTGTGGCCCTGAGTCCCCCAACCCCTCT
GGCTCTCACAGTGCCTTACCTCAAGCTTAGCTTTATTCGGCTCTCGCTGGCTCTGCC
GGGAACGTGGAGAAGCAGAAGGCTGCCACCATCCAGAGAGGCTGAGTCCCATTGCCAG
AGCCCCAGGAGATGGGAGCCAAAGCTGCCAGCTTGACGGGCTCACGAGGACCCGCGA
TGTCTCTCTCGGCCCTTCACTCTTTGGCTACACGGGTCTCTGCAGACTTGGCCAGGCC
GCAAGGAACAGCTCCAGGCCAGAGCGTGACATGCATTCTTACCAGACATGGACCCTGGC
TCTCCAGTTCTCTGGATCCCTCACTGGCTGGCTGTGGTGGTGTGGGAGCAGCGGCTCA
GGGGATGCCCACTCTTGGGACACCCTGCTCAGGAAATGGGAGCCAGTGTGCGGGACTGC
CTGCTGAGAAACCGGAGGCAGATGGAGGTAATATCCTTAAGATTAACCACTCAGAACTT
CAGGAAGATGCAGTTGAGAATGATGATTATGATAAAGCTGAGACGTTACAACAAAGATTA
GAAGACCTGGAACAAGAGAAAATCAGCCTGCACCTTCACTTCTTCAAGGCAGCCAGCT
CTTAGCAGTTTCTGGGTACCTGGCAGCACAAGTCCAGGCTGCCTTGGCCGTGGGGCC
ACTCAGCAGGCCAGCGGAGATGACACCCACACCCCACTGAGAATGGAGCCGAGGCTGTTG
GAACCCACTGCTCAGGACAGCTTGCACGTGCCATCACGAGACGAGACTGGCTTCTTCAG
GAAAAGCAGCAGCTACAGAAAGAAATCGAAGCTCTCAAGCAAGGATGTTTGTGCTGGAA
GCCAAAGATCAACAGCTGAGAAGGGAAATAGAGGAGCAAGAGCAGCAACTCCAGTGGCAG
GGCTGCGACCTGACCCCACTGGTGGGCCAGCTGTCCCTGGGTGAGCTGCAGGAGGTCAGC
AAGGCCTTGCAGGACACCCTGGCCTCAGCCGGTCAGATTCCCTTCCATGCAGAGCCACCG
GAAACCATAAGGAGCCTCCAGGAAAGAATAAAATCCCTCACTTGTCACTTAAAGAAATC
ACTACTAAGGTGTATGAGTGAGAAATCTGCAGCACCTGAGGAAGAAAGTTAACGAT
ATTGAAACCCAACTACCAGCCTTGTGTAAGCCAAAATGCATGCCATATCAGGAAACCAT
TTCTGGACGGCTAAAGACCTCACCGAGGAGATTAGATCATTAAACATCAGAGAGAGAAGGG
CTGGAGGGACTCCTCAGCAAGCTGTTGGTGTGAGTTCAGGAATGTCAAAAAGCTGGGA
AGTGTTAAGAAGATTACAACAGACTGAGAAGAGAAGTGGAGCACCAGGAGACTGCCTAT
GATGGAGTCTCACTCTGTGCCCCAGTCTGGAGTGCAGTGGTACGATCTTGCTCACTGCAA
CCTCTGCCTCCCGAGTTCAAGCAATTCTCCTGCCTCAGCCTCCGGAGTAGCTGGGATTAC
AGGTGCCCCGCCACCATGCCTGGCTAATTTGTATTCTTAGTAGAGATGGGGTTTTACCAT
GTTGACCAGACTGGTCTAAAACCTCCTGACCTCAAGTATCCACCTTCTCAGCCTCCCAA
AGTGCTGGGATTACAGACATGAGCCACTGTGCCCTGGCCCTACAA

```

Restriction Sites: Please inquire

ACCN: NM_001164539

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001164539.1</u> , <u>NP_001158011.1</u>
RefSeq Size:	2676 bp
RefSeq ORF:	2268 bp
Locus ID:	27185
UniProt ID:	<u>Q9NRI5</u>
Cytogenetics:	1q42.2
Gene Summary:	<p>This gene encodes a protein with multiple coiled coil motifs which is located in the nucleus, cytoplasm and mitochondria. The protein is involved in neurite outgrowth and cortical development through its interaction with other proteins. This gene is disrupted in a t(1;11) (q42.1;q14.3) translocation which segregates with schizophrenia and related psychiatric disorders in a large Scottish family. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (c) lacks several 3' exons but has an alternate 3' exon, as compared to variant L. The resulting isoform (c, also known as isoform 26) is shorter and has a distinct C-terminus, as compared to isoform L.</p>