

Product datasheet for **SC327225**

TXNDC16 (NM_001160047) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TXNDC16 (NM_001160047) Human Untagged Clone
Tag:	Tag Free
Symbol:	TXNDC16
Synonyms:	ERp90; KIAA1344
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001160047, the custom clone sequence may differ by one or more nucleotides

```
ATGTTTTCCGGCTTCAATGTCTTTAGAGTTGGGATCTCTTTTGTGATAATGTGCATTTTT
TACATGCCAACAGTAACTCTTTACCAGAACTGAGTCCTCAGAAATATTTTAGTACATTG
CAACCAGGAAAAGCCTCTTTAGCTTATTTTGTCAAGCTGATTCCCAAGAACATCTGTA
TTTCTTGAAGAACTGAATGAGGCTGTTAGACCTCTGCAGGACTATGGAATTTGAGTTGCC
AAGGTTAATTGTGTCAAAGAAGAAATATCAAGATACTGTGGAAAAGAAAAGGATTTGATG
AAAGCATATTTATTCAAGGGCAACATATTGCTCAGAGAATTCCTACTGACACCTTGTTT
GATGTGAATGCCATTGTGCGCCATGTTCTCTTTGAAGTGAATATATTACCAACCTGGAA
GACCTTCAGAACATAGAAAATGCTCTGAAAGGAAAAGCAAATATTATATTCTCATATGTA
AGAGCCATTGGAATACCAGAGCACAGAGCAGTCATGGAAGCCGCTTTTGTGTATGGGACT
ACATACCAATTTGTCTTAACCACAGAAATTGCCCTTTTGGAAAGTATTGGCTCTGAGGAT
GTGGAATATGCACATCTCTACTTTTTTTCATTGTAACTAGTCTTGGACTTGACCCAGCAA
TGTAAGAAGAACTAATGGAACAGCCATTGACTACACTGAACATTCACCTGTTTATTAAG
ACAATGAAAGCACCTCTGTTGACTGAAGTTGCTGAAGATCCTCAACAAGTTTCAACTGTC
CATCTCCAACCTGGGCTTACCACTGGTTTTTATTGTTAGCCAACAGGCTACTTATGAAGCT
GATAGAAGAACTGCAGAATGGGTTGCTTGGCGTCTTCTGGGAAAAGCAGGAGTTCTACTC
TTGTTAAGGGACTCTTTGGAAGTGAACATTCTCAAGATGCTAATGTGGTCTTCAAAGA
GCAGAAGAGGGAGTTCCAGTGAATTTTTGGTATTACATGATGTTGATTTAATAATATCT
CATGTGGAAAATAATATGCACATTGAGGAAATACAAGAAGATGAAGACAATGACATGGAA
GGTCCAGATATAGATGTTTCCAGTATGAAAGTGGCAGAACTGTTTTTCCAGAGATAGGAAG
AGAAAATTACCTTTGGAACCTACAGTGGAACTAACAGAAGAAACATTTAATGCAACAGTG
ATGGCTTCTGACAGCATAGTACTCTTCTATGCTGGTTGGCAAGCAGTATCCATGGCATT
TTGCAATCCTATATTGATGTGGCAGTTAACTGAAAGGCACATCTACTATGCTTCTACT
AGAATAAACTGTGCAGATTGGTCTGATGTATGTAAGCAAAATGTTACTGAATTTCTCT
ATCATAAAGATGTACAAGAAAGGCGAGAAGCCAGTATCTTATGCTGGAATGTTAGGAACC
GAAGATCTCCTAAAATTTATCCAGCTCAACAGGATTTTATATCCAGTGAATATAACATCG
ATCCAAGAAGCAGAAGAATATTTAAGTGGGGAATTATATAAGACCTCATCTTGATTCT
AGTGTGTCAGTATTGGGACTATTTAGTCCAACCATGAAAACAGCAAAAGAAGATTTTAGT
GAAGCAGGAACTACCTAAAAGGATATGTTATCACTGGAATTTATTCTGAAGAAGATGTT
TTGCTACTGTCAACCAATATGCTGCAAGTCTTCCAGCCCTGCTGCTTGCAGACACACA
GAAGGCAAAATAGAGAGCATCCCACTAGCTAGCACACATGCACAAGACATAGTTCAAATA
ATAACAGATGCACTACTGGAATGTTTCCGGAAATCACTGTGGAATCTTCCCAGTTAT
TTCAGACTTCAGAAACCATTTATTGATTTTGTTCAGTATGGCACTGTAAATCCTCAGTAT
AAAAAAGCAATATTGACACTGGTAAAGCAGAAATCTTGGATTCAATTTACTCCATGCTGG
TTAAATCTAAAGAATACTCCAGTGGGGAGAGGAATCTTGAGGGCATATTTTGATCCTCTG
CCTCCCCTTCTCTTCTTGTGTTTGGTGAATCTGCATTGAGGTGGCCAAGTATTTGCATTT
CCTTCAGACCAGGCTATAATTGAAGAAAACCTTGTATTGTGGCTGAAGAAATTAGAAGCA
GGACTAGAAAATCATATCACAATTTTACCTGCTCAAGAATGGAACCTCCTCTTCCAGCT
TATGATTTTCTAAGTATGATAGATGCCGCAACATCTCAAGCTGGCACTAGGAAAAGTTCCC
AAGTGTATGAAAGAAACAGATGTGCAGGAGAATGATAAGGAACAACATGAAGATAAATCG
GCAGTCAGAAAAGAACCATTGAACTCTGAGAAATAAGCATTGGAATAGAAGTAATTGG
TTTAAAGAAGCAGAAAAATCATTTAGACGTGATAAAGAGTTAGGATGCTCAAAAGTGAAC
```

Restriction Sites: Please inquire

ACCN: NM_001160047

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001160047.1</u> , <u>NP_001153519.1</u>
RefSeq Size:	4636 bp
RefSeq ORF:	2463 bp
Locus ID:	57544
UniProt ID:	<u>Q9P2K2</u>
Cytogenetics:	14q22.1
Protein Families:	Druggable Genome, Transmembrane