

Product datasheet for **SC327246**

HSPC302 (TBCK) (NM_001163436) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HSPC302 (TBCK) (NM_001163436) Human Untagged Clone
Tag:	Tag Free
Symbol:	TBCK
Synonyms:	HSPC302; IHPRF3; TBCKL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327246 representing NM_001163436.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGTTTCCCCTGAAGGACGCTGAAATGGGAGCCTTACCTTCTTTCGCTCGGCTCTGCCACATGATGTT
TGTGGAAGCAATGGACTTCCCTCTCACACCAAATCCATCAAATTTTAGGGCGCTTCAAATCCTTAAA
ACCATCACCCATCCAGACTCTGCCAGTATGTGGATATTTCTAGGGGAAAGCATGAACGACTAGTGGTC
GTGGCTGAACATTGTGAACGTAGTCTGGAAGACTTGCTTCGAGAAAGGAAACCTGTGAGCTGTTCAACG
GTTTTGTGTATAGCATTGAGGTTCTTCAGGGCTTGCAGTATATGAACAAACATGGTATAGTACACAGG
GCATTGTCTCCTCATAATCCTGTTGGACCGAAAGGACATATAAATTGGCTAAATTTGGACTTTAT
CACATGACAGCTCATGGTATGATGTTGATTTCCAATAGGGTATCCCTCGTACTTGGCCCTGAGGTA
ATTGCACAGGGAATTTCAAACCAGTATCACATGCCAAGTAAAAACCATTGCCTTCTGGCCCCAAA
TCAGATGATGGTCTCTTGAATCATTATTTAGCTTTGTGTGGGAAGAAAATATTTTCAGAGCTTG
GATATTTCTGAAAGACTAAAATTTTGGTACTTTGGATTGTGTAGATGACACTTAAATAGTTCTGGCT
GAAGAGCATGGTTGTTTGGACATTATAAAGGAGCTTCTGAAACTGTGATAGATCTTTTGAATAAGTGC
CTTACCTTCCATCCTTCTAAGAGGCCAACCCAGATGAATTAATGAAGGCAAAAGTATTCAGTGAGGTA
TCACCTTTATATACCCCTTTACCAAACCTGCCAGTCTGTTTTTCACTTCTCTGAGATGTGCTGATTTA
ACTCTGCCTGAGGATATCAGTCAGTTGTGTAAGATATAAATAATGATTACCTGGCAGAAAGATCTATT
GAAGAAGTGATTACCTTTGGTGTGGCTGGAGGTGACTTGGAGAAAGAGCTTGTCAACAAGGAAATC
ATTCGATCCAAACCACCTATCTGCACACTCCCCAATTTCTCTTTGAGGATGGTAAAGCTTTGGACAA
GGTCGAGATAGAAGCTCGCTTTTAGATGATACCACTGTGACATTGCGTTATGCCAGCTAAGAAATAGA
TTGAAAGATGTTGGTGGAGAAGCATTTTACCCATTACTTGAAGATGACCACTAAATTTACCTCATTCA
AACAGCAATAATGAGTTGTCTGCAGCTGCCACGCTCCCTTTAATCATCAGAGAGAAGGATACAGAGTAC
CAACTAAATAGAATTATCTCTTCGACAGGCTGCTAAAGGCTTATCCATATAAAAAAACCAAACTGG
AAAGAAGCAAGAGTTGACATTCTCCTCTTATGAGAGGTTAACCTGGGCTGCTCTTCTGGGAGTTGAG
GGAGCTATTATGCCAAGTACGATGCAATTGATAAAGACACTCCAATTCACAGATAGACAAATTGAA
GTGGATATTCCTCGTGTATCAGTACGATGAAGTGTATCATCACCAGAAGGTCATGCAAAATTTAGG
CGTGTATTAAGCCTGGGTAGTGTCTCATCCTGATCTTGTGATTGGCAAGGCTTGTACTACTTTGT
GCTCCATTCTATATCTAACTTCAATAATGAAGCCTTGGCTTATGCATGTATGTCTGCTTTTATCCC
AAATACCTGTATAACTTCTTCTTAAAAGACAACCTCACATGTAATACAAGAGTATCTGACTGTCTTCT
CAGATGATTGCATTTATGATCCAGAGCTGAGTAAATCATCTCAATGAGATTGGTTTCATTCCAGATCTC
TATGCCATCCCTTGGTTTCTTACCATGTTTACTCATGATTTCCACTACACAAAATTTCCACCTCTGG
GATACCTTACTACTTGGGAATTCCTCTTTCCATTCTGTATTGGAGTAGCAATTCCTCAGCAGCTGCGG
GACCGGCTTTTGGCTAATGGCTTAAATGAGTGTATTCTTCTTCTCCGATTTACCAGAAATTGACATT
GAACGCTGTGTGAGAGAATCTATCAACCTGTTTGTGGACTCCTAAAAGTGCTACTTACAGACAGCAT
GCTCAACCTCAAAGCCATCTTCTGACAGCAGTGGAGGCAGAAGTTCGGCACCTTATTTCTCTGTGAG
TGTCAGATCCTCAAAGACAGATCTGTCAAGAGAATCCATCCCATTAATGACCTGAAGTCAGAAGTA
TCACCACGGATTTACGACAGAGGACCTGATTGACTTGTGTGAGCTCACAGTGACAGGCCACTTCAAACA
CCCAGCAAGAAAACAAAGTCCAGTAAACCAAGCTCCTGGTGGTTGACATCCGGAATAGTGAAGACTTT
ATTCTGGTACATTTACGGAAGCATCAACATTCCATTAGTGTGCCTTCACTGAGAAGGGGAGCTT
ACCCAGGGCCCTTACTGCTATGCTCCAGAACTCAAAGGGAAGGTCATTGTGATCGTGGGGCATGTG
GCAAAACACACAGCTGAGTTTGCAGCTCACCTTGTGAAGATGAAATATCCAAGAATCTGTATTCTAGAT
GGTGGCATTAAATAAAATAAAGCCAACAGGCCTCCTACCATCCCATCTCCTCAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001163436

Insert Size: 2682 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001163436.2](#)

RefSeq Size: 5130 bp

RefSeq ORF: 2682 bp

Locus ID: 93627

UniProt ID: [Q8TEA7](#)

Cytogenetics: 4q24

Protein Families: Druggable Genome, Protein Kinase

MW: 100.7 kDa

Gene Summary: This gene encodes a protein that contains a protein kinase domain, a Rhodanase-like domain and the Tre-2/Bub2/Cdc16 (TBC) domain. The encoded protein is thought to play a role in actin organization, cell growth and cell proliferation by regulating the mammalian target of the rapamycin (mTOR) signaling pathway. This protein may also be involved in the transcriptional regulation of the components of the mTOR complex. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014]
Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.