

Product datasheet for **SC327263**

HERC6 (NM_001165136) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HERC6 (NM_001165136) Human Untagged Clone
Tag:	Tag Free
Symbol:	HERC6
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001165136, the custom clone sequence may differ by one or more nucleotides

```

ATGTACTTCTGTTGGGGCGCCGACTCCAGGGAGCTGCAGCGCCGGAGGACGGCGGGCAGC
CCCGGGCTGAGCTACTGCAGGCGGCCAGCGGGGAGCGCCACTCTCTGCTGCTGAC
AACCACAGGGTCTCTCGTGGGAGACAACAGCAGGGGTACAGTGGGCCGAGGGGCGG
CAGCGCGGGGAGCTGCCAGAACCAATTCAGGCATTGAAAACCTAATTGTTGATCTGTG
AGCTGCGGGAAGGAGCACTCCCTGGCTGTGTGCCACAAAGGAAGGTCTTCGCATGGGA
GCTGGTTCTGAAGGCGAGCTGGGATTGGAGAATTCAGGAAATAAGTTTCACACCTAAG
AAAATAATGACTCTGAATGATATAAAAAAATAACAAGTTTCTGTGGACTACCACTCC
CTGGCATTATCAAAAGATAGCCAAGTGTTCCTGGGAAAGAACAGCCATGGGCAGCTG
GGCTTGGGGAAGGAGTCCCTCCCAAGCCAGCCCGCAGAGGGTGGAGTCCCTGGAGGG
ATCCCACTGGCTCAGGTGGCTGCCGGAGGGGCTCACAGCTTTCCTGTCTCTGTGGG
ACTTCGTTTGGTGGGAAAGTAACAGTGCCGGGAGCTGGCCCTCAGTGGGCGTAATGTC
CCAGTGCAAAGCAACAAGCCTCTCTCAGTCGGTGCAGTGAAGAATCTAGGTGTGGTTTAT
ATCAGCTGTGGTGTGCACACACTGCGGTGCTTACCCAGGACGGGAAAGTTCACATTT
GGAGACAATCGCTCTGGACAGCTGGGATACAGCCCACTCTGAGAAGAGAGTCCACAA
CTTGTGGAAAGAATTGATGGCCTAGTTTCGCAGATAGATTGTGGAAGTTATCACACCTG
GCATATGTGCACACCACTGGTCAAGTGGTATCTTTTGGTCATGGACCAAGTGACACAAGC
AAGCCAATCATCCGGAGGCCCTGACAGAGAACTTTGACATTAGCTGCCTGATTTCTGCT
GAAGACTTCGTTGATGTTCAAGTCAAACACATTTTGTGGAACATATGCCAACTTTGTG
ACAACATCATCAGGATACTAGTCCACACGTGCTCCCGGAAAACCTGCCAGAAAAGC
CGAATTAGCCAGTCCATGGCAGAAAAATGGATAGCAGTAAAAGAAGAAGTACTGAACAT
GAAATGGCTAAAAGTGAATTAAGATGATTTTTTCATCTCTGCTGTGCTGACTGCAAGT
TTTTTAAAGAAAAGGGAAGTGGAGAAAAGCACTTCCATTGATGTGGACTTAGAAAAGGCA
AGAGATACCTTCAAGAAGTTAACAAAAAGGAATGGATTTCTTCCATGATAACTACGTGT
CTCGAGGATGATCTGCTCAGAGCTTCCATGCCATTCTCCACCAAGAAGCTTTATCA
GTTTTCTCTGCTCCAGAAATGCTCTGTGATGCATGATTCTAAGAACTGGAAGAAGCTG
GTGTTCCATTTGCAAAGGCTGTGTGAAATGAGTAAACAATCTTTGCAAGTCTAAAG
AAGTGTGGGCATTTTGAAGAATCTTCTCTGAATCCGCTGATCCAGATGCTTAAAGCA
GCCATCATCTCAGCTGCTTCATCAGACTAAAACCGAACAGGATCACTGTAATGTTAAA
GCTCTTTTAGGAATGATGAAAGAACTGCATAAGGTAACAAGCTAACTGTCGACTACCA
GAAAATACTTTCAACATAAATGAACTCTCAACTTATTAAGCTTTTATATAGATAGAGGA
AGACAGCTCTTTCGGGATAACCACCTGATGTCAGAAAAGAAAGCATACTGCTTATGCAT
GAAACAATTCTGCAAAAAAGGATGAATTTCTCCATCACCCAGATTTATACTTAGAGTC
AGACGAAGTGCCTGGTTAAAGATGCTCTGCGTCAATTAAGTCAAGCTGAAGCTACTGAC
TTCTGCAAAGTATTAGTGGTTGAATTTATTAATGAAATTTGCTCTGAGTCTGGAGGGTT
AGTTCAGAGTCTTCCACTGTATGTTTGAAGAGATGACCAAGCCAGAATATGGAATGTTT
ATGTATCTGAAATGGGTTCTGCATGTGGTTTCTGCCAAGCCTAAACCTGAGAAGAAA
AGATATTTCTCTTTGGAATGCTGTGGACTCTCCTTATTCAATTTAAATGTTGTAAC
CTTCTTTCCCACTGGCTCTGTATAAAAAACTTCTGGACCAAAAGCCATCATTGGAAGAT
TTAAAAGAACTCAGTCTCGGTTGGGGAAGAGTTTGAAGAAGTTCTAGATGATGCTGCT
GATGACATTGGAGATGCGCTCTGCATACGCTTTTCTATACACTGGGACCAAAAATGATGTT
GACTTAATTTCAAATGGGATCTCCATACCTGTGGACCAAAACCAAGAGAGACTATGTT
TCTAAGTATATTGATTACATTTTCAACGCTCTGTAAAAGCAGTTTATGAGGAATTTGAG
AGAGGATTTTATAGAGTCTGTGAGAAGGAGATACTTAGACATTTCTACCTGAAGAACTA
ATGACAGCAATCATTGGAATACTGATTATGACTGGAAACAGTTTGAACAGAATTTCAAAG
TATGAGCAAGGATACCAAAAATCACATCCTACTATACAGTTGTTTGGAAAGGCTTTCCAC
AAACTAACCTTGATGAAAAGAAAAAATTCCTCTTTTCTTACAGGACGTGATAGGCTG
CATGCAAGAGGCATACAGAAAATGGAAATAGTATTTGCTGTCTGAACTTTTCAAGTAA
AGAGATCACCAACATCAATAACTTGTGATAATTTCTCTCCCTCCCTAAGTATTCTACA
ATGGAAAGAAATGGAGGAAGCACTTCAAGTAGCCATCAACAACAACAGAGGATTTGTCTCA
CCCATGCTCACACAGTCA

```

Restriction Sites:	Please inquire
ACCN:	NM_001165136
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001165136.1</u> , <u>NP_001158608.1</u>
RefSeq Size:	3781 bp
RefSeq ORF:	2961 bp
Locus ID:	55008
UniProt ID:	<u>Q8IVU3</u>
Cytogenetics:	4q22.1
Protein Families:	Druggable Genome
Gene Summary:	<p>HERC6 belongs to the HERC family of ubiquitin ligases, all of which contain a HECT domain and at least 1 RCC1 (MIM 179710)-like domain (RLD). The 350-amino acid HECT domain is predicted to catalyze the formation of a thioester with ubiquitin before transferring it to a substrate, and the RLD is predicted to act as a guanine nucleotide exchange factor for small G proteins (Hochrainer et al., 2005 [PubMed 15676274]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (2) lacks an alternate in-frame exon in the central coding region, compared to variant 1. The resulting isoform (2) lacks an internal segment, compared to isoform 1.</p>