

Product datasheet for **SC327379**

MAFF (NM_001161573) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAFF (NM_001161573) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAFF
Synonyms:	hMafF; U-MAF
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001161573, the custom clone sequence may differ by one or more nucleotides ATGTCTGTGGATCCCCTATCCAGCAAAGCTCTAAAGATCAAGCGAGAGCTGAGCGAGAAC ACGCCGCACCTGTTCGGACGAGGCGCTGATGGGGCTGTCGGTGCGCGAGCTGAACCGGCAT CTGCGCGGGCTCTCCGCCGAGGAGGTGACACGGCTCAAGCAGCGCGCCGACACTCAAA AACCGTGGCTACGCCGCCAGCTGCCCGTGAAGCGCGTGTGCCAGAAGGAGGAGCTGCAG AAGCAGAAGTCGGAGCTGGAGCGCGAGGTGGACAAGCTGGCGCGGAGAACGCCGCCATG CGCCTGGAGCTCGACGCGCTGCGCGGCAAGTGGAGGCGTGCAGGGCTTCGCGCGCTCC GTGGCCCGCCCGCGGGCCCGCCACGCTCGTGGCGCCGCGCCAGCGTCATCACCATCGTC AAGTCCACCCCGGCTCGGGTCTGGCCCCGCCACGGCCCGGACCCCGCCACGGCCCCG GCCTCCTGCTCC
Restriction Sites:	Please inquire
ACCN:	NM_001161573
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001161573.1, NP_001155045.1</u>
RefSeq Size:	2234 bp
RefSeq ORF:	495 bp
Locus ID:	23764
UniProt ID:	<u>Q9ULX9</u>
Cytogenetics:	22q13.1
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	<p>The protein encoded by this gene is a basic leucine zipper (bZIP) transcription factor that lacks a transactivation domain. It is known to bind the US-2 DNA element in the promoter of the oxytocin receptor (OTR) gene and most likely heterodimerizes with other leucine zipper-containing proteins to enhance expression of the OTR gene during term pregnancy. The encoded protein can also form homodimers, and since it lacks a transactivation domain, the homodimer may act as a repressor of transcription. This gene may also be involved in the cellular stress response. Multiple transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jun 2009]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR compared to variant 3. Variants 1, 3, and 4 encode the same isoform (a).</p>