

Product datasheet for **SC327398**

ABHD14B (NM_001146314) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ABHD14B (NM_001146314) Human Untagged Clone
Tag: Tag Free
Symbol: ABHD14B
Synonyms: CIB; HEL-S-299
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >SC327398 representing NM_001146314.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCAGCAAGCGTGGAGCAGCGGAGGGCACCATCCAGGTGCAGGGCCAGGCCCTCTTCTCCGAGAG
GCCCTGCCGGCAGTGGCAGGCTCGTTCTCTGTACTGCTGCTGCATGGTATTCGCTTCTCCTCCGAG
ACCTGGCAGAACCTGGGTACACTGCACAGGCTGGCCAGGCTGGTACCGGGCTGTGGCCATTGACCTG
CCAGGTCTGGGGCACTCCAAGGAAGCAGCAGCCCTGCCCTATTGGGGAGCTGGCCCTGGCAGCTTC
CTGGCGGCTGTGGTGGATGCCTGGAGCTGGGCCCGGTTGTGATCAGTCCATCACTGAGTGGCATG
TACTCCCTGCCCTTCTCACGGCCCTGGCTCCCAGCTCCCGGGCTTTGTGCCAGTGGCCCCATCTGC
ACTGACAAAATCAATGCTGCCAACTATGCCAGTGTGAAGACTCCAGCTCTGATTGTATATGGAGACCAG
GACCCCATGGGTGAGACCAGCTTTGAGCACCTGAAGCAGCTGCCCAACCACCGGGTGTGATCATGAAG
GGGGCGGGCACCCTGTACCTGGACAAACCAGAGGAGTGGCATAACAGGGCTGTGGACTTCTGCAG
GGGCTCCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001146314

Insert Size: 633 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146314.1](#)

RefSeq Size: 2253 bp

RefSeq ORF: 633 bp

Locus ID: 84836

UniProt ID: [Q96IU4](#)

Cytogenetics: 3p21.2

MW: 22.3 kDa

Gene Summary: Has hydrolase activity towards p-nitrophenyl butyrate (in vitro). May activate transcription.
[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same protein.