

Product datasheet for **SC327417**

SFTPA1 (NM_001164647) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	SFTPA1 (NM_001164647) Human Untagged Clone
Tag:	Tag Free
Symbol:	SFTPA1
Synonyms:	COLEC4; PSAP; PSP-A; PSPA; SFTP1; SFTPA1B; SP-A; SP-A1; SP-A1 beta; SP-A1 delta; SP-A1 epsilon; SP-A1 gamma; SPA; SPA1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC327417 representing NM_001164647. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGGCTGTGCCCTCTGGCCCTCAACCTCATCTTGATGGCAGCCTCTGGTGTGTGTGCGAAGTGAAG
GACGTTTGTGTTGGAAGCCCTGGTATCCCCGGCACTCCTGGATCCCACGGCCTGCCAGGCAGGGACGGG
AGAGATGGTCTCAAAGGAGACCCTGGCCCTCCAGGCCCATGGGTCCACCTGGAGAAATGCCATGTCCT
CCTGAAATGATGGGCTGCCTGGAGCCCTGGTATCCCTGGAGAGTGTGGAGAGAAGGGGAGCCTGGC
GAGAGGGGCCCTCCAGGGCTTCCAGCTCATCTAGATGAGGAGCTCCAAGCCACACTCCACGACTTTAGA
CATCAAATCCTGCAGACAAGGGGAGCCCTCAGTCTGCAGGGCTCCATAATGACAGTAGGAGAGAAGGTC
TTCTCCAGCAATGGGCAGTCCATCACTTTTGTGCCATTGAGGAGCATGTGCCAGAGCAGGCGGCCGC
ATTGCTGTCCAAGGAATCCAGAGGAAAATGAGGCCATTGCAAGCTTCGTGAAGAAGTACAACACATAT
GCCTATGTAGGCCTGACTGAGGGTCCCAGCCCTGGAGACTTCGGTACTCAGACGGGACCCTGTAAAC
TACACCAACTGGTACCGAGGGGAGCCCGCAGGTCGGGGAAAAGAGCAGTGTGTGGAGATGTACACAGAT
GGCAGTGAATGACAGGAACCTGCTGACTCCCGACTGACCATCTGTGAGTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001164647

Insert Size: 747 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164647.1](#)

RefSeq Size: 2159 bp

RefSeq ORF: 747 bp

Locus ID: 653509

UniProt ID: [Q8IWL2](#)

Cytogenetics: 10q22.3

MW: 26.2 kDa

Gene Summary: This gene encodes a lung surfactant protein that is a member of a subfamily of C-type lectins called collectins. The encoded protein binds specific carbohydrate moieties found on lipids and on the surface of microorganisms. This protein plays an essential role in surfactant homeostasis and in the defense against respiratory pathogens. Mutations in this gene are associated with idiopathic pulmonary fibrosis. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010]

Transcript Variant: This variant (4) differs in the 5' UTR compared to variant 1. Variants 1, 3 and 4 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.