

Product datasheet for **SC327418**

CDA11 (TATDN1) (NM_001146160) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDA11 (TATDN1) (NM_001146160) Human Untagged Clone
Tag:	Tag Free
Symbol:	TATDN1
Synonyms:	CDA11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC327418 representing NM_001146160. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGATTACAGGTGAAATCTACAAGACAGTAAAGATGCACTGCATTTGGCACAACAAATGGTATGTTT
TTCAGTACAGTTGGATGTCATCCTACAAGATGTGGTGAATTTGAAAAGAATAACCCTGATCTTTACTTA
AAGGAGTTGCTAAATCTTGCTGAAAACAATAAAGGAAAGTTGTGGCAATAGGAGAATGCGGACTTGAT
TTTGACCGACTGCAGTTTGTCCCAAAGATACTCAACTCAAATATTTGAAAACAGTTTGAAGTGTCA
GAACAAACAAAATTACCAATGTTTCTTCATTGTCGAAACTCACATGCTGAATTTTGGACATAATGAAA
AGAAATAGAGATCGGTGTGTAGGGGAGTGGTGCATTTGATGGTACCAAGGAAGCAGCAGCTGCT
TTGATTGACTTGGATCTTTATATAGGATTAATGGTTGCTCACTGAAAACGAAGCTAATTTGGAAGTT
TTGAAGTCAATTCCTAGTGAAAAATTAATGATTGAGACAGATGCACCTTGGTGTGGAGTCAAAAAGTACA
CATGCTGGATCAAAATATATAAGAACTGCATTTCTACCAAAAAGAAGTGGGAAAGTGGGCACTGCTTA
AAAGACAGAAATGAACCCGCCATATAATTCAAATATTGGAGATAATGTCAGCAGTGAGAGATGAGGAT
CCACTGGAATTAGCCAATACACTATATAACAATACTATTTAAAGTATTTTTCTGGAATATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001146160

Insert Size: 753 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146160.1](#)

RefSeq Size: 930 bp

RefSeq ORF: 753 bp

Locus ID: 83940

UniProt ID: [Q6P1N9](#)

Cytogenetics: 8q24.13

MW: 28.2 kDa

Gene Summary: Putative deoxyribonuclease.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) has multiple differences compared to variant 4. These differences cause translation initiation at a downstream AUG and result in an isoform (b) that is shorter than isoform c.