

## Product datasheet for **SC327429**

### **COPS7A (NM\_001164094) Human Untagged Clone**

#### **Product data:**

<b>Product Type:</b>	Expression Plasmids
<b>Product Name:</b>	COPS7A (NM_001164094) Human Untagged Clone
<b>Tag:</b>	Tag Free
<b>Symbol:</b>	COPS7A
<b>Synonyms:</b>	CSN7; CSN7A; SGN7a
<b>Mammalian Cell Selection:</b>	None
<b>Vector:</b>	<u><a href="#">pCMV6-XL5</a></u>
<b>E. coli Selection:</b>	Ampicillin (100 ug/mL)
<b>Fully Sequenced ORF:</b>	>NCBI ORF sequence for NM_001164094, the custom clone sequence may differ by one or more nucleotides ATGAGTGC GGAAGTGAAGTGACAGGGCAGAACCAGGAGCAATTTCTGCTCCTAGCCAAG TCGGCCAAGGGGCAGCGCTGGCCACACTCATCCATCAGGTGCTGGAGGCCCTGGTGTC TACGTGTTTTGGAGAACTGCTGGACATGCCAATGTTAGAGAGCTGGCTGAGAGTGACTTT GCCTCTACCTTCCGGCTGCTCACAGTGTGGTTATGGGACATACGCTGACTACTTAGCT GAAGCCC GGAATCTTCTCCACTAACAGAGGCTCAGAAGAATAAGCTTCGACACCTCTCA GTTGTCACCTGGCTGCTAAAGTAAAGTGTATCCCATATGCAGTGTGCTGGAGGCTCTT GCCCTGCGTAATGTGCGGCAGCTGGAAGACCTTGTGATTGAGGCTGTGTATGCTGACGTG CTTCGTGGCTCCCTGGACCAGCGCAACCAGCGGCTCGAGGTTGACTACAGCATCGGGCGG GACATCCAGCGCCAGGACCTCAGTGCCATTGCCCGAACCTGCAGGAATGGTGTGTGGGC TGTGAGGTCGTGCTGTCAGGCATTGAGGAGCAGGTGAGCCGTGCCAACCAACACAAGGAG CAGCAGCTGGGCCTGAAGCAGCAGATTGAGAGTGAGGTTGCCAACCTTAAAAAACATT AAAGTTACGACGGCAGCAGCAGCCGACACATCTCAGGACCTGAGCAACACCTGACT GAGCTGAGGGAACCAGCTCCTGGCACCAACCAGCGCCAGCCAGCAAGAAAGCCTCAAAG GGCAAGGGGCTCCGAGGGAGCGCCAAGATTTGGTCCAAGTCAAT
<b>Restriction Sites:</b>	Please inquire
<b>ACCN:</b>	NM_001164094
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001164094.1</a></u> , <u><a href="#">NP_001157566.1</a></u>
<b>RefSeq Size:</b>	1919 bp
<b>RefSeq ORF:</b>	828 bp
<b>Locus ID:</b>	50813
<b>UniProt ID:</b>	<u><a href="#">Q9UBW8</a></u>
<b>Cytogenetics:</b>	12p13.31
<b>Gene Summary:</b>	<p>This gene encodes a component of the COP9 signalosome, an evolutionarily conserved multi-subunit protease that regulates the activity of the ubiquitin conjugation pathway. Alternatively spliced transcript variants that encode the same protein have been described. [provided by RefSeq, Mar 2014]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1-4 encode the same protein.</p>