

Product datasheet for **SC327480**

MTHFSD (NM_001159380) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	MTHFSD (NM_001159380) Human Untagged Clone
Tag:	Tag Free
Symbol:	MTHFSD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC327480 representing NM_001159380. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAATCACAAAATTTAGCTGACTTTCCCGACCTGTTTCATCACAGGATACCCAACCTTTAAGGGTCT
TATCTGGCTTGCCAAAACATCAAAGACCTAGACGTTTTTGCCAGAACACAGGAAGTTAAAGTGGACCT
GATAAACCACTGGAAGCGTTCGGCTGCTGGTGCAGAGCAAAAAACATTGTTGGTTCCAACACCA
CGACTGAGAACGGGATTGTTAATAAGATCACACCACCCCTGGGCAACTAAAGACATCTTGAGAAAA
TGTGCCACCTCTCAGGGTGTGAGGAACACAGTGTCCCATAGGCTTGACTCCAGAGTCCCTCGTGGAT
TTAGTTTGGTGGGATCCGTCGCCGTTTCTGAAAAAGGCTGGAGAATCGGAAGGGAGAAGGCTACGCC
GATCTGGAATATGCCATGATGGTATCCATGGGCGCCGTCAGCAAGGAGACGCCGGTGGTCAACATCGTC
CACGACTGCCAGGTCGTGGACATCCCTGAAGAGCTTGTGAGGAGCACGACATCACTGTGGACTACATC
CTCACTCCAACCAGAGTCATCGCCACAGGCTGCAAGCGCCCAAAGCCAATGGGAATCACCTGGTTCAAG
ATCAGCCTGGAGATGATGGAGAAAATCCCATACTGAGGAGCCTCCGCGCCGAGAGCAGCAGGCTGGG
AAGGATGTCACCCTCCAGGGTGAGCACCAGCACCTCCGGAACCAGGCTGCCAGCAGACAGTGCCTTG
AGTGTGGCAGGAGGCCCGGACACCCGACCAGAAACCAATTCATGGAGGCAGCCCTGGCTCC
CCACCAGGGGAGGGTGCCCGCTTGACGCCGATGTTACGTTGGGAACCTCCCGGGGACGCCCGTGTG
AGTGACCTGAAGAGACCCTGCGGGAACCTCGGCTCCGTCGCCCTGCGGCTCACCTGGCAGGGCCCGCG
CGCAGAGCCTTCCTCCATTACCCGACTCTGCCGAGCCAGCAGGCCGTCCTCTGCTTGCAGGGCCTG
CGCCTGGGACCCGACACCTGAGGGTGGCGCTGGCCAGGCAGCAGAGGGACAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001159380

Insert Size: 1092 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001159380.1](#)

RefSeq Size: 3002 bp

RefSeq ORF: 1092 bp

Locus ID: 64779
Cytogenetics: 16q24.1
MW: 39.8 kDa