

Product datasheet for **SC327489**

LRTM2 (NM_001163926) Human Untagged Clone


Product data:

Product Type: Expression Plasmids
Product Name: LRTM2 (NM_001163926) Human Untagged Clone
Tag: Tag Free
Symbol: LRTM2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >SC327489 representing NM_001163926.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGGCGCCGGGCAGCAGCCCTGGGCAGAGGGGCAGGCTCGCCCTGCAGTGGAGGCAAGTCTCCTGG
ATCACCTGCTGGATCGCCCTGTATGCTGTGGAGGCCCTCCCCACCTGCCCTTTCTCCTGCAAGTGTGAC
AGCCGCAGCCTGGAGGTGGACTGCAGTGGCCTTGGCCTCACACGGTGCCCCAGACGTGCCCGCAGCC
ACCCGAACCTCTTGCTCTTGAACAATAAGCTGAGTGCCTGCCAAGCTGGGCTTTCGCCAACCTCTCC
AGCCTGCAGCGGTTGGACCTGTCCAACAACCTTCTGGACCGGCTGCCCGCTCCATTTTCGGGGACCTG
ACGAATCTGACTGAGCTTCAGCTGCGCAATAACAGCATCAGGACCTGGACAGGGACCTGTGCGGCAC
TCGCCGCTGCTCCGCCACCTGGACCTGTCCATCAACGGCCTGGCCAGTTGCCCCCTGGTCTTTTCGAC
GGGCTCCTGGCTCTGCGCTCCCTCTCGCTTCCGCTCCAACCGTCTGCAGAATCTGGACCGGCTGACATTT
GAACCCCTAGCAAACCTGCAGCTGCTGCAGGTCGGGGATAAACCCCTGGGAGTGTGACTGTAACCTGCGT
GAGTTCAAACACTGGATGGAGTGGTTCTCTACCGAGGGGGACGCTTGGACCAGCTTGCCTGCACCCCTG
CCCAAGGAGCTGAGGGGAAGGACATGCGGATGGTCCCATGGAGATGTTCAACTACTGCTCCAGCTG
GAGGACGAGAATAGCTCAGCTGGCTGGATATTCCTGGGCCACCCTGCACCAAGGCCAGTCCAGAGCCT
GCTAAGCCCAAGCCGGGGCTGAGCCGGAGCCGGAGCCAGCACAGCCTGCCACAGAAGCAGAGGCAC
CGGCCGGCAGCGTGAGGCGAGCCATGGGCACGGTGATCATTGCAGGGGTGCTGTGCGGCGTCTGCTGC
ATCATGATGGTGGTGGCCGCTGCCTATGGCTGCATCTACGCCCTCCTCATGGCCAAGTACCACGGGAG
CTCAAAAAGCGCCAGCCCTGATGGGGACCCGAGGGCGAGCACGAGGACCAGAAGCAGATCTCTTCT
GTGGCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



Plasmid Map:


ACCN: NM_001163926

Insert Size: 1113 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001163926.1](#)

RefSeq Size: 3953 bp

RefSeq ORF: 1113 bp

Locus ID: 654429

UniProt ID: [Q8N967](#)
Cytogenetics: 12p13.33
Protein Families: Transmembrane
MW: 41.2 kDa