

## Product datasheet for SC327521

### BZW2 (NM\_001159767) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	BZW2 (NM_001159767) Human Untagged Clone
Tag:	Tag Free
Symbol:	BZW2
Synonyms:	HSPC028; MST017; MSTP017
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC327521 representing NM_001159767. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC CGATCGCC
ATGAATAAGCATCAGAAGCCAGTGCTAACAGGCCAGCGTTCAAACCTCGGAAAAGGGATGAAAAAGAG
AAATTCGAACCCACAGTCTTCAGGGATACACTTGTCCAGGGGCTTAATGAGGCTGGTGATGACCTTGAA
GCTGTAGCCTAAATTTCTGGACTCTACAGGCTCAAGATTAGATTATCGTCGCTATGCAGACACACTTTC
GATATCCTGGTGGCTGGCAGTATGCTTGCCCTGGAGGAACGCGCATAGATGATGGTGACAAGACCAAG
ATGACCAACCACTGTGTGTTTTAGCAAAATGAAGATCATGAAACCATCCGAAACTATGCTCAGGTCTTC
AATAAACTCATCAGGAGATATAAGTATTTGGAGAAGGCATTTGAAGATGAAATGAAAAAGCTTCTCCTC
TTCCTTAAAGCCTTTTCCGAAACAGAGCAGACAAAGTTGGCGATGCTGTGGGGATTCTGCTGGGCAAT
GGCACCTGCCCGCCACCATCCTCACCAGTCTCTTACCAGACAGCTTAGTCAAAGAAGGCATTGGCGCC
TCATTTGCTGTCAAGCTTTTCAAAGCATGGATGGCAGAAAAAGATGCCAACTCTGTTACCTCGTCTTTG
AGAAAAGCCAACCTTAGACAAGAGGCTGCTTGAACCTCTTCCAGTTAACAGACAGAGTGTGGATCATTTT
GCTAAATACTTCACTGACGCAGGTCTTAAGGAGCTTCCGACTTCTCCGAGTCCAGCAGTCCCTGGGC
ACCAGGAAGGAAGTGCAGAAGGAGCTCCAGGAGCGTCTTCTCAGGAATGCCCGATCAAGGAGGTGGTG
CTTTATGTCAAAGAAGAAATGAAGAGGAATGATCTTCCAGAAACAGCAGTATTGGTCTTCTGTGGACA
TGTATAATGAACGCTGTTGAGTGAACAAGAAGGAAGAACTTGTTCAGAGCAGGCTCTGAAGCACCTG
AAGCAATATGCTCCCCTGCTGGCCGTGTTGAGCTCCCAAGGCCAGTCAGAGCTGATCCTCCTCCAGAAG
GTTCAGGAATACTGCTACGACAACATCCATTTTCATGAAAGCCTTTTCAAAGATTGTGGTTCTCTTTTAT
AAAGCTGATGTTCTGAGCGAAGAAGCAATACTGAAATGGTATAAGGAAGCACATGTTGCTAAAGGCAAA
AGTGTCTTTCTTGACCAGATGAAGAAATTTGTTGAGTGGTTACAAAATGCAGAAGAAGAATCCGAATCG
GAAGGTGAGGAAAAT TAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM\_001159767

Insert Size: 1260 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001159767.1](#)

RefSeq Size: 1883 bp

RefSeq ORF: 1260 bp

Locus ID: 28969

UniProt ID: [Q9Y6E2](#)

Cytogenetics: 7p21.1

MW: 48.2 kDa

Gene Summary: May be involved in neuronal differentiation.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (1) represents the longest transcript. Variants 1, 2 and 4 encode the same isoform (1).