

Product datasheet for **SC327528**

SLC37A4 (NM_001164280) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC37A4 (NM_001164280) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC37A4
Synonyms:	G6PT1; G6PT2; G6PT3; GSD1b; GSD1c; GSD1d; PRO0685; TRG-19; TRG19
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC327528 representing NM_001164280. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAGCCAGGGCTATGGCTATTATCGCACTGTGATCTTCTCAGCCATGTTTGGGGGCTACAGCCTG
TATTACTTCAATCGAAGACCTTCTCCTTTGTCATGCCATCATTGGTGGAAAGAGATCCCTTTGGACAAG
GATGATTTGGGGTTCATCACCAGCAGCCAGTCGGCAGCTTATGCTATCAGCAAGTTTGTGAGTGGGGT
CTGTCTGACCAGATGAGTGTCTGCTGGCTCTTCTCTTCTGGGCTGCTCCTGGTTGGCCTGGTCAACATA
TTCTTTGCTGGAGCTCCACAGTACCTGTCTTTGCTGCCCTCTGGTTCTTAATGGCCTGGCCAGGGG
CTGGGCTGGCCCCATGTGGGAAGTCTCGGGAAGTGGTTTGGCCATCTCAGTTTGGCACTTGGTGG
GCCATCTGTCAACCAGCATGAACCTGGCTGGAGGGCTGGGCCCTATCCTGGCAACCATCCTTGCCAG
AGCTACAGCTGGCGCAGCACGCTGGCCCTATCTGGGGCACTGTGTGGTGTCTCCTTCTCTGTCTC
CTGCTCATCCACAATGAACCTGTGATGTTGGACTCCGCAACCTGGACCCATGCCCTCTGAGGGCAAG
AAGGGCTCCTTGAAGGAGGAGAGCACCTGCAGGAGCTGCTGTGCCCTTACCTGTGGGTGCTCTCC
ACTGGTTACCTTGTGGTGTGGAGTAAAGACCTGCTGTACTGACTGGGGCCAGTTCTTCTTATCCAG
GAGAAAGGACAGTCAGCCCTGTAGGTAGCTCTACATGAGTGCCTGGAAGTTGGGGCCTTGTAGGC
AGCATCGCAGCTGGCTACCTGTGACACCGGCCATGGCAAAGCGGGACTGTCCAACACGGGAACCT
CGCCATGGCCTGTTGCTGTTGATGATGGCTGGCATGACAGTGTCCATGTACCTCTCCGGGTAACAGTG
ACCAGTACTCCCCAAGCTCTGGATCCTGGTATTGGGAGCTGATTTGGTTTCTCCTCGTATGGCCCC
ATTGCCCTGTTTGGAGTCATAGCCAACGAGAGTGGCCCTCCCAACTTGTGTGGCACCTCCACGCCATT
GTGGGACTCATGGCCAATGTGGCGGCTTTCTGGCTGGGCTGCCCTTCCAGCACCATTGCCAAGCACTAC
AGTTGGAGCACAGCCTTCTGGGTGGCTGAAGTATTTGTGCGGCCAGCACGGCTGCCTTCTCTCCTA
CGAAACATCCGCACCAAGATGGCCGAGTGTCCAAGAAGGCTGAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001164280
Insert Size:	1290 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001164280.1</u>
RefSeq Size:	2063 bp
RefSeq ORF:	1290 bp
Locus ID:	2542
UniProt ID:	<u>O43826</u>
Cytogenetics:	11q23.3
Protein Families:	Transmembrane
MW:	46.4 kDa
Gene Summary:	<p>This gene regulates glucose-6-phosphate transport from the cytoplasm to the lumen of the endoplasmic reticulum, in order to maintain glucose homeostasis. It also plays a role in ATP-mediated calcium sequestration in the lumen of the endoplasmic reticulum. Mutations in this gene have been associated with various forms of glycogen storage disease. Alternative splicing in this gene results in multiple transcript variants.[provided by RefSeq, Aug 2009]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR, compared to variant 1. Variants 1, 4, and 5 encode the same protein (isoform 1).</p>