

Product datasheet for **SC327536**

FBXW12 (NM_001159929) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FBXW12 (NM_001159929) Human Untagged Clone
Tag:	Tag Free
Symbol:	FBXW12
Synonyms:	FBW12; FBXO12; FBXO35
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001159929, the custom clone sequence may differ by one or more nucleotides

```

ATGAGCCCCACTCACCAGATTCAAGATCCCAAGCATTGGAATAGGATTGCAGACAGTGAT
TACCTGTGGAGGTCACTATCTCTGCAGAGATGGGACTGCAGCAACTTCACCAATCAACAC
CTGGGCACACACACATGGAAGCAATTTTTCTGCATCAAAGAAGGAAGGAGCTTCGGTTG
GCATTGGCACAGCCGATAAATTATCTACAAAGTAAGTAAAGAACATCGCATTGAGACG
GAGTTGGCTTATCTCTCAGGAAATAGACTTACAGTGGATGAACAGGAGAAATCAATCATT
TGTAAGTATCTCCAAGCAAGAGCTTTGTGCCTGGGATGTGCAAGAGGGTACCATGATC
TGGTCAAGCCCAGTCCAGGAGTCCATTCTCAAATCTGGTAACCCTCCCTCAGATGCAT
CTCGCCATCACTATGGATCGGAAAAAACTATCAAAGTGTGGAAGTGCAGGACAGGGAC
GCTCTGGCTGTTCTCCCATGCCACAGCCCTGTTATTGCATGGAAGCCTATCTTACAAG
GATGGCCATTCCCTGATGGTTGGCGATGCTGCAGGTGACATCTACACATTTACACTGCCT
GGGTTAAGAGATGTTTCTAAAGTACTGCATTTCAATATGGTATTGTACTTCTACTACTGC
TCTCTGACAAGAAATGGGTATTTGCATGTGGGACATACAGTTCGTACCTTGCCACAGGTA
TTCTCTACAGAGTCCCTTACTGAGACCATCAGAAGGCAAGTTCCTCTGTCTACCTTTCTC
CCACATAAATTATGTGCCAGCGCTGCTGGACCCAAAGGTGAAAAACAGGATAACACTG
ATGTCCCAAAGTAGCACTGGAAAAAGACAGAATTTATCACCTTTGATCTAACAACCAAG
AAGACTGGAGGCCAAACAGTCAATCAAGCATATGAGATCGCAAGTTTCCAGGTGGCAGCT
CATCTGAAGTGCCTATCTGGATGGGAGCCAGTATGGATATGATTGTCTTTACCAGT
GGACCATACTTGTACTCTTCAAGCATCACTGGCTTCTCTGCTGCAACGATTTGAGGACCAT
CAGGCAGCCATCAACAATTCTGGGTGGATCCTTGCTATGTGCTCACCACATCCGAGAAC
TCTGTGCACGTGTACATGTGGGAAGAAGGAGGCCCATCCATACCTCAGGAGCTGCTGT
CACCTGGAAAACACGTGGCATGATCACACAACAGACAGCTGCATCTCCAGTGTGATGTGT
GATAATGCAAGCATAGTACTTAGGGTGAGGAAAGTAAGTGACTCCAGCATTCTGGTGATG
TATTCTTTGAATACG

```

Restriction Sites: Please inquire



[View online »](#)

ACCN:	NM_001159929
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001159929.1, NP_001153401.1</u>
RefSeq Size:	1473 bp
RefSeq ORF:	1338 bp
Locus ID:	285231
UniProt ID:	<u>Q6X9E4</u>
Cytogenetics:	3p21.31
Protein Families:	Druggable Genome
Gene Summary:	<p>Members of the F-box protein family, such as FBXW12, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603034), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (3) has a distinct N-terminus and is shorter than isoform 1.</p>