

Product datasheet for **SC327545**

GABA A Receptor alpha 5 (GABRA5) (NM_001165037) Human Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | GABA A Receptor alpha 5 (GABRA5) (NM_001165037) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | GABRA5 |
| Synonyms: | DEE79; EIEE79 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

Fully Sequenced ORF: >SC327545 representing NM_001165037.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGACAATGGAATGTTCTCTGGTTTTATCATGATCAAAAACCTCCTTCTCTTTTGTATTCCATGAAC
TTATCCAGTCACTTGGCTTTTACAGATGCCAACCAAGTTCAGTAAAAGATGAGACCAATGACAACATC
ACGATATTTACCAGGATCTTGGATGGGCTCTTGGATGGCTACGACAACAGACTTCGGCCCCGGCTGGGA
GAGCGCATCACTCAGGTGAGGACCGACATCTACGTCACCAGCTTCGGCCCCGGTGTCCGACACGAAATG
GAGTACACCATAGACGTGTTTTTCCGACAAAGCTGAAAAGATGAAAGGCTTCGGTTAAGGGGCCCATG
CAGCGCTCCCTCTCAACAACCTCCTTCCAGCAAGATCTGGACCCAGACACGTTCTTCCACAACGGG
AAGAAGTCCATCGCTCACAACATGACCACGCCCAACAAGCTGCTGCGGCTGGAGGACGACGGCACCTG
CTCTACACCATGCGCTTGACCATCTCTGCAGAGTGCCCATGCAGCTTGAGGACTTCCCGATGGATGCG
CAGCTTGGCCTCTGAAATTTGGCAGCTATGCGTACCCTAATTCTGAAGTCGTTTACGTCTGGACCAAC
GGCTCCACCAAGTCGGTGGTGGTGGCGGAAGATGGCTCCAGACTGAACCAGTACCACCTGATGGGGCAG
ACGGTGGGCACTGAGAACATCAGCACCAGCACAGGCGAAATACACAATCATGACAGCTCACTTCCACCTG
AAAAGGAAGATTGGCTACTTTGTATCCAGACCTACCTTCCCTGCATAATGACCGTGATCTTATCAG
GTGTCCTTTTGGCTGAACCGGAATCAGTCCCAGCCAGGACAGTTTTTGGGGTACCACGGTGTGACC
ATGACGACCTCAGCATCAGCGCCAGGAATCTCTGCCAAAAGTGGCCTACGCCACCGCCATGGACTGG
TTCATAGCCGTGTGCTATGCCTTCGTCTTCTCGGCGCTGATAGAGTTTCCACGGTCAATTACTTTACC
AAGAGAGGCTGGGCTGGGATGGCAAAAAGCCTTGAAGCAGCCAAGATCAAGAAAAGCGTGAAGTC
ATACTAAATAAGTCAACAACGCTTTTACAACCTGGGAAGATGTCTACCCCCAAACATTCGAAGGAA
CAGACCCAGCAGGGACGTCGAATACAACCTCAGTCTCAGTAAAACCTCTGAAGAGAAGACTTGTAA
AGCAAAAAGACTTACAACAGTATCAGCAAAAATTGACAAAATGTCCCGAATCGTATTCCCAGTCTTGTC
GGCACTTTCAACTTAGTTTACTGGCAACGATTTTGAATAGGGAGCCGGTGATAAAAGGAGCCGCTCT
CCAAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_001165037

Insert Size: 1389 bp

| | |
|-------------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| OTI Annotation: | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001165037.1</u> |
| RefSeq Size: | 2776 bp |
| RefSeq ORF: | 1389 bp |
| Locus ID: | 2558 |
| UniProt ID: | <u>P31644</u> |
| Cytogenetics: | 15q12 |
| Protein Families: | Druggable Genome, Ion Channels: Cys-loop Receptors, Transmembrane |
| Protein Pathways: | Neuroactive ligand-receptor interaction |
| MW: | 52.1 kDa |
| Gene Summary: | <p>GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA-A receptors, which are ligand-gated chloride channels. Chloride conductance of these channels can be modulated by agents such as benzodiazepines that bind to the GABA-A receptor. At least 16 distinct subunits of GABA-A receptors have been identified. Transcript variants utilizing three different alternative non-coding first exons have been described. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (2) represents use of an alternate promoter and 5' UTR, compared to variant 1. Both variants 1 and 2 encode the same protein.</p> |