

## Product datasheet for **SC327551**

### Cannabinoid Receptor I (CNR1) (NM\_001160258) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cannabinoid Receptor I (CNR1) (NM_001160258) Human Untagged Clone
Tag:	Tag Free
Symbol:	CNR1
Synonyms:	CANN6; CB-R; CB1; CB1A; CB1K5; CB1R; CNR
Mammalian Cell Selection:	None
Vector:	<u><a href="#">pCMV6-XL5</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001160258, the custom clone sequence may differ by one or more nucleotides

```

ATGAAGTCGATCCTAGATGGCCTTGCAGATACCACCTTCCGCACCATCACCCTGACCTC
CTGTACGTGGGCTCAAATGACATTACAGTACGAAGACATCAAAGGTGACATGGCATCCAAA
TTAGGGTACTTCCACAGAAATCCCTTTAACTTCCTTTAGGGGAAGTCCCTTCCAAAG
AAGATGACTGCGGGAGACAACCCAGCTAGTCCAGCAGACCAGGTGAACATTACAGAA
TTTTACAACAAGTCTCTCTCGTCTTCAAGGAGAATGAGGAGAACATCCAGTGTGGGGAG
AACTTCATGGACATAGAGTGTTCATGGTCTGAACCCAGCCAGCAGCTGGCCATTGCA
GTCCTGTCCCTCACGCTGGGCACCTTACGGTCTGGAGAACCTCCTGGTGTGTGCGTC
ATCCTCCACTCCCGCAGCCTCCGCTGCAGGCCTTCTACCACTTCATCGGCAGCCTGGCG
GTGGCAGACCTCCTGGGAGTGCATTTTTGTCTACAGCTTCACTGACTCCACGTGTTCC
CACCGCAAAGATAGCCGCAACGTGTTTCTGTTCAAAGTGGGTGGGGTACGGCCTCCTTC
ACTGCCTCCGTGGGCAGCCTGTTCTCACAGCCATCGACAGGTACATATCCATTCACAGG
CCCCTGGCCTATAAGAGGATTGTACCAGGCCAAGGCCGTGGTGGCGTTTTGCCTGATG
TGGACCATAGCCATTGTGATCGCCGTGCTGCCTCTCCTGGGCTGGAAGTCCGAGAAACTG
CAATCTGTTTGTCTCAGACATTTCCACACATTGATGAAACCTACCTGATGTTCTGGATC
GGGGTACCAGCGTACTGCTTCTGTTTCATCGTGTATGCGTACATGTATATTTCTCTGGAAG
GCTCACAGCCACGCCGTCCGCATGATTCAGCGTGGCACCAGCAAGCCGCATGGCATTAGG
TTAGCCAAGACCCTGGTCTGATCCTGGTGGTGTGATCATCTGCTGGGGCCCTCTGCTT
GCAATCATGGTGTATGATGTCTTTGGGAAGATGAACAAGCTCATTAAAGACGGTGTGTTGCA
TTCTGCAGTATGCTCTGCCTGCTGAACTCCACCGTGAACCCATCATCTATGCTCTGAGG
AGTAAGGACCTGCGACACGCTTTCCGGAGCATGTTTCCCTCTTGTGAAGGCACTGCGCAG
CCTCTGGATAACAGCATGGGGACTCGGACTGCCTGCACAAACACGCAAACAATGCAGCC
AGTGTTCACAGGGCCGAGAAAGCTGCATCAAGAGCACGGTCAAGATTGCCAAGGTAACC
ATGTCTGTGCCACAGACAGTCTGCCGAGGCTCTG

```



[View online >](#)

<b>Restriction Sites:</b>	Please inquire
<b>ACCN:</b>	NM_001160258
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001160258.1</a></u> , <u><a href="#">NP_001153730.1</a></u>
<b>RefSeq Size:</b>	5901 bp
<b>RefSeq ORF:</b>	1419 bp
<b>Locus ID:</b>	1268
<b>UniProt ID:</b>	<u><a href="#">P21554</a></u>
<b>Cytogenetics:</b>	6q15
<b>Protein Families:</b>	Druggable Genome, GPCR, Transmembrane
<b>Protein Pathways:</b>	Neuroactive ligand-receptor interaction
<b>Gene Summary:</b>	<p>This gene encodes one of two cannabinoid receptors. The cannabinoids, principally delta-9-tetrahydrocannabinol and synthetic analogs, are psychoactive ingredients of marijuana. The cannabinoid receptors are members of the guanine-nucleotide-binding protein (G-protein) coupled receptor family, which inhibit adenylyate cyclase activity in a dose-dependent, stereoselective and pertussis toxin-sensitive manner. The two receptors have been found to be involved in the cannabinoid-induced CNS effects (including alterations in mood and cognition) experienced by users of marijuana. Multiple transcript variants encoding two different protein isoforms have been described for this gene. [provided by RefSeq, May 2009]</p> <p>Transcript Variant: This variant (4) contains two alternate exons in the 5' UTR, compared to variant 1. Variants 1, 3, 4 and 5 all encode isoform a. This variant was designated CB1C by PubMed ID: 15289816.</p>