

Product datasheet for **SC327560**

STAU2 (NM_001164384) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	STAU2 (NM_001164384) Human Untagged Clone
Tag:	Tag Free
Symbol:	STAU2
Synonyms:	39K2; 39K3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC327560 representing NM_001164384.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTTCAAATAAATCAGATGTTCTCAGTGCAGCTGAGTCTTGGTGAAGCAGACATGGGAATCCGAAGGC
AGCAGTATAAAGAAGGCTCAGCAGGCTGTTGCCAATAAAGCTTTGACTGAATCTACGCTTCCCAAACCA
GTTCAGAAGCCACCCAAAAGTAATGTTAACAATAACCCAGGCAGTATAACTCCAAGTGTGGAAGTGAAT
GGGCTTGCTATGAAAAGGGGAGAGCCTGCCATCTACAGGCCATTAGATCCAAAGCCATTCCAAATTAT
AGAGCTAATTACAACCTTCGGGGCATGTACAATCAGAGGTATCATTGCCAGTGCCTAAGATCTTTTAT
GTTCAGCTCACTGTAGGAAATAATGAATTTTTGGGGAAGGAAAGACTCGACAAGCTGCTAGACACAAT
GCTGCAATGAAAGCCCTCAAGCACTGCAGAATGAACCTATTCCAGAAAGATCTCCTCAGAATGGTGAA
TCAGGAAAGGATGTGGATGATGACAAAGATGCAAATAAGTCTGAGATCAGCTTAGTGTGTTGAAATTGCT
CTGAAGCGAAATATGCCTGTCAGTTTTGAGGTTATTAAGAAAGTGGACCACCACATATGAAAAGCTTT
GTTACTCGAGTGTGAGTAGGAGTTCCTGCAGAAGGAGAAGGAAATAGCAAAAACTCTCCAAGAAG
CGCGCTGCGACCACCGTCTTACAGGAGCTTAAAAAAGTCCACCTCTTCCGTGGTGAAAGGCCAAAA
CTATTTTTTAAAAACGCCCTAAAAACAATAGTAAAGGCCGGACCAGAATATGGCCAAGGGATGAACCT
ATTAGCCGCTGGCGCAAATTAACAGGCCAAAAAGGAAAAGGAGCCGGATTATGTTTTGCTTTCAGAA
AGAGGAATGCCTCGAGCTCGAGAATTTGTGATGCAGGTGAAGGTAGGCAATGAAGTTGCTACAGGAACA
GGACCTAATAAAAAGATAGCCAAAAAATGCTGCAGAAGCAATGCTGTTACAACCTGGTTATAAAGCA
TCCACTAATCTTCAGGATCAACTTGAGAAGACAGGGGAAAACAAGGATGGAGTGGTCCAAGCCTGGG
TTTCTGAACCAACAAATAACTCCTGCACTACTCTAGGCTATTTGTACCCAAAGATGAACCAA
CCTTCAAGCTCTTCTTCAAGTATATCTCCACATCGAATAGTTCAGCTACAATTGCCAGGAACTCCTT
ATGAATGGAACATCTTCTACAGCTGAAGCCATAGGTTTAAAAGGAAGTTCTCCTACTCCCTTGTCT
CCAGTACAACCTTCAAAACAAGTGAATATTTAGCAAGGATTCAAGGCTTTCAGGTATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001164384

Insert Size: 1440 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001164384.1
RefSeq Size:	4225 bp
RefSeq ORF:	1440 bp
Locus ID:	27067
UniProt ID:	Q9NUL3
Cytogenetics:	8q21.11
Protein Families:	Transcription Factors
MW:	52.7 kDa
Gene Summary:	<p>Staufen homolog 2 is a member of the family of double-stranded RNA (dsRNA)-binding proteins involved in the transport and/or localization of mRNAs to different subcellular compartments and/or organelles. These proteins are characterized by the presence of multiple dsRNA-binding domains which are required to bind RNAs having double-stranded secondary structures. Staufen homolog 2 shares 48.5% and 59.9% similarity with drosophila and human staufen, respectively. The exact function of Staufen homolog 2 is not known, but since it contains 3 copies of conserved dsRNA binding domain, it could be involved in double-stranded RNA binding events. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2009]</p> <p>Transcript Variant: This variant (6) differs in the 5' and 3' UTRs and coding sequence compared to variant 1. The resulting isoform (e) has shorter and distinct N- and C-termini compared to isoform a. Variants 5 and 6 both encode the same isoform (e).</p>