

Product datasheet for **SC327593**

TJAP1 (NM_001146016) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TJAP1 (NM_001146016) Human Untagged Clone
Tag:	Tag Free
Symbol:	TJAP1
Synonyms:	PILT; TJP4
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF:	<p>>NCBI ORF sequence for NM_001146016, the custom clone sequence may differ by one or more nucleotides</p> <pre> ATGACTAGTGCCGCCCTGCTAAGAAACCTACCGTAAGGCACCACCAGAGCATCGGGAG CTGCGTTTGAAATTCCTGGATCCCGGCTTGAGCAGGAGGAACCCCTGACTGATGCAGAA AGGATGAAGCTCTTACAGGAGGAGAATGAAGAGCTTCGCCGGCGCCTGGCCTCCGCCACC AGACGCACTGAGGCCCTGGAACGTGAGCTGGAATTTGGCAGGACTGCCTGGAGCTGGAG CTGGGCCAGAGCCGCGAGGAGCTGGACAAATTTAAGGATAAGTTCCGCAGGCTGCAGAAC AGCTACACGGCTTCCCAGCGGACCAACCAGGAGTTGGAGGACAAGCTGCACACACTGGCC TCTCTTAGCCACAGCTGGATTTTGAATCAAGAAGGCTGAGATGGATAGGAAGACGCTG GACTGGGAGATTGTGGAGCTGACCAACAAGCTGTGGATGCCAAGAACACCATCAACAAG CTGGAAGAGCTCAATGAGCGGTACCGGCTGGACTGCAACCTGGTGTACAGCTCCTCAAG TGCAACAAGTCCCCTCCGAAACCACAAGTTTGGCGATCTGCCCTGTGAGCTACAGGAC ATGTTTCGGAAACATTTGCACAGTGGTCAAGAGGCCGCCAGCCAGGTCTGCTCCCAGC CTAGCCCCAGGGCTGTGGTGCCTACCTCAGTCATTGCCCGAGTGTAGAGAAGCCGGAG TCTCTACTGCTCAATTCAGCCCAGTCAGGCAGCGCCGGGCGCCCTTGGCTGAGGATGTC TTTGTGCATGTGGACATGAGTGAGGGTGTCCCAGGTGATCCAGCCAGTCCCCCGGCCCT GGCAGCCCCACCCACAACCCAATGGGGAGTGCCACTCTCTGGTACTGCCAGGGGCTCC CCGGAGGAAGAGCTGCCCTGCCAGCCTTTGAGAAGCTGAACCCCTACCCAACCCCGTCT CCACCACCCACTGTATCCTGGCCGAGGGTAATAGAGTTCTCTGAGGATAAGGTTTCGG ATCCCCCGCAACAGCCCCCTGCCAACTGCACTTACGCTACCCGCCAGGCCATTTCCCTG AGCCTGGTAGAGGAGGGGAGTGAGCGGGCCCGCCAGCCAGTCCAGCAGCACCCTGCC TCAGCCCAGGCCACCCACCAGCCAGCCAGCACCCTAACACTCAGTGCCCCA GCTAGCTCTGCCAGCTCTGAAGAGGACCTGCTGGTCAGCTGGCAGCGGGCATTGTGGAC CGTACTCCACCACCTGCTGTGGCCAGCGCACAGCCTTTGGACGCGATGCCCTCCCT GAGCTGCAGCGCCATTTGCCCATAGCCCCGCTGACAGAGATGAGGTGGTCCAGGCACCT TCTGCCGACCCGAAGAGAGTGAGCTTTTGTACCCACAGAACCTGACTCTGGCTTTCC AGGGAGGAAGAAGAGCTGAACCTGCCTATCAGTCTGAGGAAGAGCGCCAGAGCTGCTG CCCATTAACAGGGGCACAGAGGAGGGGCCAGGCACTTCCCACACCAGGGCAGGGCTGG CCACTCCCAGCTCCAGTCGCCCCAGCGCAGCCCAAGAGGATGGGGTTTACCACCTG CACCGCAAGGACAGCCTGACCCAGGCCAGGAGCAGGGCAACCTGCTCAAC </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001146016
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146016.1](#), [NP_001139488.1](#)

RefSeq Size: 2761 bp

RefSeq ORF: 1674 bp

Locus ID: 93643

UniProt ID: [Q5JTD0](#)

Cytogenetics: 6p21.1

Protein Pathways: Tight junction

Gene Summary: This gene encodes a tight junction-associated protein. Incorporation of the encoded protein into tight junctions occurs at a late stage of formation of the junctions. The encoded protein localizes to the Golgi and may function in vesicle trafficking. Alternatively spliced transcript variants have been described. A related pseudogene exists on the X chromosome. [provided by RefSeq, Mar 2009]
Transcript Variant: This variant (1) encodes the longer isoform (a). Variants 1, 2, and 8-15 all encode the same isoform (a).