

Product datasheet for **SC327598**

ZNF583 (NM_001159861) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF583 (NM_001159861) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF583
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC327598 representing NM_001159861.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCCAAGGATTTGGTGACATTTGGGGATGTGGCTGTAATTTCTCTCAAGAGGAATGGGAATGGCTG
AACCTGCTCAGAGGAATTTGTACAGGAAAGTGATGTTGGAGAACTACAGGAGCTTGGTATCATTGGGA
GTTTCTGTTTCTAAGCCAGATGTGATCTCATTATTGGAGCAAGGAAAAGGCCCTGGATGGTGAAGAAG
GAGGGAACAAGAGGCCCATGCCCTGATTGGGAGTATGTATTTAAAAACAGTGAATTTTCATCAAAGCAA
GAGACATATGAAGAATCATCCAAAGTTGTGACAGTGGGAGCAAGACATCTTAGTTATAGCCTTGACTAT
CCCAGTTTGAGAGAAGACTGTCAAAGTGAGGACTGGTATAAGAACCAGCTGGGAAGTCAAGAGGTACAT
CTTAGTCAATTAATCATCACTCATAAAGAAATCCTTCCAGAAGTCAAATAAAGAATATAACAAATCT
TGGCAAACATTCCACCAGGATACAATCTTTGATATACAACAGAGTTTTCCACCAAAGAAAAGGCACAT
AAGCATGAACCACAAAAGAAAAGTTACCGAAAAAATCTGTTGAAATGAAACATAGGAAAGTCTATGTA
GAAAAGAACTTTGAAATGTAATGATTGTGAGAAAGTCTCAACCAGAGCTCATCCCTACTCTTCAT
CAGAGAATTCATACTGGAGAGAAAACCTATGCATGTGTTGAATGTGGGAAAACGTTACGCCAGAGTGCA
AACTTGGCGCAACATAAGAGAATACATACTGGAGAGAAAACCTATGAATGTAAGAATGTAGGAAAAGCC
TTCAGCCAGAATGCACACCTGGCCCAACATCAGAGAGTTCATACTGGAGAGAAAACCTTATCAGTGTAAA
GAATGTAAAAAGCCTTCAGCCAGATTGCACACCTGACTCAGCATCAGAGAGTTCATACTGGAGAGAGA
CCTTTCGAATGTATTGAATGTGGAAAGGCCTTTAGTAATGGTTCATTTCTTGCTCAGCATCAGAGAATT
CATAACAGGAGAGAAAACCTTATGTGTGAATGTGTGTGGGAAAAGCCTTTAGCCATCGTGGATACCTAAT
GTACATCAGAGAATTCATACTGGAGAGAGACCCTACGAATGTAAAGAAATGTAGGAAAAGCCTTCAGCCAG
TATGCACACCTTGCTCAACATCAGAGAGTTCATACTGGAGAAAACCTTATGAATGTAAGATGTAGG
AAAGCCTTCAGCCAAATTGCATACCTTGATCAACATCAGAGGGTTCATACTGGAGAGAAAACCTATGAA
TGTATTGAATGTGGGAAGGCCTTTAGCAATAGTTCATCACTTGCACAACATCAGAGAAGTCACTAGGA
GAAAAACCTATATGTGAAGGAATGTAGGAAAACATTTAGCCAGAATGCAGGCCTTGCTCAACATCAG
AGAATTCATACTGGAGAGAAAACCTTATGAATGTAATGTTTGTGGGAAAAGCATTAGCTATAGTGGATCT
CTTACTCTACATCAGAGAATTCATACTGGAGAAAAGACCCTATGAATGTAAGATTGCAGGAAATCTTTC
AGGCAGCGTGACATCTTGCTCATCATGAGAGAATTCATACTATGGAGTCATTCTTGACTCTTTCCTCT
CCCTCACCTCCACATCAATCAGTTGCCAAGACCTGTAGGTTTCATCTCCTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN:	NM_001159861
Insert Size:	1710 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001159861.1</u>
RefSeq Size:	2543 bp
RefSeq ORF:	1710 bp
Locus ID:	147949
UniProt ID:	<u>Q96ND8</u>
Cytogenetics:	19q13.43
Protein Families:	Transcription Factors
MW:	66 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same protein.