

## Product datasheet for **SC327606**

### **BEGAIN (NM\_001159531) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	BEGAIN (NM_001159531) Human Untagged Clone
Tag:	Tag Free
Symbol:	BEGAIN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC327606 representing NM\_001159531.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGGAGAACTCAGCGCGCTGCAGGAGCAGAAGGGCGAGCTGCGCAAGCGGCTGTCTACACCACACAC
AAGCTCGAGAAGCTCGAGACCGAGTTCGACTCCACGCGCCACTACCTGGAGATCGAGCTGCGGCGCGCG
CAGGAGGAACTGGAGAAGGTCACGGAGAAGCTGCGCAGGATTCAGAGCAACTACATGGCACTGCAGAGG
ATCAACCAGGAGCTGGAGGACAAGCTGTACCGCATGGGCCAGCACTATGAGGAGGAGAAGCGTGCGCTG
AGCCACGAGATTGTTGCCCTCAACAGCCATCTGCTAGAAGCCAAGGTGACCATCGACAAGCTGTGAGAG
GACAACGAGCTCTATAGGAAGGACTGCAATCTAGCGGCCAGCTGCTGCAGTGCAGCCAGACCTACGGC
AGGGTCCACAAGGTGTCTGAGCTGCCCTCGATTTCAGGAGCGCGTGAGCCTGCACATGGAGAAGCAC
GGCTGCAGCTGCCATCCCCGCTCTGCCACCCGGCCTACGCCGACAGCGTCCCCACCTGCGTCATTGCC
AAGGTGCTGGAGAAGCCGACCCCGCCAGCCTGTCTCCCGCCTGTCGATGCCTCCGCCCGCGACCTG
GCCTTCTGCGACGGGGTGGAGAAACCAGGCCCGCGGCCCCCTACAAGGGAGACATCTACTGCAGTGAC
ACAGCCCTCTACTGCCCGGAGGAGCGGCCGCGAGACCGCGCGCCTAGCGTGGACGCGCCCGTGACCGAC
GTGGGCTTCTGCGGGCCAGAACTCCACTGACAGCGCGGCCGAGGAGGAGGAGGCCAGGGCGCGC
GCCTTCCCGCGGGCTTCCAGCATGAGGCCTTCCCCAGCTACGCAGGCTCACTGCCACGTCCAGCTCC
TACTCCAGCTTACGCGCCACGTCGGAGGAGAAGGAGCAGCGCAGGCCAGCACGCTGACCGCGTCGAG
CAGGCCATCTACCTGAACAGCCGCGACGAGCTCTCGACCGCAAGCCACCCGCCACCACCTACGAGGGC
AGCCCTCGCTTTGCCAAGGCCACGGCCGCGGTGGCGGCCCGCTGGAGGCCGAAGTGGCCCCAGGCTTC
GGGCGGACCATGTCACCGTACCCGGCCGAGACCTCCGCTTCCCGCCTCTCCGGTCCCAGCAGGCC
CTGATGCCCCAAACCTGTGGAGCTGCGGGCCAAGCCGGGGACCGCCCGCTCCCAGGAGGACATG
AGGGGCCAGTGGCGTCCCCTGAGCGTGGAGGACATCGGCCCTACTCCTACCCCGTGAGCGCTGCCGCG
CGCGCCTCACCTGCAGTCTCTGAACGCTACTACGGCGGGCCGGGGCAGCCCGGGCAAGAAGGCC
GACGGCCGCGCCAGCCGCTCTACGCCAGCTACAAGGCCGACAGCTTCTCCGAGGGGACGACCTCTCC
CAGGGCCACCTGGCAGAGCCCTGCTTCTGCGGGCGGGCGGACCTGAGCCTCAGTCCCGGCCGCTCG
GCTGACCCACTGCCCGGCTATGCACCCAGCGAGGGGGGGACGGGGACAGGCTCGGGGTGCAGCTGTGT
GGGACCGCCAGCAGCCCTGAGCCCGAGCAGGGTCCAGGGACTCCTTGGAGCCGAGCTCCATGGAGGCC
TCCCCGAAATGCATCCTGCCGCCGCTCAGCCCCAGCAGGCCTTCCCGGACTGGTGGCTCGGGG
CTGAGCCGAAGGACAGCCTACCAAGGCCAGCTCTACGGAACCTTGCTCAACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001159531

**Insert Size:** 1782 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001159531.1](#)

**RefSeq Size:** 2763 bp

**RefSeq ORF:** 1782 bp

**Locus ID:** 57596

**UniProt ID:** [Q9BUH8](#)

**Cytogenetics:** 14q32.2

**Protein Families:** Druggable Genome

**MW:** 64.8 kDa

**Gene Summary:** May sustain the structure of the postsynaptic density (PSD).[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (1) represents the longer transcript. Variants 1 and 2 encode the same protein.