

Product datasheet for **SC327622**

RASGRP 4 (RASGRP4) (NM_001146202) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RASGRP 4 (RASGRP4) (NM_001146202) Human Untagged Clone
Tag:	Tag Free
Symbol:	RASGRP4
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF:	<p>>NCBI ORF sequence for NM_001146202, the custom clone sequence may differ by one or more nucleotides</p> <pre> ATGAACAGAAAAGACAGTAAGAGGAAGTCCCACCAGGAATGCACCGGAAAAATAGGAGGG CGAGGCCGGCCCCGCAAGTGCGCCGCCACAAGACATGCCCCAGCCCTCGGGAATCAGC AAGGTCATGGCTTCCATGAACCTGGGCCTGCTGAGTGAGGGCGGCTGCAGCGAAGATGAG CTGCTGGAGAAATGCATCCAGTCCTTCGATTAGCTGGCAGCCTGTGCCACGAGGACCAC ATGCTCAACATGGTGCTGGCCATGCACAGCTGGGTGCTGCCGTCGCGCCGACCTGGCTGCC CGCTGTGACCTCATACCAGAAGGCCACAGGGGACACCCAGGAGCTGAGACGGCTGCAG ATCTGTACCTGGTCAGGTAAGTGGCTGATGCGACACCCTGAGGTGATGCACCAGGATCCC CAGCTAGAAGAAGTCATAGGTCGTTTCTGGGCCACCGTGGCCCGGGAGGGCAACTCAGCC CAGAGAAGACTGGGAGACTCTTCTGACCTCCCAGGCTGGGCAAAAAGCGCAAAGTGTC TTGCTTTTCGACCACTTGAGACGGGGAGCTGGCTCAGCACCTCACCTACCTGGAGTTC CGGTCCTCCAGGCTATCAGCCCCAGGACCTGCGGAGCTACGTTTTGCAGGGCTCAGTA CGAGGCTGCCCGCCCTGGAGGGCTCCGTAGGTCTCAGCAACAGCGTGTCCCGCTGGGTG CAGGTGATGGTGCTGAGCCGTCGCGGGCCCTACAGCGTGCACAGGTGCTGGACAAGTTC ATTCACGTGGCAGAGGCTCCACCAGCTGCAGATTTCAACACGCTGATGGCAGTCACA GGGGCCCTGTGTACAGTGCCATCTCCAGACTCAAGGACTCCCATGCCACCTGAGCCCT GACAGCACCAAGGCCCTCCTGGAGCTCACTGAGCTCCTTGCCCTCCACAACAACACTACGCC CGTACCGCCGCACCTGGGCTGGCTGCGCGGGTTTCCGGCTGCCTGTACTGGGCGTGCAC CTCAAGGACCTGGTGCCCTGCATGAGGCACAGCCCGACAGGTTGCTGACGGCCGCTG CACCTACCAAGCTGAACAACCTCTACCTGCGGCTGCAGGAGCTGGTGGCCCTCCAAGGG CAGCATCCACCCTGCAGCGCCAATGAGGATCTGCTGCACCTGCTCAGCTCTCCCTGGAC CTCTTCTACACGGAAGACGAGATCTATGAGCTTTCTTATGCCCGGGAGCCGCTTGTCCC AAGAGCCTGCCACCCTCCCCTTCAATGCACCTCTGGTGGTGGAGTGGGCCCTGGTGTG ACACCAAGCCGACAGGGTCACTGCGGTGCGCATGTGGAGCAGCTGGTGGAGTCTGTG TTCAAGAATTATGACCCTGAAGGCCGAGGAACAATCTCTCAGGAGGACTTTGAGCGACTC TCGGGCAATTTCCCTTCGCTGCCATGGGCTTCAACCCACCCACGCCAGGGGAGAGGA TCCTTCAGCAGAGAGGAGCTGACAGGGTACCTGCTCCGGGCCAGCGCCATCTGCTCCAAG TTGGCCCTGGCCTTCTGCACACCTTCCATGAGGTACCTTCCGAAAGCCTACCTTCTGC GACAGCTGCAGTGGCTTCTCTGGGTGTCACCAAGCAAGGCTACCGCTGTCGGGAGTGC GGGCTGTGTTGCCACAACACTGCAGAGACCAGGTGAAGGTAGAATGTAAGAAGAGGCCA GGGGCCAAGGGCGATGCAGGACCCCGGAGCTCCTGTCCATCCACACCAGCTCCCAT GCCAGCTGTGGTCCGAGGAAAATCACTCTACACGCTATCCCTGGAGCCTGAGACTGGG TGCCAGCTTCGCCATGCCTGGACCCAGACTGAATCCCCACACCCTTCTGGGAAACAGAT ACGGTCCCCTGCCGGTATGGACCCACCATCAACTGCATCCTCCAAGCTGGATTCC </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001146202
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146202.1](#), [NP_001139674.1](#)

RefSeq Size: 3193 bp

RefSeq ORF: 1980 bp

Locus ID: 115727

UniProt ID: [Q8TDF6](#)

Cytogenetics: 19q13.2

Protein Families: Druggable Genome

Protein Pathways: MAPK signaling pathway

Gene Summary: The protein encoded by this gene is a member of the Ras guanyl nucleotide-releasing protein (RasGRP) family of Ras guanine nucleotide exchange factors. It contains a Ras exchange motif, a diacylglycerol-binding domain, and two calcium-binding EF hands. This protein was shown to activate H-Ras in a cation-dependent manner in vitro. Expression of this protein in myeloid cell lines was found to be correlated with elevated level of activated RAS protein, and the RAS activation can be greatly enhanced by phorbol ester treatment, which suggested a role of this protein in diacylglycerol regulated cell signaling pathways. Studies of a mast cell leukemia cell line expressing substantial amounts of abnormal transcripts of this gene indicated that this gene may play an important role in the final stages of mast cell development. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2009]

Transcript Variant: This variant (b) uses an alternate in-frame splice junction compared to variant a. The resulting isoform (b) has the same N- and C-termini but is shorter compared to isoform a.