

Product datasheet for **SC327661**

Prominin 2 (PROM2) (NM_001165977) Human Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Prominin 2 (PROM2) (NM_001165977) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | PROM2 |
| Synonyms: | PROML2 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |




[View online »](#)

Fully Sequenced ORF: >SC327661 representing NM_001165977.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAAGCACACACTGGCTCTGCTGGCTCCCCTGCTGGCCTGGGCTGGGCTGGCCTGAGTCAGCTG
GCTGCAGGGGCCACAGACTGCAAGTTCTTGGCCCGGCAGAGCACCTGACATTACCCCAGCAGCCAGG
GCCCGGTGGCTGGCCCTCGAGTTCGTGCGCCAGGACTCCTGGACTCCCTCTATGGCACCCTGCGCCGC
TTCTCTCGGTGGTGCAGCTCAATCCTTTCCCTTCAGAGTTGGTAAAGGCCCTACTGAATGAGCTGGCC
TCCGTGAAGGTGAATGAGGTGGTGCAGTACGAGGCGGGCTACGTGGTATGCGCTGTGATCGCGGGCCTC
TACCTGCTGCTGGTGCCTACTGCCGGGCTTTGCTTCTGCTGCTGCCGCTGCCACCGGCGCTGCGGGGA
CGAGTGAAGACAGAGCACAAAGGCGCTGGCCTGTGAGCGCGCGGCCCTCATGGTCTTCTGCTGCTGACC
ACCCTCTGCTGCTGATTGGTGTGGTCTGTGCCTTTGTCACCAACCAGCGCACGCATGAACAGATGGGC
CCCAGCATCGAGGCCATGCCTGAGACCTGCTCAGCCTCTGGGGCTGGTCTCTGATGTCCCAAGAG
CTGCAGGCCGTGGCACAGCAATTCTCCCTGCCCCAGGAGCAAGTCTCAGAGGAGCTGGATGGTGTGGT
GTGAGCATTGGGAGCGCATCCACACTCAGCTCAGGAGCTCCGTGTACCCCTTGCTGGCGCCGTGGGC
AGTTTGGGCCAGGTCCTGCAGGTCTCCGTGCACCACCTGCAAACCTTGAATGCTACAGTGGTAGAGCTG
CAGGCCGGGCAGCAGGACCTGGAGCCAGCCATCCGGGAACACCGGGACCGCCTCCTTGAGCTGCTGCAG
GAGGCCAGGTGCCAGGGAGATTGTGAGGGGCCCTGAGCTGGGCCCGCACCTGGAGCTGGGTGCTGAC
TTCAGCCAGGTGCCCTCTGTGGACCATGTCTGCACCAGCTAAAAGGTGTCCCGAGGCCAACTTCTCC
AGCATGGTCCAGGAGGAGAACAGCACCTTCAACGCCCTTCCAGCCCTGGTGGCATGCAGACATCCAGC
GTGGTGAAGAGCTGAAGAAGGCAGTGGCCAGCAGCCGGAAGGGGTGAGGACTGGCTGAAGAGTTC
CCGGCTTGGAGGCAGCTTCCCGCTGGGCCAGGCACTGCAGGAGGTGGAGGAGACAGCCGCCCTAC
CTGCAGGAGGTGCAGAGATACGAGACTACAGGTGGATCGTGGGCTGCGTGTGCTCCGTGGTCCCTA
TTCGTGGTGTCTGCAACCTGCTGGGCCTCAATCTGGGCATCTGGGCCTGTCTGCCAGGACGACCCC
AGCCACCCAGAAGCCAAGGGCGAGGCTGGAGCCCGCTTCTCATGGCAGGTGTGGGCCTCAGCTTCTCT
TTTGCTGCACCCCTCATCCTCCTGGTGTTCGCCACCTTCTGGTGGTGGCAACGTGCAGACGCTGGTG
TGCCAGAGCTGGGAGAACGGCGAGCTCTTTGAGTTTGCAGACACCCAGGGAACCTGCCCCCGTCCATG
AACCTGTCGCAACTTCTTGGCCTGAGGAAGAACATCAGCATCCACCAAGCCTATCAGCAGTGAAGGAA
GGGCAGCGCTCTGGACAGTCTGCAGCTCAACGACTCCTACGACTGGAGGAGCACCTGGATATCAAC
CAGTATACCAACAAGCTACGGCAGGAGTTGCAGAGCCTGAAAGTAGACACACAGAGCCTGGACCTGCTG
AGCTCAGCCGCCCGCCGGACCTGGAGGCCCTGCAGAGCAGTGGGCTTACGCGCATCCACTACCCCGAC
TTCTCTGTTAGATCCAGAGGCCGTGGTGAAGACCAGCATGGAGCAGCTGGCCCAGGAGCTGCAAGGA
CTGGCCCAGGCCAAGACAATTCTGTGCTGGGGCAGCGGCTGCAGGAGGAGGCCAAGGACTCAGAAAC
CTTACCAGGAGAAGGTGCTCCCCAGCAGAGCCTTGTGGCAAAGCTCAACCTCAGCGTCAGGGCCCTG
GAGTCTCTGCCCGAATCTCCAGCTGGAGACCTCAGATGTCTAGCCAAATGTCACCTACCTGAAAGGA
GAGCTGCCTGCCTGGGCAGCCAGGATCCTGAGGAATGTGAGTGTGTTTCTGGCCCGGAGATGGGC
TACTTCTCCAGTACGTGGCCTGGGTGAGAGAGGAGGTGACTCAGCGCATTGCCACCTGCCAGCCCTC
TCCGGAGCCCTGGACAACAGCCGTGTGATCCTGTGTGACATGATGGCTGACCCCTGGAATGCCTTCTGG
TTCTGCCTGGCATGGTGCACCTTCTTCTGATCCCCAGCATCATCTTGGCGTCAAGACCTCCAATAC
TTCCGTCTATCCGAAACGCCTCAGCTCCACCAGCTCTGAGGAGACTCAGCTCTTCCACATCCCCCGG
GTTACCTCCCTGAAGCTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001165977

Insert Size: 2505 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001165977.2](#)

RefSeq Size: 2935 bp

RefSeq ORF: 2505 bp

Locus ID: 150696

UniProt ID: [Q8N271](#)

Cytogenetics: 2q11.1

Protein Families: Transmembrane

MW: 91.9 kDa

Gene Summary: This gene encodes a member of the prominin family of pentaspan membrane glycoproteins. The encoded protein localizes to basal epithelial cells and may be involved in the organization of plasma membrane microdomains. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]
Transcript Variant: This variant (2) differs in the 3' UTR compared to variant 1. Variants 1, 2 and 3 encode the same protein.