

## Product datasheet for **SC327697**

### Neurabin 1 (PPP1R9A) (NM\_001166163) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Neurabin 1 (PPP1R9A) (NM_001166163) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPP1R9A
Synonyms:	Neurabin-I; NRB1; NRBI
Mammalian Cell Selection:	None
Vector:	<u><a href="#">pCMV6-XL5</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001166163, the custom clone sequence may differ by one or more nucleotides

```

ATGTTGAAAACCTGAGTCTTCAGGTGAACGAACCACTCTCAGAAGTGCCTCTCCTCACAGG
AATGCATATCGAACTGAGTTTCAGGCACTGAAAAGTACCTTTGACAAACCAAGTCAGAT
GGGGAACAAAAACAAAAGAAGGTGAGGGCTCCCAGCAGAGCAGGGGGAGGAAATATGGC
TCCAATGTCAACAGAATTAACCTATTTATGCAGATGGGTATGGAACCAACGAGAAT
GCTGCAGTCATTGCCAAAACAAGGGGAAAGGTGGACATTCATCTCCTCAGAGAAGAATG
AAGCCCAAGAATTTCTGGAAAAACAGATGGCTCAGTTGTTAAGTTGGAGTCTTCTGTT
TCTGAACGAATTAGTAGATTTGACACTATGTACGATGGCCCTTCATATTCCAAGTCACT
GAGACTCGAAAGATGTTTGAGAGAAGTGTGCATGAATCAGGACAGAACAACCGCTATTCC
CCAAAGAAAGAGAAAGCTGGAGGGAGTGAACCTCAGGATGAATGGGGAGGTTCCAAGTCC
AACAGAGGCAGTACTGATTCCTTGACAGCCTTAGCTCCCGAACTGAGGCTGTCTCCCA
ACTGTGAGTCAACTGAGTGCAGTATTTGAGAACACTGATTCTCCCACTGCCATCATTTCT
GAGAAGGCTGAAAACAATGAATACTCAGTGACTGGGCATTATCCCTTGAATTTACCATCT
GTTACTGTTACAAATCTTGACACATTTGGTCACCTGAAGGATTCTAATTCCTGGCCTCCT
TCAAACAAGCGAGGTGTTGATACAGAGGATGCTCACAAGAGTAATGCAACTCCAGTACCA
GAAGTGGCTTCTAAAAGTACCTCTCTAGCTTCGATACCTGGTGAAGAGATCCAGCAGAGC
AAGGAACCCGAGGACTCCACATCTAATCAACAGACTCCCGACAGCATTGACAAAGATGGT
CCTGAAGAACCTTGTGCTGAAAGTAAGGCAATGCCAAAGTCCGAAATCCCTTCACCACAA
AGCCAACCTGTTAGAAGATGCTGAAGCTAATTTGGTTGGAAGGGAGGCAGCAAAGCAACAG
AGGAAAGAACCTTGACAGGTGGTATTTACCTCTCCTGATGCTTCTGCATCCAGTTGTGGA
AAAGAAGTACCTGAAGATTCAAATAATTTGATGGTTCCCATGTGTACATGCACAGTGAC
TATAATGTGTATAGGGTGAGATCCAGGTATAATTCAGACTGGGGAGAGACAGGCACTGAG
CAGGATGAGGAGGAAGATAGTGATGAGAACAGTACTATCAGCCTGATATGGAGTACTCG
GAAATTGTTGGATTGCCAGAAGAAGAAGAAATCCAGCAATAGGAAAATTAAGTTTAGT
AGTGCTCCTATTAAGGTTTTCAACACATACTCCAATGAAGACTATGACAGGAGAAATGAC
GAAGTTGACCCTGTGGCTGCTTCAGCTGAGTATGAACTGAAAAACGTGTAGAAAAGCTG

```



[View online »](#)

GAACTTTTCCAGTGGAGCTAGAGAAAGATGAGGATGGTCTTGGTATAAGTATTATTGGA  
 ATGGGTGTTGGAGCAGATGCTGGACTTAAAAAGCTGGGAATATTCGTCAAGACAGTAACA  
 GAAGGTGGTGTCTCAACGGGATGGCAGAATACAAGTCAATGACCAGATTGTGGAAGTG  
 GATGGAATCAGCTTGGTGGTGTGACACAGAATTTGCAGCAACAGTTCTCAGAAACACC  
 AAGGGCAACGTCAGATTTGTTATTGGGCGGGAAAAACCAGGACAAGTGAGCGAGGTTGCC  
 CAGTTGATAAGCCAGACACTGGAACAGGAGAGGGCGCCAGAGAGCTGCTGGAACAGCAC  
 TATGCCAGTATGATGCCGACGATGACGAGACAGGAGAATATGCCACAGATGAAGAAGAA  
 GATGAGGTAGGACCTGTCTTCTGGCAGCGACATGGCCATTGAAGTCTTTGAGCTGCCT  
 GAGAATGAGGACATGTTTTCCCATCAGAACTGGACACAAGCAAGCTCAGTCACAAGTTC  
 AAAGAGTTGCAAATCAAACATGCAGTTACAGAAGCAGAGATTCAAAAAATTGAAGACCAAG  
 CTGCAGGCAGCAGAAAATGAGAAAGTGAGGTGGGAACTAGAAAAACCCAACTCCAACAA  
 AACATAGAAGAGAATAAGGAAAGAATGTTGAAGTTGAAAGCTACTGGATTGAGGCCAA  
 ACATTATGCCACACAGTGAATGAGCATCTCAAAGAGACTCAAAGCCAGTATCAGGCCCTG  
 GAAAAGAAATACAACAAGGCAAAGAAGTTGATCAAGATTTTCAACAAAAAGAGCTTGAT  
 TTCATCAAAAGACAGGAAGCAGAAAGAAAGAAAATAGAAGATTTGAAAAAGCTCATCTT  
 GTGGAAGTGCAAGGCCTCCAAGTGCGGATTAGAGATTTGGAAGCTGAGGTATTCAGGCTA  
 CTGAAGCAAAATGGGACTCAAGTTAACAATAATAACAACATCTTTGAGAGAAGAATCT  
 CTTGGTGAAGTCTCTAAAGGGGATACCATGGAGAACTTGGATGGCAAGCAGACATCTTGC  
 CAAGATGGCCTAAGTCAAGACTTGAATGAAGCAGTCCCAGAGACAGAGCGCCTGGATTCA  
 AAAGCACTGAAAACCGAGCCAGCTCTCTGTGAAGAACAGACGCCAGAGACCCTCTAGG  
 ACAAGACTGTATGATAGTGTAGTTCCACAGATGGGGAGGACAGTCTAGAGAGAAAAGAA  
 TTTACCTTCAATGATGACTTCAGTCCCAGCAGTACCAGTTCAGCAGACCTCAGCGGCTTA  
 GGAGCAGAACCTAAAACACCAGGGCTCTCTCAGTCTTAGCACTGTCATCAGATGAGATC  
 CTTGATGATGGACAGTCTCCCAAACACAGTCAGTGTGAGAATCGGCGCCTTCAGGAATGG  
 AGTGTGCAGCAGGTTTCTCACTGGTTAATGAGCCTAAATCTGGAGCAGTATGTATCTGAA  
 TTCAGTGCCTAAAACATCACTGGAGAACAGCTCCTGCAGTTGGATGGAAATAAACTTAAG  
 GCTCTTGGAAATGACAGCATCCCAGGACCGAGCAGTGGTCAAAAAGAACTCAAGGAAATG  
 AAGATGTCTCTAGAGAAGGCTCGGAAGGCCAAGAGAAAATGGAAAAACAAGAGAAAAG  
 CTAAGGAGAAAGGAGCAAGAGCAAATGCAGAGGAAGTCCAAAAAGACAGAAAAGATGACG  
 TCAACTACAGCCGAGGTGCTGGTGAGCAG

- Restriction Sites:** Please inquire
- ACCN:** NM\_001166163
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001166163.1](#), [NP\\_001159635.1](#)

**RefSeq Size:** 9601 bp

**RefSeq ORF:** 3273 bp

**Locus ID:** 55607

**UniProt ID:** [Q9ULJ8](#)

**Cytogenetics:** 7q21.3

**Protein Families:** Druggable Genome

**Gene Summary:** This gene is imprinted, and located in a cluster of imprinted genes on chromosome 7q12. This gene is transcribed in both neuronal and multiple embryonic tissues, and it is maternally expressed mainly in embryonic skeletal muscle tissues and biallelically expressed in other embryonic tissues. The protein encoded by this gene includes a PDZ domain and a sterile alpha motif (SAM). It is a regulatory subunit of protein phosphatase I, and controls actin cytoskeleton reorganization. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]

Transcript Variant: This variant (5) lacks an exon in the middle region and four consecutive exons in the 3' coding region, as compared to variant 1. The reading frame is not changed, but the resulting isoform (5) is shorter and lacks two internal segments, as compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.