

## Product datasheet for **SC327831**

### ZNF323 (ZSCAN31) (NM\_030899) Human Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	ZNF323 (ZSCAN31) (NM_030899) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZSCAN31
Synonyms:	ZNF20-Lp; ZNF310P; ZNF323
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC327831 representing NM_030899. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGGCTTCAACAGAGGAACAGTACGATCTTAAGATTGTGAAAGTGGAGGAAGACCCTATCTGGGACCAA
GAAACCCACCTTCGAGGGAACTTTTCTGGCCAAGAAGCCTCCCGACAACCTTTTAGGCAGTTTTGT
TACCAAGAGACTCCTGGTCCCGGAGAAGCTCTGAGCCGGCTCCGAGAATCTGTCATCAGTGGCTAAGG
CCAGAAATCCACACCAAGAGCAAATCTTGAGCTGCTGGTGTGGAGCAATTCCTGACTATCCTGCCT
GAGGAGCTCCAGGCTGGGTGCGGGAGCACCATCCGAGAGTGGGAGGAGGCTGTGGCTGTAGTTGAA
GATCTGGAACAAGAGCTTAGTGAGCCAGGAACCAGGCTCCAGACCATGAACATGGACATTCTGAAGTG
CTTTGGAGGATGTGGAACATCTGAAGGTCAAGCAGGAACCAACAGACATACAGCTTCAGCCTATGGTG
ACACAGCTCAGATATGAATCTTTTTGCCTCCACCAATTTCAAGAACAAGATGGTAAAGTATACCTGAG
AACCAGGAGTTGGCATCAAAGCAAGAAATCTTAAAAGAAATGGAACATTTGGGGATAGCAAACCTCAA
AGAGATGTATCTTTGGATTCTAAGTACAGAGAACTTGTAAACGAGACAGCAAGGCAGAAAAGCAGCAG
GCATTTCCACTGGAGAGAGACGCCACAGGTGCAATGAATGTGGGAAAAGCTTCACTAAGAGTTCAGTA
CTCATTGAGCACCAGAGAATCCACACTGGGAGAAGCCATATGAATGTGAAGAATGTGGGAAAGCCCTT
AGCCGGAGGTCAAGCCTGAATGAACATCGGCGGAGCCACACTGGAGAGAAAACCTATCAATGTAAGGAG
TGTGGGAAAGCCTTCAGTGCCAGCAATGGCCTCACTCGACAGAAAGAAATCCACAGGGGAAAACCA
TATGAATGCAAAGTGTGTGGGAAGGCTTTCCTCCTCAGCTCATGCCTTGTTCAGCATCAGAGGATACAC
ACTGGAGAGAAGCGCTATCAGTGTCTGAGTGTGGCAAAGCCTTCATTCAGAATGCAGGGCTTTTCCAG
CATCTCCGAGTCCACACTGGTGAGAAAACCTATCAGTGCAGTCACTGAGTAACTCTTTAGTAAGCGG
ACACTTCTTAAGAAACATCAGAAAATCCACACTGGAGAGAGACCATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_030899

**Insert Size:** 1221 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_030899.4](#)

**RefSeq Size:** 3039 bp

**RefSeq ORF:** 1221 bp

**Locus ID:** 64288

**UniProt ID:** [Q96LW9](#)

**Cytogenetics:** 6p22.3-p22.1

<b>Domains:</b>	LER, zf-C2H2
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	47.3 kDa
<b>Gene Summary:</b>	<p>This gene encodes a protein containing multiple C2H2-type zinc finger motifs. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016]</p> <p>Transcript Variant: This variant (1) encodes isoform 1. Variants 1, 2, 3, 4 and 7 encode the same isoform (1).</p>