

Product datasheet for **SC327838**

ATAD3C (NM_001039211) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ATAD3C (NM_001039211) Human Untagged Clone
Tag:	Tag Free
Symbol:	ATAD3C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC327838 representing NM_001039211. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCAAAGGACGCCCTGAATCTGGCGCAGATGCAGGAGCAGACGCTGCAGTTGGAGCAACAGTCCAAG
CTCAAACAACCTGTCAATGAGGATTTACGGAAGCAGGAGGAGTCCGTGCAGAAGCACCATCAGACCTTC
TTGGAGTCCATCAGGGCGGCTGGCACCTTGTGGGAAGGATCCGTGCCTTTGTGACAGACCGGGAC
AAAGTGACAGCCACGGTGGCTGGCTGACGCTGCTGGCTGTCGGGGTCTACTCAGCCAAGAATGCGACA
GCCGTCACTGGCCGCTACATCGAGGCTCGGCTGGGAAGCCGTCCTAGTGAGGGAGACGTCCCGCATC
ACGGTCTTGAGGCGCTGCGGCACCCATCCAGCAGGTCAGCGCGGCTCCTCAGTCGACCCAGGAC
GTGCTGGAGGGTGTGTGCTTAGTCCCAGCCTGGAAGCACGGGTGCGCGACATCGCCATAATGACAAGG
AACATCAAGAAGAACCAGGGCCTGTACAGGCACATCCTGCTGTACGGGCCACCAGGCACCGGAAGACG
CTGTTTGCCAAGAACTCGCCCTGCACTCAGGCATGGACTACGCCATCATGACAGGCGGGGACGTGGCC
CCCATGGGGCGGAAGCGTGACCGCCATGCACAAGCTCTTGACTGGGCAATACCAGCCGGCGCGGC
CTCCTGCTCTTTGTGGATGAAGCGGACGCCTTCTTCGGAAGCGAGCCACTGAGAAGATAAGCGAGGAC
CTCAGGGCCACACTGAACGCCTTCTGTACCGCACGGGCCAGCACAGCAAAATTCATGCTGATCCTG
GCCAGTGCACCCCGAGCAGTTCGACTGGCCATCAATGCCATCGACGTGATGGTCCACTTCGAC
CTGCCAGGGCAGGAGGAGCGGGCGCCTGGTGAGAATGTATCTTAACGAGTATGTTCTTAAGCCGGCC
ACAGAAGGAAAGCGGCTCTGAAGCTGGCCAGTTTGACTACGGGAGGAAGTGCTTAGAGATCGCTCGG
CTGACAGAGGGCATGTCATGCCGAAGATCGCACAGCTGGCCGTGCTTGGCAGGCCACGGCGTATGCC
TCCAAGGACGGGGTCTGACCGAGGCCATGATGGACGCCTGCGTGCAAGACTTTGTCCAGCAGCACCAG
CAGATGATGCGCTGGCTGAAGGGGGAGAGGCCTGGGCCGAGGACGAGCAACCCTCATCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001039211



[View online »](#)

Insert Size:	1236 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001039211.2</u>
RefSeq Size:	3859 bp
RefSeq ORF:	1236 bp
Locus ID:	219293
UniProt ID:	<u>Q5T2N8</u>
Cytogenetics:	1p36.33
MW:	46.4 kDa