

Product datasheet for SC327842

DUSP7 (NM_001947) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DUSP7 (NM_001947) Human Untagged Clone
Tag:	Tag Free
Symbol:	DUSP7
Synonyms:	MKPX; PYST2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC327842 representing NM_001947. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAAAACAGCTCCGCGGCCCGCCAGCGCGGGCGCACATGTCGACTTCGGGGCGGGCGGGCTGGG
GGCACCCGGGGGGTCCGAGCCGGTGC GGCGTCCGGCGCAGGCACCGGGCGGGCGCGGGC
ACGGGGCAGGGGCCATGCCCTGCAAGAGCGCCGAGTGGCTGCAGGAGGAGCTGGAGGCGCGGGC
GCGTCTTGCTGCTGCTGACTGCCGGCCGACGAGCTCTTCGAGTCGTCGCACATCGAGACGGCCATC
AACCTGGCCATCCCGGGCTCATGTTGCGCCGCTGCGCAAGGGCAACCTGCCCATCCGCTCCATCATC
CCCAACACGCGGACAAGGAGCGCTTCGCCACGCGCTGCAAGGGCGCCACCGTGTGCTCTACGACGAG
GCCACGGCCGAGTGGCAGCCGAGCCGGCGCTCCCGCTCCGTGCTCGGCTGCTCCTACAGAAGCTG
CGCGACGACGGCTGCCAGGCCTACTACCTCCAAGGTGGTTTCAACAAGTTTCAAACAGAGTACTCTGAG
CACTGCGAGACCAACGTGGACAGCTCTTCTCGCCGAGCAGCTCGCCACCCACCTCAGTGTGGGCTG
GGGGGCTGCGCATCAGCTCTGACTGCTCCGACGGCGAGTCGGACCGAGAGCTGCCACGAGTGGCCACC
GAGTCAGACGGCAGCCCTGTGCCATCCAGCAACAGCCTTCCCTGTCCAGATCCTGCCCTACCTTAC
CTCGGCTGCGCAAGGACTCCACCAACCTGGACGTGCTCGCAAGTATGGCATCAAGTATATCCTCAAT
GTCACACCAACCTACCAACGCCTTCGAGCACGGCGGAGTTCACCTACAAGCAGATCCCCATCTCT
GACCACTGGAGCCAGAACCTCTCCAGTCTTCCCTGAGGCCATCAGCTTATTGACGAAGCCCGCTCC
AAGAAGTGTGGTGTCTGGTGCCTGCTGGCAGGCATCAGCCGCTCAGTGACGGTCACTGTGGCCTAT
CTGATGCAGAAGATGAACCTGTCACTCAACGACGCTACGACTTTGTCAAGAGGAAAAAGTCCAACATC
TCGCCAACTTCAACTTCATGGGCAGCTGCTGGACTTTGAGCGGACGCTGGGGCTAAGCAGCCCGTGC
GACAACCACGCGTGCAGTGAGCAGCTCTACTTTCCACGCCACCAACCAACCTGTTCCCACTCAAT
ACGCTGGAGTCCACGTGA
ACGGGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAG
GATGACGACGATAAGGTTTAAACGGCCGGCCGGT
```



[View online »](#)

Restriction Sites: SgfI-NotI

Plasmid Map:



ACCN: NM_001947

Insert Size: 1260 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001947.3](#)

RefSeq Size: 3248 bp

RefSeq ORF: 1260 bp

Locus ID: 1849

UniProt ID: [Q16829](#)

Cytogenetics: 3p21.2

Protein Families: Druggable Genome, Phosphatase

Protein Pathways: MAPK signaling pathway

MW: 45 kDa

Gene Summary: Dual-specificity phosphatases (DUSPs) constitute a large heterogeneous subgroup of the type I cysteine-based protein-tyrosine phosphatase superfamily. DUSPs are characterized by their ability to dephosphorylate both tyrosine and serine/threonine residues. DUSP7 belongs to a class of DUSPs, designated MKPs, that dephosphorylate MAPK (mitogen-activated protein kinase) proteins ERK (see MIM 601795), JNK (see MIM 601158), and p38 (see MIM 600289) with specificity distinct from that of individual MKP proteins. MKPs contain a highly conserved C-terminal catalytic domain and an N-terminal Cdc25 (see MIM 116947)-like (CH2) domain. MAPK activation cascades mediate various physiologic processes, including cellular proliferation, apoptosis, differentiation, and stress responses (summary by Patterson et al., 2009 [PubMed 19228121]).[supplied by OMIM, Dec 2009]