

Product datasheet for **SC327851**

DNAJA4 (NM_018602) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DNAJA4 (NM_018602) Human Untagged Clone
Tag:	Tag Free
Symbol:	DNAJA4
Synonyms:	MST104; MSTP104; PRO1472
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC327851 representing NM_018602. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
 ATGGCCCGGGCGGCAGTCAGAGCTGGAGCTCCGGGGAATCAGACGGGCAGCCAAAGGAGCAGACGCC
 GAGAAGCCAGACACAAGATGGTGAAGGAGACCCAGTACTATGACATCCTGGCGTGAAGCCAGCGC
 TCCCGGAGGAGATCAAGAAGGCCTATCGGAAGCTGGCGCTCAAGTACCACCCGACAGAAGCCGGAT
 GAGGGCGAGAAGTTAACTCATATCCAGGCATATGAAGTGCTTCAGATCCAAAGAAAAGGGATGT
 TATGACCAAGGCGAGAGCAGGCAATTAAGAAGGAGGCTCAGGCAGCCCGAGCTTCTCTCACCCATG
 GACATCTTTGACATGTTCTTTGGTGGTGGACGGATGGCTAGAGAGAGAAGAGCAAGAATGTTGTA
 CACCAGTTATCTGTAACCTTGAAGATCTATATAATGGAGTCACGAAGAAATTGGCCCTCCAGAAAAT
 GTAATTTGTGAGAAATGTGAAGGTGTTGGTGGGAAGAAGGGATCGGTGGAGAAGTGCCCGCTGTGCAAG
 GGGCGGGGATGCAGATCCACATCCAGCAGATCGGGCCGGGCATGGTACAGCAGATCCAGACCGTGTGC
 ATCGAGTGCAAGGGCCAGGGTGAGCGCATCAACCCCAAGGACCGCTGCGAGAGCTGCAGCGGGGCCAAG
 GTGATCCGTGAGAAGAAGATTATCGAGGTACATGTTGAAAAGGTATGAAAGATGGGCAAAAGATACTA
 TTTCATGGAGAAGGAGATCAGGAGCCTGAGCTGGAGCCTGGTGTGTCATAATTGTGCTTGATCAGAAG
 GATCATAGTGTCTTTCAGAGACGAGGCCATGACTTGATCATGAAAATGAAAATTCAGCTTTCTGAAGCT
 CTTTGTGGCTTCAAGAAGACGATAAAAACATTGGACAATCGAATTCCTGTTATTACATCCAAAGCAGGT
 GAGGTGATAAAGCACGGGACCTGAGATGCGTGCGGATGAAGGAATGCCATCTACAAAGCACCCCTG
 GAAAAAGGGATTCTGATCATACAGTTTTTGTAAATCTTTCCTGAAAAACACTGGCTTTCTCTGAAAAAG
 CTTCTCAGCTGGAAGCTTTACTCCCTCCTCGACAGAAAGTGAGGATTACAGATGACATGGATCAGGTG
 GAGCTGAAGGAGTTTGTCCCAATGAGCAGAACTGGCGTCAGCACAGGGAGGCCTACGAGGAGGACGAA
 GACGGGCCCCAGGCTGGAGTGAGTGCCAGACGGCATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC


[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_018602
Insert Size:	1281 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_018602.3</u>
RefSeq Size:	3203 bp
RefSeq ORF:	1281 bp
Locus ID:	55466
UniProt ID:	<u>Q8WW22</u>
Cytogenetics:	15q25.1
Domains:	DnaJ_CXXCXGXG
MW:	48 kDa