

Product datasheet for **SC327864**

SPRYD5 (TRIM51) (NM_032681) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPRYD5 (TRIM51) (NM_032681) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPRYD5
Synonyms:	SPRYD5; TRIM51A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC327864 representing NM_032681. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
ATGAATTCTGGAATCTTGCAAGTCTTCCAGAGGGCACTCACCTGTCCCATCTGCATGAACACTTCTCTCA
GACCCAGTCACCATAGACTGTGGGCACAGCTTTTGCCGGCCCTGTTGTACCTCACTGGCAAGACACG
GCAGTTCTTGCTCAGTGTCTGAATGCAAGAAGACAACGCGGCAGAGAAACCTCAACACTGACATTTGT
TTGAAGAACATGGCTTTTCATTGCCAGAAAAGCCAGCCTCCGGCAATTCCTTAGCTCTGAGGAGCAAATA
TGTGGGATGCACAGAGACAAAGAAGATGTTCTGTGAAGTGGACAAGAGCCTGCTCTGTTTGCCGTGC
TCCAACCTCTCAGGAGCACCGGAATCACATACACTGTCCCATTGAGTGGGCTGCTGAGGAACGCCGGGAG
GAGCTCCTAAAAAAATGCAGTCTTTATGGGAAAAGCTTGTGAAAATCTCAGAAATCTGAACATGGAA
ACCACAAGAACCAGATGCTGGAAGGATTATGTGAGTTTAAAGATAGAAGCAATCAGAGCTGAATATCAG
AAGATGCCTGCATTTCTCCATGAAGAAGAGCAACATCACTTGGAAAGGCTGCGAAAGGAGGGCGAGGAC
ATTTTTCAGCAACTCAATGAAAGCAAAGCCAGAATGGAACATTCAGGGAGCTTTTAAAGAGGAATGTAT
GAGGATCTGAAGCAAATGTGCCATAAAGCAGATGTGGAGCTACTCCAGGCTTTTGAGACATATTACAC
AGGTATGAGTCTCTGCTGCTGCAAGTGTCTGAGCCTGTGAATCCAGAGCTCAGTGCAGGGCCCATCACT
GGACTGCTGGACAGCCTCAGTGGATTGAGAGTTGATTTTACTCTGCAGCCTGAAAGAGCCAATAGTCAT
ATCTTCTGTGTGAGATTTGAGAAGCATGAATGTTGGATGTGACCTCAAGATGATCCCGATATCACT
GGAAAATCTGAATGTTTTCTTGTATGGGGGCTCAGGCTTTCACATCTGGCAAATATTATTGGGAGGTT
CATATGGGGGACTCTTGAATTGGGCTTTTGGTGTCTGTAACAATTATTGGAAAGAGAAGAGACAGAAT
GACAAGATAGATGGAGAGGAGGACTCTTTCTTCTGGATGTGTTAAGGAGGACACTCACTGCAGTCTC
TTTACCACCTCCCACCTGTGGTGCAATATGTTCCAAGACCTACCAGCACAGTAGGATTATCTCTGGAT
TGTGAAGGTAGAACCGTGAGCTTTGTTGATGTTGATCAAAGTTCCTGATATACCATCCCAATTGC
TCCTTCTACCTCCTCTCAGGCCTATCTTTGCTGTAGTCACTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_032681
Insert Size:	1359 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_032681.3
RefSeq Size:	1648 bp
RefSeq ORF:	1359 bp
Locus ID:	84767
UniProt ID:	Q9BSJ1
Cytogenetics:	11q12.1
Protein Families:	Druggable Genome
MW:	52.3 kDa