

Product datasheet for **SC327874**

ZNF577 (NM_032679) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF577 (NM_032679) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF577
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327874 representing NM_032679.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAAAAATGCCACGATTGTAATGTCTGTGAGGAGAGAGCAAGGCAGTTCTTCAGGGGAGGGTCATTG
TCATTCGAAGATGTGGCTGTGGCTTCACCAGGGAGGAGTGGCAGTTTTTGGACCACTCTCAGAAGGTC
TTGTACAAGGAAGTAATGTTGGAGAACTACATCAACCTAGTATCAATAGGGTATCGAGGCACCAAGCCA
GATTTCGCTCTTCAAGTTGGAGCAAGGAGAACCCCAAGGATAGCAGAAGGAGCAGCCACAGTCAATC
TGTCCAGGTTTTGTAATCCAGAGTAGAAGATATGCAGGAAAAGATTCTGATGCATTTGGTGGATATGGG
AGATCATGCCTCCATATCAAGCGTGACAAAACCTTACTGGAGTTAAATACCATAGATGTGTTAAACCC
AGCAGCCCTAAATCACAGCTCAATGACCTACAAAAATTTGTGCAGGAGGAAACCATGAATGCAGT
GTGTGCGGGAGAGCCTTCTCCAGGAAAGCACAGCTTATTCAACATCAGAGAACTGAAAGAGGAGAGAAA
CCCCATGGATGTGGTGAATGTGGGAAAACATTCATGAGGAAGATTCAGCTCACTGAGCATCAGAGAACT
CACACAGGAGAGAAGCCCATGAATGTAGTGAATGTGGAAAAGCCTTCTCCAGAAAGTCACAGCTCATG
GTCCATCAGAGAACCCATACAGGAGAGAAAACCTACAGATGCAGCAAATGCGGAAAAGCCTTCAGCCGG
AAGTGCCGGCTCAATAGACATCAGCGATCACATACTGGAGAGAAAACCTCTATGGGTGCAGTGTGTGGG
AAAGCCTTTTCTCAGAAGGCATACCTCACTGCACACCAGAGACTTCACACAGGAGATAAGCCTTATAAA
TGCAGTGATTGTGGAAGAACCTTCTATTTTAAGTCAGACCTGACCAGACATCAGAGGATTCATACGGGA
GAGAAACCTTATGAATGTAGTGAGTGTGAAAAAGCCTTTAGAAGCAAGTCGAAGCTCATTGAGCATCAG
CGTACTCACACTGGAGAGAGACCATATTCATGCAGAGAATGTGGCAAAGCCTTTGCCACATGTCAGTC
CTCATTAACATGAGAAAACTCACATAAGAGAGACAGCCATAAATTCAGTGACGGTGGAGAAACCTTCC
TCAAGGAGTCATACCTCCTTATACATGAGCGAACTCATACAAGAGCAAAAGACTGTGAACACAGTACCT
ATAGAAATGCCTTCCTCAGGAACCCCGCCATTGTTAAACAAGAGTGAGCGCCTAGTGGGCAGAAATGTA
GTGATTGTGGAACAACCTTTTCCAAGAAATCAAGCCTTTGTAGTTAATCAGGAATTTGAACAGAGAATA
AGCCTCACAAATGAAGTGAATGTGGCCCATCAGTAATAAATTATATCTTGTATCTTACAGATATTGTA
TCAGAAATA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_032679

Insert Size: 1458 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_032679.2</u>
RefSeq Size:	3096 bp
RefSeq ORF:	1458 bp
Locus ID:	84765
UniProt ID:	<u>Q9BSK1</u>
Cytogenetics:	19q13.41
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors
MW:	54.9 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (a).