

Product datasheet for **SC327882**

ZSCAN5B (NM_001080456) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZSCAN5B (NM_001080456) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZSCAN5B
Synonyms:	ZNF371; ZNF495B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327882 representing NM_001080456.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTGCAAATTGGACACTCTCATGGGTGAGGAGGACCCCTGCAACAGCCCTGGGTGAGACACTCCA
CGGTCTGTGGCGTCCCCAGAACTCAACTTGAAATCACGACAGGAACCCCTGAGACTTGGCACATGAAT
TTCAGGATGTTTCAGCTGCCAGAGGAGTCGGACCCATCCAGGCTCTGAGGAACTCACTGAGCTGTGC
CATCTGTGGCTGAGGCCGACCTCCACACCAAGAGCAGATCCTGGACATGCTGGTGATGGAGCAGTTC
ATGATCTCCATGCCCCAGGAGCTCCAGGTCTTAGTCAAGGTGAACGGTGTGCAGAGCTGCAAGACCTG
GAGGACCTGCTACGAAATACAGAAGACCCAAGAAATGGTCTATAGTCAACTTGTCTGGCAAGAATAT
CTTATGCTGAACCTCAGATGTCGAGATGGCTGAAGCCCCGCCAGTGTGAGAGATGATCCGAGAGACGTG
TCCAGCCAGTGGGCTCCTCTGTGAACCAGATGCATCCGGGACAGGCCAGGCCGCCGAGAGCAGCAG
ATCCTGCCAGGTGCGTGCCTGTGACCTGTCCAGGAGGAGGAGGAGGACTTTCTGCTACACAAGAGTATTGAC
GTAAACAGGTGACCAAACTCTCCGAGACCAAGCAGACCTTGAGAAGGATCTGAAGGAAACAGGGAA
GAGAACCAGGACTGTATCCCAAGAGCTCAGCTTCAAAGAGTCCCAATCTGGTGAGAGCAAAGGAG
GGGAAGGAACCCAGAAAAGAGCCTCTGTGAAAATGTGGATGCTGACACACCTTCTGCCTGCGTTGTG
GAGAGAGAAGCTTTGACTCACAGCGGGAACAGAGGAGATGCTCTGAATCTGAGCAGTCCCAAGAGC
AAACCAGATGCCTCCTCCATTTCCCAAGAGAGCTCAAGGAGAAGCCACACCTGTGGGCAACAGAGAA
TCCCCAGGACAAGCTGAGATCAATCCAGTTCATTTCCAGGCCCTGCGGGCCCGGTGAGTACCCAGAT
GGCCAAGAAGCAAGGCCCTGCCGCCCTTTCATGTGACGTGTGCAATAAATCATTTAAGTATTTTCC
CAGCTAAGCATCCACAGGAGTCCACACAGGAGACAGACCTTTCAATGTGATCTCTGTGGAAGCGC
TTCTTGACGCCCTCAGACCTCCGAGTTCACGAGGAGTCCACACTGGCGAGAGGCCCTACATGTGTGAC
GTCTGCCAAAAGCGGTTTCGCCACGAGTCCACCTTACAGGCCCAAGAGGATCCACACCGGGGAGAGG
CCGTTTAAATGTAATACTGCAGCAAAGTTTTTCAGCCACAAGGGGAACCTGAACGTTACACAGCGCACC
CACTCCGAGAGAAGCCCTACAAATGTCCACCTGTCAAAGGCCCTCCGTGAGTGGGACATTCAAG
CGTCACCTGAAAACACACCGGGAACCACTCCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001080456

Insert Size: 1488 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080456.2</u>
RefSeq Size:	1626 bp
RefSeq ORF:	1488 bp
Locus ID:	342933
UniProt ID:	<u>A6NJL1</u>
Cytogenetics:	19q13.43
MW:	55.8 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]