

Product datasheet for **SC327889**

THUMPD2 (NM_025264) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	THUMPD2 (NM_025264) Human Untagged Clone
Tag:	Tag Free
Symbol:	THUMPD2
Synonyms:	C2orf8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327889 representing NM_025264.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCGGAGGCGCGTGGAGAGCCAGGGTCCGGGCTGAGGCTGGCGCCCGATTCTTCTGCACTGCGGGT
CGCGGCTGGAGCCGTTTCGTAATGCGAGAGGTGCGGGCGCGGCTGGCGGCCACGCAGTTGAATATATT
TCAGGAAAGGTTTTTTTACCACCTGTTCTGATTTGAATATGTTGAAGAAATTTAAATCTGCAGAAAGA
TTATTTTGTGATTAATAAGCAGTTTCCACTTATTATTTCTTCTGTAAGTAAAGGAAAAATATTTAAT
GAAATGCAAAGACTTATAAATGAAGATCCAGGAAGTTGGTTGAATGCCATTTCAATTTGGAAAAATCTT
CTTGAACCTTGATGCAAAAAAGGAAAACTTTCTCAGAGAGATGATAACCAACTAAAAAGAAAAGTGGGA
GAAATGAAATCATTGCAAAGAAATTTAAATAGAACAAATGCAAAGATAGAAGAGAATAGGGACTGC
CAGCTGAAAAACAAATAAAGAAGAACTCTGGAGCAAAGAGATTTTACCCTAAAAGCGAAAAGTTT
CAAGAAGAAGAATTTAGAATGACATAGAGAAAGCAATTGATACTCATAATCAGAATGACTTGACTTTC
AGAGTATCTTGTGCTGCAGTGGAACTATTGGAAAGGCCCTCACTGCACAGGAGGTAGGAAAAGTAATT
GGAATTGCTATTATGAAACACTTTGGATGGAAAGCAGACTTGAGGAATCCACAATTAGAGATCTTTATA
CATCTAAATGACATTTACTCTGTGGTGGGATTCTGTGTTTCAGGGTTTCCCTAGCCAGCAGAGCTTAC
ATCAAGACAGCTGGACTGCGATCTACAATAGCGTGGGCAATGGCATCTCTGGCTGACATTAAGGCTGGT
GCATTTGTTTTAGATCCAATGTGTGGACTTGGAACTACTTTTGGAAAGCTGCTAAAGAATGGCCAGAT
GTGTATTATGTAGGTGCTGATGTCAGCGACTCACAGTTACTAGGTTACTTTGGACAATCTGAAAGCTGCA
GGCCTTGAGGATAAAATGAATTAATCTCTGTTATAGAATTGCCATTGCCTTCAGAAAGTGT
GATATTATTTCTGACATTCATTTGGGAAAAAGTTAAGTTAGGAAAAGACATCAAAAGCATTCTA
CAAGAAATGGAAAGAGTGCTTCATGTTGGCGGAACCATTTGATTGTTGCTTAGTGAAGATCACCACAGG
CGCCTTACAGATTGTAAGAGAGCAACATCCCTTTCAATTCCAAGGACAGTCACACAGATGAACCTGGA
ATTAATAAGTGCCTGAATCCTGAAGAAAAACTGGTGCATTCAAGACAGCGTCAACTTCATTGAAAGCC
AGTAACCACAAATTTAGACAGAATGTCACCATTTGGCTCCTTGGTACCAGTGGAAATGCTACAAAGTT
AGCCTTGGAAAGACAGATGCGTTCATATGTAATATAAGAAGTCGCACTCTTCTGGACTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_025264

Insert Size: 1512 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_025264.4](#)

RefSeq Size: 2224 bp

RefSeq ORF: 1512 bp

Locus ID: 80745

UniProt ID: [Q9BTF0](#)

Cytogenetics: 2p22-p21

Protein Families: Druggable Genome

MW: 56.4 kDa