

Product datasheet for **SC327900**

ZNF776 (NM_173632) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF776 (NM_173632) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF776
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_173632, the custom clone sequence may differ by one or more nucleotides

```

ATGGCGGCGGCCGCTGAGGCCCGGCTCAGGGCACTGTGACCTTTGAAGATGTGGCT
GTGAACTTTTCCAGGAGGAATGGAGTCTCCTTAGTGAGGCTCAGAGATGCCTTTATCAT
GACGTGATGCTGGAGAACCTGACACTTATATCTTCTCTAGGTTGTTGGTATGGAGCAAAA
GACGAGACACCTTCTAAGCAGACCCTTTCTATACAACAGGAGTCCCCACTCAGGACACAT
TGGACAGGTGTATGTACCAAGAAGGTCCACCTCTGGGGAATGTGTGGCCCTCTCTGGGA
GATATCTTACACCAGGGAACACAACACAATCAGAAATTGAATGGGTTTGGGGCATATGAA
AAAAAATTGGATGACGATGCAAAACCATCATCAAGACCAGAAGCAGCACATTGGAGAGAAA
TCGTACAGAAGCAATGCCAAGGGAACATCTTTTGTAAGAAGTGTAAATCCATATGTCA
CATGAGCCATTTATCTTTTATGAGGTTGGGAAAGACTTTTTGTCCAGCTTGAGATTACT
CAACAAGAGGACATTCACACTTCAGGGAAGTCAAACCTTTGAAACTAAGCATGGGATACCC
CTTCAGGGTGGAAAACTCATTACATCTGTGGAGAGTCCACAATACCGTTTAGCAACAAA
CACTCACTTGTCTTCACCAGAGACTTCTCCCTAGAGAAGGACCTTATGTATGCAGTGAT
TCTGGGAAATTCAGTAGCAAAAGTAATAGTTTTAATAATCATCAGGGAGTTCGCACTGGA
AAAAGACCTTATCAGTGTGGACAATGTGATGAATCATTTTGGTATAAGGCCACCTCACT
GAACACCAGAGAGTTCACACTGGAGAAAGACCTTATGAGTGTGGAGAATGTGATAAATCT
TTTAGTCATAAGCACAGTCTTGTGACCATCAGCGAGTTCACACTGGAGAAAGACCTTAT
GAATGTGACGAATGTGGGAAATCTTTAGCCATAAGCGCAGCCTTGTTTACCACCAGCGA
GTTTCACTGGAGAAAGACCTTATCAGTGTGGAGAATGTGGGAAATCGTTTAATCACAAG
TGCAACCTCATTACAGCATCAGCGAGTTCACACTGGAGAAAGACCTTTGAGTGTACGGCA
TGTGGGAAGTTATTTAGGAGCAACTCCACCTAAAGGAACACCAGAGAGTTCACACTGGA
GAAAGACCCTATGAGTGTAAAGAATGTAGGAAATCATTTAGGTACAAGTACACCTCACT
GAACACCAGAGAGTTCACACTGGAGAAAGGCCATATGAGTGTAGAGAATGTGGGAAATGT
TTTCATCAAAAGGCGAGTCTCATTCAACATCAGCAGATCCACTCTGGAGAAAGGCCACAT
GAGTGTGGAGAATGTGGGAAATGTTTTATCAAAAGGCGAGTCTCATTGACATCAGCAG
ATTCACCTCTGGAGAAAGGCCACATGAGTGTGGAGAATGTGGAAAGTGTTCGTCAAAAG
GGAAACCTCATTAAACATCAACGAGTTCACACGGGAGAAAGACATCATGAATGT
  
```


[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_173632
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_173632.3</u> , <u>NP_775903.3</u>
RefSeq Size:	5292 bp
RefSeq ORF:	1557 bp
Locus ID:	284309
UniProt ID:	<u>Q68DI1</u>
Cytogenetics:	19q13.43
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longer isoform (1).