

## Product datasheet for **SC327902**

### **GADL1 (NM\_207359) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	GADL1 (NM_207359) Human Untagged Clone
Tag:	Tag Free
Symbol:	GADL1
Synonyms:	ADC; CSADC; HuADC; HuCSADC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC327902 representing NM\_207359.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGCAGCGACTCGGACCGCCAGTGCCTGTGGACGGAGATTTGATCAACAAGAGATGATTCCAAGT
AAGAAGAATGCTGTTCTGTGGATGGGGTTGTGCTGAATGGTCCACAACAGATGCAAAAGCTGGAGAA
AAATTTGTTGAAGAGCCCTGTAGGCTAATAATGGAAGAGGTGGTTTTGAAAGCTACAGATGCAATGAG
AAGGTGTGTGAATGGAGGCCCTCTGAACAACGAAACAGCTTCTTGATTTGGAGATGAGAGACTCAGGC
GAGCCACCCATAAACTATTGGAACCTGTCTGGGATGTCATACACTACAGTGTCAAACTAACCACCCA
AGATTTTTCAACCAATTGTATGCTGGACTTGATTACTCCTTGGTGGCCCGATTTATGACCGAAGCA
TTGAATCCAAGTGTTTATACGTATGAGGTGTCCCAGTGTCTGTTAGTGAAGAAGCGGTTCTGAAG
AAAATGATTGAATTTATTGGCTGAAAAGAAGGGGATGGAATTTAACCCAGGTGGCTCAGTGTCCAAT
ATGTATGCAATGAATTTAGCTAGATACAAATATTGTCCTGATTAAGGAAAAGGGGCTGTCTGGTTCG
CCAAGATTAATCCTTTTACATCTGCAGAGTGCATTACTCTATGAAGAAGGCAGCCTCTTTTCTGGG
ATTGGCACTGAGAATGTTTGCTTTGTGGAAACAGATGGAAGAGGTAAAATGATACCTGAGGAACTGGAG
AAGCAAGTCTGGCAAGCCAGAAAAGAGGGGGCAGCACCGTTTCTTGTCTGTGCCACTTCTGGTACAAC
GTGTTGGGAGCTTTTGACCCTCTGGATGAAATAGCAGACATCTGCGAGAGGCACAGCCTCTGGCTTCAT
GTAGATGCTTCTTGGGTGGCTCAGCTTGTATGTCGAGGAAGCACCGCAAGCTTCTGCATGGCATCCAC
AGGGCTGACTCTGTGGCCTGGAACCCACACAAGATGCTGATGGCTGGGATCCAGTGTCTGTCTCCTT
GTGAAAGACAAATCTGATCTTCTAAAAATGCTACTCTGCCAAGGCATTTACCTTCCAGCAGGAT
AAATCTATGATGTGAGCTATGACACAGGAGACAAGTCTATCCAGTGTAGCAGAAGACCAGATGCATTC
AAGTTCTGGATGACCTGGAAGGCCCTGGGTACATTAGGCCTTGAAGAAAGAGTTAATCGTCTCTGCT
TTATCTAGGTACCTAGTAGATGAAATCAAGAAAAGAGAAGGATTCAAGTTACTGATGGAACCTGAATAT
GCCAATATTTGCTTTTGGTACATTCCACCAGCCTCAGAGAGATGGAAGAAGACCAGGTTCTGGGCA
AACTTAATTTGGTGGCCCGCCAGCCATTAAGGAGAGGATGATGAAGAAGGGAAAGCTTGATGCTGGGCTAC
CAGCCGCACCGGGAAAGGTCAACTTCTCCGCCAGGTGGTATCAGCCCTCAAGTGAGCCGGGAGGAC
ATGGACTTCTCCTGGATGAGATAGACTTACTGGGTAAGACATGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** Sgfl-Mlul

**ACCN:** NM\_207359

**Insert Size:** 1566 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_207359.2</a>
<b>RefSeq Size:</b>	3655 bp
<b>RefSeq ORF:</b>	1566 bp
<b>Locus ID:</b>	339896
<b>UniProt ID:</b>	<a href="#">Q6ZQY3</a>
<b>Cytogenetics:</b>	3p24.1-p23
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	59.2 kDa
<b>Gene Summary:</b>	May catalyze the decarboxylation of L-aspartate, 3-sulfino-L-alanine (cysteine sulfinic acid), and L-cysteate to beta-alanine, hypotaurine and taurine, respectively. Does not exhibit any decarboxylation activity toward glutamate.[UniProtKB/Swiss-Prot Function]