

Product datasheet for **SC327903**

ACBD5 (NM_145698) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ACBD5 (NM_145698) Human Untagged Clone
Tag:	Tag Free
Symbol:	ACBD5
Synonyms:	RDLKD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC327903 representing NM_145698.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGCTCTTCTCTCGTTTCATGCAGGCTCTTGGGAAAGCTGGTGTGCTGCTGCCTGATTCCCGCCGAC
AGACCTTGGGACCGGGCCAACACTGGCAGCTGGAGATGGCGGACACGAGATCCGTGCACGAGACTAGG
TTTGAGGCGGCCGTGAAGGTGATCCAGAGTTTGC CGAAGAATGGTTCAATCCAGCCAACAAATGAAATG
ATGCTTAAATTTTATAGCTTCTATAAGCAGGCAACTGAAGGACCCGTAAACTTTCAAGGCTGGATTT
TGGGATCCTATTGGAAGATATAAATGGGATGCTTGGAGTTCACCTGGGTGATATGACCAAAGAGGAAGCC
ATGATTGCATATGTTGAAGAAATGAAAAAGATTATTGAACTATGCCAATGACTGAGAAAGTTGAAGAA
TTGCTGCGTGTACATAGGTCATTTTATGAAATGTGCGAGGACAAAAGAGTGGCAGGAGTTCTGATATA
ACCTCAGATCTTGGTAAATGTTCTCACTTCTACTCCAACGCCAAAACCGTTAATGGTAAAGCTGAAAGC
AGTGACAGTGGAGCCGAGTCTGAGGAAGAAGAGGCCCAAGAAGAAGTAAAGGAGCAGAACAAGTGAT
AATGATAAGAAAATGATGAAGAAGTCAGCAGACCATAAGAATTTGGAAGTCATTGCTCACTAATGGCTAT
GATAAAGATGGCTTTGTTCCAGGATATACAGAATGACATTCATGCCAGTCTTCCCTGAATGGCAGAAGC
ACTGAAGAAGTAAAGCCATTGATGAAAACCTGGGGCAAACCTGGAAAATCTGCTGTTTGCATTCACCAA
GATATAAATGATGATCATGTTGAAGATGTTACAGGAATTCAGCATTTGACAAGCGATTGACAGAGTGAA
GTTTACTGTGATTCTATGGAACAATTTGGACAAGAAGAGTCTTTAGACAGCTTTACGTCCAACAATGGA
CCATTTCAAGTATTACTTGGGTGGTCATTCCAGTCAACCCATGGAAAATCTGGATTTCTGTAAGATATT
CAAGTACCTCTGGAAAATGGCAACATTTGGGAATATGCAGGTGGTTCAGTTGAAGAAAAGGTGAAGTC
AAGCATGGAGGAGAAGATGGCAGGAATAACAGCGGAGCACCACACCGGGAGAAGCGAGGCCGAGAAAAT
GACGAATTCTCTAATGTTAGAAGAGGAAGAGACATAGGATGCAACACTTGAGCGAAGGAACCAAGGCC
CGGCAGGTGGGAAGTGGAGGTGATGGGGAGCGCTGGGGCTCCGACAGAGGGTCCCGAGGCAGCCTCAAT
GAGCAGATCGCCCTCGTCTGATGAGACTGCAGGAGGACATGCAGAATGTCTTCAGAGACTGCAGAAA
CTGAAAACGCTGACTGCTTTCAGGCAAATCATCAACATCAACATTGCAGACTGCTCCTCAGCCCACC
TCACAGAGACCATCTTGGTGGCCCTTCGAGATGTCTCTGGTGTGCTAACGTTTGCCATCATATGGCCT
TTTATTGCACAGTGGTTGGTGTATTTATACTATCAAAGAAGGAGAAGAAAATGAACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_145698

Insert Size:	1578 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_145698.4</u>
RefSeq Size:	3996 bp
RefSeq ORF:	1578 bp
Locus ID:	91452
UniProt ID:	<u>Q5T8D3</u>
Cytogenetics:	10p12.1
Domains:	ACBP
Protein Families:	Transmembrane
MW:	59 kDa
Gene Summary:	<p>This gene encodes a member of the acyl-Coenzyme A binding protein family, known to function in the transport and distribution of long chain acyl-Coenzyme A in cells. This gene may play a role in the differentiation of megakaryocytes and formation of platelets. A related protein in yeast is involved in autophagy of peroxisomes. A mutation in this gene has been associated with autosomal dominant thrombocytopenia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]</p> <p>Transcript Variant: This variant (1) encodes isoform 1.</p>