

Product datasheet for **SC327907**

ZNF480 (NM_144684) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF480 (NM_144684) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF480
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327907 representing NM_144684.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGTGTGATGAAAAAGCCCAGAAGAGAAGGAAGAAAGCAAAAGGAGTCAGGGATGGCTCTTCTC
CAGGGACACTTAACATTCAGGGACGTGGCCATAGAATTCTCTCAGCGGAGTGGAATGCCTGGACCCT
GCACAGAGGGCTTTATACAAGGATGTGATGTTGGAGAACTACAGGAACCTGGTCTCCCTGGGAATCTCT
CTTCTGACCTGAATATTAACCTCATGTTGAGACAAAGGAGGGAGCCCTGGTCTGGTGAAGTGAAGTG
AAAATAGCAAAAAATTCAGATGGGAGGGAGTGCATCAAAGGTGTGAACACAGGGAGCAGCTATGCATTG
GGAAGCAATGCAGAAGACAAACCAATTAACAAACAACCTGGAGTATCCTTTCACCTACATCTGTCTGAA
CTGGAGCTATTTCCAGATGAAAGGGTAATAAATGGATGTAATCAAGTTGAAAACCTTATCAACCACAGT
TCCTCTGTTTCTGTCTTCAAGAAATGTCTTCCAGTGTCAAACCCCATTTTTAATAGGAATGATTTT
GATGATTTCTTCAATTTCTCCACAAGAACAGAAAGTACACCTTAGAGAAAAACCTTATGAATGTAATGAG
CATAGCAAAGTCTTTAGAGTATCTTCCAGCCTTACTAAACATCAAGTAATCCATACTGTAGAGAAACCT
TACAAATGTAATTCATGCGGCAAGGTCTTTAGTCGCAATTCACACCTTGCAGAACATTGTAGAATTCAT
ACTGGAGAGAAACCTTACAAATGTAATGTCTGTGGCAAGGTTTTAGTTACAATTCAAAACCTTGCACGA
CATCAAAGAAATTCATACCAGAGAGAAGCCGTATGAATGTAATGAATGTGGTAAAGTCTTCAGTAATAAT
TCTTACCTTGCACGACATCAAAGAAATTCATGCTGAAGAGAAACCTTACAAATGTAATGAATGTGGTAAA
GGCTTCAGTCATAAGTCATCTCTAGCAAATCATTGGAGAATTTATACTGGAGAGAAGCCTTACAAATGT
GATGAATGTGGCAAGGCCTTCTATAGGATTGCGCTCCTTGTACGACATCAGAAAAATTCATACTGGAGAG
AAACCTTACAAATGTAATGAATGTGAAAGGTCTTTATTCAAAATTCGCACCTAGCACAAACATTGGAGA
ATTCATACAGGAGAGAAACCTTACAAATGTAATGAATGTGAAAGTATTTAATCAACTTCAAATCTT
GCACGACATCGAAGAATTCATACTGGAGAGAAGCCTTACAAATGTAATGAATGTGGTAAAGCATTAGT
GAGTATTCAGGCCTTTCAGCCATCTTGTAAATCCACACTGGAGAGAAGCCTTACAAATGTAATGTAATGT
GGCAAGGCATTACAGACACAAGTTATCACTAACCAATCATCAGAGAATCCATACTGGAGAAAGACCTTAC
AAATGTAATGAATGTGGCAAGGTCTTCAATCGAATTGCACACCTTGCACGACATCGGAAAATTCATACT
GGAGAGAAACCTTACAAATGTAATGAGTGTGGCAAGGCCTTTCAGTCATTTACATACCTAGCACAAACAT
TGGACAATTCATATGGGATAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:	NM_144684
Insert Size:	1608 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_144684.3
RefSeq Size:	4767 bp
RefSeq ORF:	1608 bp
Locus ID:	147657
UniProt ID:	Q8WV37
Cytogenetics:	19q13.41
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors
MW:	61.7 kDa
Gene Summary:	Involved in transcriptional regulation as an activator.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longest protein. Sequence Note:.