

Product datasheet for **SC327910**

NAPRT1 (NAPRT) (NM_145201) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NAPRT1 (NAPRT) (NM_145201) Human Untagged Clone
Tag:	Tag Free
Symbol:	NAPRT
Synonyms:	NAPRT1; PP3856
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327910 representing NM_145201.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGCGGAGCAGGACCCCGAGGCGCGCGCGCGCGCGCGCTGCTCACTGACCTCTACCAGGCC
ACCATGGCGTTGGCTATTGGCGCGCGGGCCGGGCGGGACCGCGCGAGTTCGAGCTCTTCTCCGC
CGCTGCCCGTTCCGGCGCGCCTTCGCCTTGGCCCGCGGCTTGCAGACTGTGTGCCCTTCCGCGCGCC
TTCGCGCTGCGGGACCGCGACGTGCAGTTCCTGGCCTCGGTGCTGCCCCAGACACGGATCCTGCGTTC
TTCGAGCACCTTCGGGCCCTCGACTGCTCCGAGGTGACGGTGCGAGCCCTGCCGAGGGCTCCCTCGCC
TCCCCGGAGTGCCGCTCCTGCAGGTGTCCGGGCCGCTCCTGGTGGTGCAGCTGCTGGAGACACCCTG
CTCTGCCTGGTCACTACGCCAGCCTGGTGGCCACCAACGCAGCGCGGCTTCGCTTGATCGCAGGGCCA
GAGAAGCGGCTGTAGAGATGGCCTGAGGCGGGCTCAGGGCCCCGATGGGGGCTGACAGCCTCCACC
TACAGTACCTGGGCGGCTTCGACAGCAGCAGCAACGTGTAGCGGGCCAGCTGCGAGGTGTGCCGTG
GCCGGGACCTGGCCACTCCTTCGTCACTTCTTTTCAGGCAGCGAGGTGCCCCCTGACCCGATGTTG
GCCCAGCAGCTGGTGGGGCCCTGGGGTGGACCTGGCGGCCAAAGCCAGGTGTGGCTGGAGCAGGTG
TGTGCCAACCTGGGGCTGGGGGTGCAGGAGCCGATCCAGGCGAGCGGGCAGCCTTTGTGGCCTATGCC
TTGGCTTTTCCCGGGCCTTCCAGGGCCTCCTGGACACTACAGCGTGTGGAGGAGTGGTCTCCCAAC
TTCCTAGCAGTCGCCCTGGCCCTGGGAGAGCTGGGCTACCGGGCAGTGGGCGTGGAGCTGGACAGTGGT
GACCTGCTACAGCAGGCTCAGGAGATCCGCAAGTCTTCCGAGCTGTGCAGCCAGTTCAGGTGCC
TGGCTGGAGTCAGTCTCATCGTAGTCAGCAACAACATTGACGAGGAGGCGCTGGCCGACTGGCCAG
GAGGGCAGTGAGGTGAATGTATTGGCATTGGCACCAGTGTGGTCACTGCCCAACAGCCTTCCCTG
GGTGGCGTCTAAGCTGGTGGCCGTGGGGGCCAGCCACGAATGAAGCTGACCGAGGACCCCGAGAAG
CAGACGTTGCCTGGGAGCAAGGCTGCTTCCGGCTCCTGGGCTCTGACGGGTCTCCACTCATGGACATG
CTGCAGTTAGCAGAAGAGCCAGTGCCACAGGCTGGGCAGGAGCTGAGGGTGTGGCCTCCAGGGGCCAG
GAGCCCTGCACCGTGAGGCCAGCCAGGTGGAGCCACTACTGCGGCTCTGCCTCCAGCAGGGACAGCTG
TGTGAGCCGCTCCATCCCTGGCAGAGTCTAGAGCCTTGGCCAGCTGTCCCTGAGCCGACTCAGCCCT
GAGCACAGGCGGCTGCGGAGCCCTGCACAGTACCAGGTGGTGTGTCCGAGAGGCTGCAGGCCCTGGT
AACAGTCTGTGTGCGGGCAGTCCCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_145201
Insert Size:	1617 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_145201.5</u>
RefSeq Size:	1723 bp
RefSeq ORF:	1617 bp
Locus ID:	93100
UniProt ID:	<u>Q6XQN6</u>
Cytogenetics:	8q24.3
MW:	57.6 kDa
Gene Summary:	<p>Nicotinic acid (NA; niacin) is converted by nicotinic acid phosphoribosyltransferase (NAPRT; EC 2.4.2.11) to NA mononucleotide (NaMN), which is then converted to NA adenine dinucleotide (NaAD), and finally to nicotinamide adenine dinucleotide (NAD), which serves as a coenzyme in cellular redox reactions and is an essential component of a variety of processes in cellular metabolism including response to stress (Hara et al., 2007).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).</p>