

Product datasheet for **SC327913**

SLC39A5 (NM_173596) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC39A5 (NM_173596) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC39A5
Synonyms:	LZT-Hs7; MYP24; ZIP5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327913 representing NM_173596.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATGGGGTCCCAGTGAATCATCTGCTGGCCGGCTTCTGTGTGTGGTCTCTGGGCTGGGTAGGG
GGCTCAGTCCCCAACCTGGGCCCTGCTGACGAGGAGCAGAACCATTACCTGGCCAGCTGTTTGGCCTG
TACGGCGAGAATGGGACGCTGACTGACAGGGGGCTTGGCGGGCTTCTCCACAGCCTGGGGCTAGGCCGA
GTTCAGGGGCTTCGCTGGGACAGCATGGGCTCTGACTGGACGGGCTGCATCCCCAGCTGCAGACAAT
TCCACACACAGGCCACAGAACCCTGAGCTGAGTGTGGATGTCTGGGACAGGGATGCCTCTGGGTCCCTCA
GGGTGGGGTGAACGGAAGAGTCAAAGGCCCTCACCTACCCCGTGGGCCAGCCCCCTCGGGCTGGAC
CTCCTTACAGGCTTCTGTTGCTGGACCACTATTGGCTGACCACCTGAATGAGGATTGTCTGAACGGC
TCCCAGCTGCTGGTCAATTTGGCTTGGCCCCGCTGCTCCTGACCCCTCGTCAGTTTGTCTGCTG
TGCCAGCCCTGCTTATCAGATCGACAGCCGCTGTCATCGGCGCTCCGGCCCTGCACCCCAAGGG
GATCTACTATCTGCCCTGCTTACAGTGCCTGGCAGTCTGTTGCTCAGCCTCCCTTCTCCCTATCC
CTGCTGCTGCTGGGCTCCTGGGACCTGCTACTACGGCCCTTCTGGGCTTCTGGGGCCCTGGCG
GTGGGCACTCTTTGTGGGGATGCACTGCTACATCTGCTACCGCATGCACAAGAAGGGCGGCACGCAGGA
CCTGGCGGACTACCAGAGAAGGACTGGGCCCGGGGCTGTGAGTCTCGGAGGCCTTCTCTGCTCTTT
GTGCTGGAGAACATGCTGGGGCTTTTGGCGCACCGAGGGCTCAGGCCAAGATGCTGCAGGCGAAAACGA
AGGAATCTCGAAACACGCAACTTGGATCCGGAGAATGGCAGTGGGATGGCCCTTACAGCCCTACAGGCA
GCTCCAGAGCCAGGGGCTCAGGGCCAGAGGGAGAAGAACAGCCAGCACCACAGCTCTGGCCCCCTCCT
GGGCACCAAGGCCACAGTCATGGGACCCAGGGTGGCACTGATATCACGTGGATGGTCTCCTGGGAGAT
GGTCTACACAACCTCACTGATGGGCTGGCCATAGGTGCTGCCTTCTGATGGCTTCTCCAGCGGCCTC
AGTACCACCTTAGCGGTCTTCTGCCATGAGCTGCCCCACGAACTGGGTGACTTTGCCATGCTGCTCCAG
TCAGGGCTGTCTTTTGGCGGCTGCTGCTGCTGAGCCTCGTGTCTGGAGCCCTGGGATTGGGGGTGCA
GTCCTGGGGTGGGGCTCAGCCTGGGCCCTGTCCCCCTCACTCCCTGGGTGTTTGGGGTCACTGCTGGG
GTCTTCTCTATGTGGCCCTTGTGGACATGCTACCAGCCCTGCTTCTGCTCCTCCGGAGCCCTGCCTACG
CCCCATGTGCTCCTGCAGGGGCTGGGGCTGCTGCTGGGGGGCGGCCTCATGCTTGCATAACCCTGCTG
GAGGAGCGGCTACTGCCCGTACCCTGAGGGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_173596
Insert Size:	1623 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_173596.2
RefSeq Size:	2015 bp
RefSeq ORF:	1623 bp
Locus ID:	283375
UniProt ID:	Q6ZMH5
Cytogenetics:	12q13.3
Protein Families:	Transmembrane
MW:	56.5 kDa
Gene Summary:	<p>The protein encoded by this gene belongs to the ZIP family of zinc transporters that transport zinc into cells from outside, and play a crucial role in controlling intracellular zinc levels. Zinc is an essential cofactor for many enzymes and proteins involved in gene transcription, growth, development and differentiation. Mutations in this gene have been associated with autosomal dominant high myopia (MYP24). Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Sep 2014]</p> <p>Transcript Variant: This variant (1) represents the longer transcript. Both variants 1 and 2 encode the same protein.</p>